


**SUPPLIES & ACCESSORIES**  
ECG LEAD WIRES & TELEMETRY  
AAMI Color Coding



### ECG Combiner Lead Wires & Cables

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and disposable lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

#### ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



##### Disposable, DIN Shielded, Multi-Pinch, AAMI

700-0006-30	3 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-32	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-33	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets



##### Reusable, DIN Shielded, Multi-Pinch, AAMI

700-0006-35	3 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-37	5 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-38	5 Lead (V-Lead)	Multi-Pinch	74cm /29 in	1 Set

Multi-pinch

700-0010-00	3 Lead	Multi-Pinch	130cm/51 in	1 Set
700-0010-01	5 Lead	Multi-Pinch	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set



Snap

##### Reusable, DIN Shielded, Snap, AAMI

700-0007-35	3 Lead	Snap	74cm /29 in	1 Set
700-0007-37	5 Lead	Snap	74cm /29 in	1 Set
700-0010-18	3-Lead	Snap	130cm /51 in	1 Set
700-0010-19	5-Lead	Snap	130cm /51 in	1 Set

#### ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter



##### Radiolucent, Reusable, DIN Unshielded, Multi-Pinch, AAMI

700-0010-03	3 Lead, Radiolucent	Multi-Pinch	74cm /29 in	1 Set
700-0010-04	5 Lead, Radiolucent	Multi-Pinch	74cm /29 in	1 Set
700-0010-05	3 Lead, Radiolucent	Multi-Pinch	122cm/48 in	1 Set
700-0010-06	5 Lead, Radiolucent	Multi-Pinch	122cm/48 in	1 Set
700-0010-07	3 Lead, Radiolucent	Multi-Pinch	183cm/72 in	1 Set
700-0010-08	5 Lead, Radiolucent	Multi-Pinch	183cm/72 in	1 Set

## ECG Combiner Lead Wires & Cables

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and disposable lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Combiner Cables

#### Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496



700-0008-55	3 Lead	Spacelabs	305 cm/10 ft	1 Cable
700-0008-57	5 Lead	Spacelabs	305 cm/10 ft	1 Cable
700-0008-32	5 Lead	Spacelabs	457 cm/15 ft	1 Cable
700-0008-56	10 Lead (5+5)	Spacelabs	305 cm/10 ft	1 Cable

### ECG Combiner Adapters

#### Spacelabs Telemetry, AAMI



700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)		1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)		1 Adapter

## ECG Combiner Lead Wires & Cables

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and disposable lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Combiner Cables for Standardization

#### Other manufacturers

700-0008-35	6 lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable
700-0008-38	10 Lead (5+5)	GE	305 cm/10 ft	1 Cable
700-0008-39	6 Lead	GE	305 cm/10 ft	1 Cable
700-0008-47	3 Lead	Philips	305 cm/10 ft	1 Cable
700-0008-48	10 Lead (5+5)	Philips	305 cm/10 ft	1 Cable
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable
700-0008-43	6 lead	Medtronic LifePak	305 cm/10 ft	1 Cable
700-0008-52	3 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable
700-0008-54	6 lead	Zoll, 90° angle	305 cm/10 ft	1 Cable



### ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Shielded Lead Wire



91496

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

#### ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

Snap Sets			
700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16*	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18*	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

\*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220



Pinch

Pinch Sets			
700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set

\*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220

## ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Shielded Lead Wire



91496

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI



24" Snap Individual Color			
012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-34	Black	61 cm/24 in	5/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

### ECG Cables, Shielded, AAMI



With 17 PIN Adapter			
700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

## ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Shielded  
Lead Wire



91496

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Accessories

#### Electrodes



690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
-------------	--	-------	--------



690-0262-00	ECG electrode, transparent tape, general/stress, adult, single use	Adult	30/Box
-------------	--	-------	--------

### Neonatal ECG Supplies & Accessories

#### DIN Unshielded, Neonatal ECG Cable



700-0008-12	3 Lead	Neonatal	1 Cable
-------------	--------	----------	---------

#### DIN Unshielded, Neonatal Electrodes



685-0037-00	Pre-wired neonatal electrodes, Pre-gelled, silver/silver chloride, single use	3-Lead	100 Sets/Box
-------------	---	--------	--------------

#### Electrode



690-0290-00	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	30/Box
-------------	---	---------	--------

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap

#### Snap Sets

012-0283-11	3 Lead	61 cm/24 in	1 Set
012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0605-10	5 Lead	64 cm/25.2 in	1 Set



Pinch

#### Pinch Sets

012-0284-11	3 Lead	61 cm/24 in	1 Set
012-0286-11	5 Lead	61 cm/24 in	1 Set

### ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Snap

#### Snap Individual Color

012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack

#### Snap, Individual Color

*For use with Spacelabs Ultraview Telemetry 9034X/9134X only*



Snap

700-0010-12	Red	74 cm/29 in	5/Pack
700-0010-09	Black	74 cm/29 in	5/Pack
700-0010-10	Brown	74 cm/29 in	5/Pack
700-0010-13	White	74 cm/29 in	5/Pack
700-0010-11	Green	74 cm/29 in	5/Pack

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch

#### Pinch Individual Color

012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack

### Telemetry Transmitter Adapters



#### For use with ECG Combiner Lead Wires

700-0008-88	Aria Telemetry Combiner Adapter (96281)	1 Adapter
-------------	---	-----------



#### For use with DIN Unshielded Lead Wires

706-0191-00	Aria Telemetry Grouper (96281)	5/Bag
-------------	--------------------------------	-------

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Telemetry Accessories

#### Disposable Pouches and Electrodes



015-0500-02	Ultraview Telemetry pouch (90347)	100/Box	100/Box
690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
690-0262-00	ECG transparent electrode, general/stress	Adult	30/Box

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### SpO<sub>2</sub> Cable & Sensors



#### SpO<sub>2</sub> Cable

700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------



#### Reusable SpO<sub>2</sub> Sensors

015-0660-00	TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

#### Disposable SpO<sub>2</sub> Sensors



015-0662-00	TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO <sub>2</sub> vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------