

XPREZZON.

Introducing Xprezzon.<sup>®</sup>  
Connecting clinicians  
with patient information  
starts here.




# XPREZZON.



A powerful expression  
of performance and design.

Xprezzon provides easy access to patient data, creating information at the point of care for the clinician. The optional dual-display allows a dedicated full screen view of physiologic information plus a secondary screen for an alternative presentation of the patient's clinical data.



## STYLISH

A modern and innovative design provides an appealing and comfortable experience for both patients and their families. The sleek, frameless display makes cleaning easier to minimize the risk of infection.

## VERSATILE

Ergonomic and versatile mounting solutions enhance the overall appearance of the care environment. Various solutions for both monitor and cable management ensure a more organized and efficient workflow at the bedside.

---

## SENSITIVE

Xprezzon is sensitive to a patient's need for a good night's rest. It intuitively dims the display in low ambient light.

## TRENDY

A high-resolution screen displays crisp and visually rich Custom Trends which enable faster and better informed patient assessment. With a single touch, clinicians can access different trend views based on observation needs, type of procedure, or clinical protocols. Trends can be displayed continuously for real-time updates of critical patient information. They can also be sized and arranged as preferred, providing clinicians with their own uniquely customized views.



## ANYWHERE

Your patient's information in the palm of your hand. Whether inside the hospital or on the go, clinicians can stay connected with their patients. Spacelabs' innovative ICS provides the ultimate freedom and confidence - access to patient information on a clinician's mobile device.



## ENLIGHTENED

The sleek, frameless display integrates highly visible LED alarm lights. These illuminate in the front and rear, enhancing the clinician's ability to identify which monitor is alarming in the busy care environment.



## CONNECTED

Xprezzon's connectivity protects our customers' investments in previous generations of Spacelabs patient monitoring and information systems. Xprezzon's backward compatibility presents a seamless transition between platforms.

---

# SHAPING THE FUTURE OF PATIENT MONITORING

## A Quick Look at the Xprezzon High Acuity Patient Monitor Highly Visible, Easy to Use, Heightened Connectivity

- 3, 5, or Diagnostic 12 lead ECG, Resp, SPO2, NIBP, Temp
- 4 invasive (expandable to 8) and Cardiac Output
- Up to 8 waveforms, 24 parameters
- Highly visible 19" display and integrated LED alarm lights (Front and Rear)
- Secondary 17" (43.18 cm) or optional 19" (48.3 cm) interactive touchscreen display
- Full Bed Review of Remote Bedside Monitors
- Vital Signs Calculations
- EMR connectivity with ICS
- 24 hour Patient Data Shuttle
- USB Bar Code Scanner Ready
- Varitrend® 4 OCRG
- Optional Exergen USB Spot Check Thermometer
- Optional Printer/Recorder, Capnography, Multi-Gas, BIS, SvO2 and Dual SpO2 Modules
- LifeSync compatible

Features may vary depending on purchased configuration.

