

SUPPLIES & ACCESSORIES TEMPERATURE



Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

Continuous Temperature

Cables

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



| | | | |
|-------------|--|--------------|--------|
| 700-0031-00 | TruLink Dual Temperature Cable | 30 cm/12 in. | 1 Each |
| 700-4000-00 | Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors | 305 cm/10 ft | 1 Each |

Reusable Probes/Sensors

Skin Surface



| | | | |
|-------------|---|--------------|--------|
| 690-0028-01 | YSI 700 series compatible with 6.5 mm DIA phone jack type connections | 305 cm/10 ft | 1 Each |
|-------------|---|--------------|--------|

Esophageal / Rectal








| | | | |
|-------------|--------------------------------------|--------------|--------|
| 690-0272-00 | 12 French, YSI 400 series compatible | 305 cm/10 ft | 1 Each |
| 118-0256-01 | 12 French, YSI 700 series compatible | 305 cm/10 ft | 1 Each |

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



| Order Code | Description | Length/Size | Qty |
|---|--------------------------------------|-------------|--------|
| Single use Probes/Sensors - Requires YSI 400 series cable (700-4000-00) | | | |
| Skin | | | |
|  701-4000-00 | Skin Temperature Sensor | | 20/Box |
| Esophageal/ Stethoscope | | | |
|  701-4009-00 | 9 French, YSI 400 series compatible | | 20/Box |
| 701-4012-00 | 12 French, YSI 400 series compatible | | 20/Box |
| 701-4018-00 | 18 French, YSI 400 series compatible | | 20/Box |
| Esophageal/ Rectal | | | |
|  701-4109-00 | 9 French, YSI 400 series compatible | | 20/Box |
| 701-4112-00 | 12 French, YSI 400 series compatible | | 20/Box |
| Foley catheter with temperature probe | | | |
|  701-4214-00 | 14 French, YSI 400 series compatible | | 20/Box |
| 701-4216-00 | 16 French, YSI 400 series compatible | | 20/Box |
| 701-4218-00 | 18 French, YSI 400 series compatible | | 20/Box |
| Tympanic temperature sensor | | | |
|  701-4300-00 | Adult, YSI 400 series compatible | | 20/Box |