

SUPPLIES & ACCESSORIES


ELECTROCARDIOGRAPHY

IEC Color Coding



ECG Combiner Lead Wires & Cables

Hospitals are equipped with numerous monitors and other devices that use ECG lead wires, each with their own unique connectors which complicates patient transfers, undermines efficiency, and makes purchasing and maintenance more complex. Our ECG combiner lead wire system can standardize most equipment to one lead wire. The reusable or disposable lead wire sets are available in a variety of leads and lengths with either a snap or a multi-pin electrode interface. Validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry systems.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Multi-pin

Disposable, DIN Shielded, Multi-Pinch, IEC

700-0006-44	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-53	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets



Multi-pin

Reusable, DIN Shielded, Multi-Pinch, IEC

700-0006-41	3 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set



Snap

Reusable, DIN Shielded, Snap, IEC

700-0007-41	3 Lead	Snap	74cm /29 in	1 Set
700-0007-44	5 Lead	Snap	74cm /29 in	1 Set

Indicates discontinued part, available while supplies last

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)

700-0008-81	3 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------



700-0008-83	5 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------

ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



91496

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, Shielded, Reusable, IEC



Snap

Snap Sets

700-0007-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0007-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0007-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0007-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0007-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0007-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set



Pinch

Pinch Sets

700-0006-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0006-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0006-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0006-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0006-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0006-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set

ECG Lead Wires, Individual Colors, Shielded, Reusable, IEC

ECG Cables, Shielded, IEC



With 17 PIN Adapter

700-0008-09	3 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-07	5 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-01	10 Lead, shielded	305 cm/10 ft	1 Cable

Indicates discontinued part, available while supplies last

ECG Accessories



Electrodes

690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
-------------	--	-------	--------

ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



91496

Order Code	Description	Length	Qty
------------	-------------	--------	-----

Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable

700-0008-13	3 Lead		1 Cable
-------------	--------	--	---------



DIN Unshielded, Neonatal Electrodes

717-0080-00	Pre-wired neonatal electrodes, pre-gelled, silver/silver chloride, CE	3-Lead	100 Sets/Box
-------------	---	--------	--------------

Electrode

690-0290-00	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	30/Box
-------------	---	---------	--------




Indicates discontinued part, available while supplies last

DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded

Order Code	Description	Length/Size	Qty
ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC			
	Snap Sets		
012-0328-10	5 Lead	61 cm/24 in	1 Set
	Pinch Sets		
012-0325-10	3 Lead	61 cm/24 in	1 Set
012-0327-10	5 Lead	61 cm/24 in	1 Set
	For use with DIN Unshielded Lead Wires		
706-0190-00	Aria Telemetry Grouper (96281), IEC		5/Bag

Telemetry Accessories



Disposable Pouches and Electrodes			
015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box

Indicates discontinued part, available while supplies last

DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



DIN Un-Shielded

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

SpO₂ Cable & Sensors

SpO₂ Cable



700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

Note: Only the below genuine Spacelabs and Nellcor branded SpO₂ sensors are compatible with Spacelabs Telemetry systems

Reusable SpO₂ Sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

Disposable SpO₂ Sensors



015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last