


# SUPPLIES & ACCESSORIES








## INVASIVE BLOOD PRESSURE



## Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description	Length/Size	Qty
<b>Integrated Adapter Cables (connects directly to transducer)</b>			
<b>Edward Lifesciences</b>			
 700-0293-01	Single pressure integrated cable	 396 cm/13 ft	1 Each
700-0292-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
<b>Abbott Laboratories/Hospira/ICU Medical</b>			
 700-0295-01	Single pressure integrated cable	 396 cm/13 ft	1 Each
700-0294-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
<b>Interface Cable &amp; Adapters</b>			
<b>Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)</b>			
690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
<b>Interface Adapter- Transducer interface cable required from transducer company</b>			
 700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
<b>Transducer Kits &amp; Accessories</b>			
<b>Pressure Transducer Kits, Single use</b>			
 690-0008-00	Transducer Kit, Patient Mount, DTX/Plus	152.5 cm/60 in	5/Box
 690-0009-00	Transducer Kit, Pole Mount, DTX/Plus	152.5 cm/60 in	5/Box