


## 3-lead Cable Colors and Lengths – 3 Lengths

700-0378-00 Short  
700-0373-00 Medium  
700-0372-00 Long



Snap	Color	Short		Medium		Long	
		mm	Inch	mm	Inch	mm	Inch
C5	Red	205	8	510	20	865	34
M	Yellow	255	10	560	22	915	36
C5R	Green	205	8	510	20	865	34

## 4-lead Cable Colors and Lengths – 3 Lengths

700-0379-00 Short  
700-0375-00 Medium  
700-0374-00 Long



Snap	Color	Short		Medium		Long	
		mm	Inch	mm	Inch	mm	Inch
C5	Red	230	9	510	20	865	34
C3	Orange	255	10	560	22	915	36
M	White	255	10	560	22	915	36
C1	Brown	230	9	510	20	865	34

## IEC 12-lead Cable Colors and Lengths - 1 Length

700-0376-00



IEC Snap	Color	Length	
		mm	Inch
L	Yellow	610	24
F	Green	585	23
C6	Violet	635	25
C5	Black	535	21
C4	Brown	460	18
C3	Green	355	14
C2	Yellow	330	13
C1	Red	330	13
N	Black	380	15
R	Red	610	24

## AHA 12-lead Cable Colors and Lengths – 1 Length

700-0377-00



AHA Snap	Color	Length	
		mm	Inch
LA	Black	610	24
LL	Red	585	23
V6	Violet	635	25
V5	Orange	535	21
V4	Blue	460	18
V3	Green	355	14
V2	Yellow	330	13
V1	Red	330	13
RL	Green	380	15
RA	White	610	24

All dimensions and colors are approximate.

35301 SE Center St, Snoqualmie, WA 98065 | T +1 425 396 3300 | F +1 425 396 3301  
 Unit B, Foxholes Centre, John Tate Road, Hertford, SG13 7DT, UK | T +44 (0)1992 507700 | F +44 (0)1992 501213