

SUPPLIES & ACCESSORIES

CAPNOGRAPHY & MULTIGAS



NomoLine Sampling Lines

NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|



NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients



| NomoLine LH Nasal CO2 Cannula | | | |
|-------------------------------|----------------|-----------|--------|
| 704-3817-00 | Adult | 2 m/7 ft | 25/Box |
| 704-3819-00 | Infant | 2 m/7 ft | 25/Box |
| 704-3820-00 | Adult w/O2 | 2 m/7 ft | 25/Box |
| 704-4624-00 | Adult w/O2 | 4 m/13 ft | 25/Box |
| 704-3821-00 | Pediatric w/O2 | 2 m/7 ft | 25/Box |



| NomoLine LH Nasal/Oral CO2 Cannula | | | |
|------------------------------------|----------------|-----------|--------|
| 704-3824-00 | Adult w/O2 | 2 m/7 ft | 25/Box |
| 704-4628-00 | Adult w/O2 | 4 m/13 ft | 25/Box |
| 704-3825-00 | Pediatric w/O2 | 2 m/7 ft | 25/Box |



| NomoLine LH Single Nasal Prong CO2 Cannula | | | |
|--|-------|----------|--------|
| 704-3826-00 | Adult | 2 m/7 ft | 25/Box |

NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



| NomoLine LH Airway Adapter Set | | | |
|--------------------------------|-----------------|-----------|--------|
| 704-3814-00 | Adult/Pediatric | 2 m/7 ft | 25/Box |
| 704-4620-00 | Adult/Pediatric | 4 m/13 ft | 25/Box |
| 704-4369-00 | Infant/Neonatal | 2 m/7 ft | 25/Box |
| 704-4621-00 | Infant/Neonatal | 4 m/13 ft | 25/Box |

Indicates discontinued part, available while supplies last

NomoLine Sampling Lines

NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|



NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

NomoLine HH Nasal CO2 Cannula



| | | | |
|-------------|----------------|----------|--------|
| 704-3830-00 | Adult | 2 m/7 ft | 25/Box |
| 704-3832-00 | Infant | 2 m/7 ft | 25/Box |
| 704-3833-00 | Adult w/O2 | 2 m/7 ft | 25/Box |
| 704-3834-00 | Pediatric w/O2 | 2 m/7 ft | 25/Box |



NomoLine HH Nasal/Oral CO2 Cannula


| | | | |
|-------------|----------------|-----------|--------|
| 704-3837-00 | Adult w/O2 | 2 m/7 ft | 25/Box |
| 704-3838-00 | Pediatric w/O2 | 3 m/10 ft | 25/Box |



NomoLine HH Single Nasal Prong CO2 Cannula

| | | | |
|-------------|-------|----------|--------|
| 704-3839-00 | Adult | 2 m/7 ft | 25/Box |
|-------------|-------|----------|--------|

NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine HH Airway Adapter Set

| | | | |
|-------------|-----------------|-----------|--------|
| 704-3827-00 | Adult/Pediatric | 2 m/7 ft | 25/Box |
| 704-3828-00 | Adult/Pediatric | 3 m/10 ft | 25/Box |
| 704-3829-00 | Infant | 2 m/7 ft | 25/Box |
| 704-4367-00 | Infant/Neonatal | 2 m/7 ft | 25/Box |

Indicates discontinued part, available while supplies last

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

Sidestream Monitoring Supplies

Please note: the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.

CO2 Nasal Cannula with male luer connector



| | | | | |
|----|-------------|-----------|----------|--------|
| * | 704-0176-00 | Adult | 2 m/7 ft | 25/Box |
| CE | 704-0176-01 | Adult | 2 m/7 ft | 25/Box |
| ** | 704-0177-00 | Pediatric | 2 m/7 ft | 25/Box |
| ** | 704-0178-00 | Infant | 2 m/7 ft | 25/Box |

CO2 Nasal Cannula w/O2 Delivery and male luer connector



| | | | | |
|----|-------------|-----------|-----------|--------|
| * | 704-0181-00 | Adult | 2 m/7 ft | 25/Box |
| CE | 704-0181-01 | Adult | 2 m/7 ft | 25/Box |
| * | 704-0180-00 | Pediatric | 2 m/7 ft | 25/Box |
| CE | 704-0180-01 | Pediatric | 2 m/7 ft | 25/Box |
| ** | 704-0179-00 | Infant | 2 m/7 ft | 25/Box |
| ** | 704-0186-00 | Adult | 3 m/10 ft | 25/Box |

CO2 Oral/Nasal Cannula w/O2 Delivery and male luer connector



| | | | | |
|----|-------------|-------|-----------|--------|
| * | 704-0183-00 | Adult | 2 m/7 ft | 25/Box |
| CE | 704-0183-01 | Adult | 2 m/7 ft | 25/Box |
| ** | 704-0187-00 | Adult | 3 m/10 ft | 25/Box |

Indicates discontinued part, available while supplies last

Sidestream Monitoring Supplies - Intubated patients

Sampling Lines & Airway Adapters, single use



| | | | | |
|----|-------------|------------------------------------|-----------|--------|
| ** | 704-0191-00 | Adult CO2 Sample Line w/ T-Adapter | 2 m/7 ft | 25/Box |
| ** | 704-0192-00 | Adult CO2 Sample Line w/ T-Adapter | 3 m/10 ft | 25/Box |



| | | | | |
|----|-------------|---|-----------|--------|
| * | 166-7121-00 | Sidestream sampling line with male luer lock on both ends | 3 m/10 ft | 50/Box |
| CE | 166-7121-02 | Sidestream sampling line with male luer lock on both ends | 3 m/10 ft | 50/Box |



| | | | | |
|---|-------------|--------------------------------------|--|--------|
| * | 704-0175-00 | Sidestream airway adapter/T-Adapter. | | 25/Box |
|---|-------------|--------------------------------------|--|--------|

* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients



Sensor & Airway Adapters - required for mainstream monitoring



| | | | |
|-------------|---------------------------------------|-----------------|--------|
| 010-1980-00 | Mainstream sensor, reusable | | 1 Each |
| 704-0173-00 | Mainstream airway adapter, single use | Adult/Pediatric | 25/Box |
| 704-0174-00 | Mainstream airway adapter, single use | Infant | 10/Box |

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module

| Order Code | Description | Size | Qty |
|---|--|-----------|--------|
| Sidestream, Non-intubated, single patient use | | | |
| CO₂ Nasal Cannula | | | |
|  * 704-0010-00 | Adult | | 10/Box |
| CE 704-0010-10 | Adult | | 10/Box |
| * 704-0012-00 | Pediatric | | 10/Box |
| CO₂ Nasal Cannula w/O₂ Delivery | | | |
| * 704-0014-00 | Adult | | 10/Box |
| CE 704-0014-10 | Adult, (for sale into EU/UK only) | | 10/Box |
| * 704-0016-00 | Pediatric | | 10/Box |
| CE 704-0016-10 | Pediatric, (for sale into EU/UK only) | | 10/Box |
| * 704-0015-00 | Infant/Neonatal | | 10/Box |
| CO₂ Oral/Nasal Cannula | | | |
| * 704-0017-00 | Adult | | 10/Box |
| CO₂ Oral/Nasal Cannula w/O₂ Delivery | | | |
| * 704-0020-00 | Adult | | 10/Box |
| CE 704-0020-10 | Adult, (for sale into EU/UK only) | | 10/Box |
| * 704-0022-00 | Pediatric | | 10/Box |
| CE 704-0022-10 | Pediatric, (for sale into EU/UK only) | | 10/Box |
| Sidestream, Intubated, single patient use | | | |
| Sidestream - short term up to 12 hours (patient non-humidified) | | | |
|  * 704-0023-00 | Airway adapter + sampling line, ET tube > 4.0 mm | Adult/Ped | 10/Box |
| CE 704-0023-10 | Airway adapter + sampling line, ET tube > 4.0 mm | Adult/Ped | 10/Box |
| * 704-0027-00 | CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i> | | 10/Box |
| CE 704-0027-10 | CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i> | | 10/Box |

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
 Dual Capnography Module

| Order Code | Description | Size | Qty |
|--|---|------------|--------|
| Sidestream - long term up to 120 hours (patient humidified) | | | |
| * 704-0025-00 | Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm | Adult/Ped | 10/Box |
| CE 704-0025-10 | Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm | Adult/Ped | 10/Box |
| * 704-0026-00 | Airway adapter + sampling line with Nafion | Ped/Infant | 10/Box |
| * 704-0028-00 | CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i> | | 10/Box |
| CE 704-0028-10 | CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i> | | 10/Box |

Mainstream, for Intubated Patients Only

Mainstream Airway Adapters, for use with CO₂ Sensor 704-0001-00

| | | | |
|----------------|--|---------|--------|
| * 704-0005-00 | CO ₂ airway adapter, single use | Adult | 24/Box |
| CE 704-0005-10 | CO ₂ airway adapter, single use | Adult | 10/Box |
| * 704-0004-00 | CO ₂ airway adapter, single use | Neonate | 24/Box |
| CE 704-0004-10 | CO ₂ airway adapter, single use | Neonate | 10/Box |
| * 704-0002-00 | CO ₂ airway adapter, reusable | Adult | 1 Each |
| CE 704-0002-10 | CO ₂ airway adapter, reusable | Adult | 1 Each |
| CE 704-0003-10 | CO ₂ airway adapter, reusable | Neonate | 1 Each |

* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)




(91518)

| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

Multigas Supplies & Accessories

Anesthetic Agent Samples Lines



| | | | |
|-------------|---|--|--------|
| 015-0128-00 | SAMPLE LINE,MALE TO MALE,LUER,SPU,7'X.06",25/BX | | 25/Box |
|-------------|---|--|--------|

| | | | |
|-------------|--|--------------|--------|
| 015-0312-01 | Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518) | 213 cm/7 ft. | 50/Box |
|-------------|--|--------------|--------|



| | | | |
|-------------|--|---------------|--------|
| 015-0313-01 | Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518) | 305 cm/10 ft. | 50/Box |
|-------------|--|---------------|--------|



| | | | |
|-------------|--|---------------|--------|
| 162-0046-01 | Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518) | 305 cm/10 ft. | 50/Box |
|-------------|--|---------------|--------|



| | | | |
|-------------|--|--------------|--------|
| 162-0047-01 | Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518) | 213 cm/7 ft. | 50/Box |
|-------------|--|--------------|--------|

Indicates discontinued part, available while supplies last

Water Traps, Reusable



| | | | |
|-------------|--|--|---------|
| 717-0447-00 | D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518) | | 10/Pack |
|-------------|--|--|---------|