





Qube<sup>®</sup> facilitates at-a-glance monitoring with an ultra-bright display featuring large numerics. Multicolor, severity-based lights illuminate on both the front and back of the monitor for enhanced alarm surveillance.





## Inspired by simplicity and style, Qube breaks with tradition to offer a fresh, new perspective in compact monitoring.

### INTRODUCING QUBE®

Small and lightweight with a battery life that goes the distance, Qube offers portability, accessibility and heightened connectivity—perfect for emergency-, general-, intensive- and post-anesthesia care units. Caregivers can access current, critical patient information— from bedside throughout transport—to provide the ultimate level of patient care.

### VERSATILE

Qube is mountable anywhere – wall, roll-stand, bed rail, anesthesia system – or can quickly be detached via a quick-release feature for immediate mobility. Long battery life and dual battery capability keep Qube on the go everywhere your patient goes. Impact resistant, chemical resistant, and protected from liquid ingress, Qube endures even with heavy use. Full functionality in a compact package.





### INSIGHTFUL

Data from a variety of sources informs and supports clinical determinations at the bedside or remotely. Graphical assessment tools and the ability to access up to 96 hours of physiologic data empower clinicians and help to facilitate decision making.

### CONFIDENT

Clinicians can customize the screen view to quickly access a multitude of graphical and tabular trend presentations based on individual patient needs and procedural guidelines. Interactive view, control, and recording for parameters, trends, and calculations can be engaged instantaneously, linking clinicians with a range of real-time patient physiologic data. Qube enables point-of-care decisions that provide ongoing reassurance and confidence to patients and their loved ones.




### ENLIGHTENED

Qube facilitates at-a-glance monitoring with an ultra-bright display featuring large numerics. Multicolor, severity-based lights illuminate on both the front and back of the monitor for enhanced alarm surveillance.

### ATTRACTIVE

The sleek, contemporary details of Qube, including ergonomic user controls and a comfortable grip handle, create a positive impact for caregiver and patient. Art meets technology.

The small, compact form of the capnography pod makes it ideal for use with Qube. Designed to attach seamlessly onto the rear of the monitor, it's ready in a snap.



### CONNECTED

With its modular architecture, Qube easily connects to any Spacelabs physiologic module or Flexport for continuity and ease when networking across different platforms. Qube also utilizes industry-standard communication protocols to allow inter-monitor communication, enhanced alarm surveillance, and broader productivity through both wired and wireless channels. Clinicians can choose how they stay connected - single or multiparameter view - for any monitored patient on the network.



## A Quick Look at the Qube Compact Patient Monitor. Portability, Accessibility, Heighted Connectivity.

- 3, 5, 12 lead ECG, Resp, SPO2, NIBP
- 4 invasive pressures and Cardiac Output
- 6 waveforms, 18 parameters
- Secondary 15" (38 cm) or optional 19" (48.3 cm) interactive touchscreen display
- Standard 8 hour battery using 2 hot swappable batteries
- 24 hour Patient Data Shuttle
- One button quick release docking station
- Wireless connectivity
- Internal strip recorder
- Optional capnography

Features may vary depending on purchased configuration.

