

# SUPPLIES & ACCESSORIES

## C50 VITAL SIGNS MONITOR













## Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
<b>ECG leadwire systems</b>				
 040-000908-00	CABLE,INTEGRATED 3-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000909-00	CABLE,INTEGRATED 5-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000910-00	CABLE,INTEGRATED 10-WIRE CLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000911-00	CABLE,INTEGRATED 3-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000912-00	CABLE,INTEGRATED 5-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000913-00	CABLE,INTEGRATED 10-WIRE CLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
<b>Pulse Oximetry Supplies</b>				
 040-000869-00	CABLE,INTEGRATED ADT FINGER CLIP SENSOR,C50	3 meters	Reusable	1 Each
 040-000769-00	CABLE,SPO2,12PIN TO DB9F,C50 <i>Sensors sold seperately</i>	3 meters	Reusable	1 Each
 040-000726-00	SENSOR,SPO2,ADT,FINGER CLIP,C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each
 040-000730-00	SENSOR,SPO2,ADT/NEO/PED,WRAP, 3',C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each









## Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
<b>Non-Invasive Blood Pressure supplies</b>				
 040-000808-00	ADULT NIBP TUBE (WHITE)2M	2 meters	Reusable	1 Each
	Reusable or Disposable Cuffs See NIBP section under TruLink Reusable or Disposable cuffs. These are compatible with the C50. <i>Requires use of the 040-00808-00 air hose (above).</i>			
 <b>714-0086-00</b>	<b>NIBP,HOSE ASSY,NEONATE,C50</b> Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	<b>2.5 meters</b>	<b>Reusable</b>	<b>1 Each</b>
	<b>Bold font indicates a new part number</b>			
<b>Capnography Supplies: requires external Phasein ISA module, 099-000007-00</b>				
 040-000400-00	<b>CABLE,CO2/AG,C50</b> Connects the Phasein ISA module to the C50	<b>1 meter</b>	<b>Reusable</b>	<b>1 Each</b>
	For airway adapter, cannula and sampling line sets, please see the capnography supply section under Nomoline SPU supplies. These are compatible with the C50.			
<b>Temperature sensors</b>				
 040-000246-00	<b>PROBE,TEMP,ADT,SKIN,10K,C50</b>	<b>3 meters</b>	<b>Reusable</b>	<b>1 Each</b>
 040-000387-00	PROBE,TEMP,NEO SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000385-00	PROBE,TEMP,ADT,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each
 040-000386-00	<b>PROBE,TEMP,NEO,RECTAL/ESOPHAGEAL,10K,C50</b>	<b>3 meters</b>	<b>Reusable</b>	<b>1 Each</b>









## Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
<b>Invasive Blood Pressure Cables and Transducer kits</b>				
	040-000454-00 INTERFACE CABLE,IBP,EDWARD CONNCTR,C50 <i>Requires 040-000211-00 transducers</i>	 3 meters	Reusable	1 Each
	040-000360-00 INTERFACE CABLE,IBP,ABBOTT CONNCTR,C50 <i>Requires 040-000870-00 transducers</i>	 3 meters	Reusable	1 Each
	040-000870-00 DISPOSABLE PRESSURE TRANSDUCER KIT(COPPER),QTY2 <i>Abbott compatible, requires 040-000360-00 cable</i>	sterile kit	Disposable	2 Each
	040-000456-00 INTERFACE CABLE,IBP,BD CONNCTR,C50 <i>Requires 040-000403-00 transducers</i>	 3 meters	Reusable	1 Each
<b>Recording Paper</b>				
	040-000116-00 RECORDER PAPER,THERMAL,C50	50mm X 20 meters		1 roll Each