

SUPPLIES & ACCESSORIES ECLIPSE MINI AND ECLIPSE PRO



Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.

Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code	Description	Length/Style	Qty																																														
Eclipse Pro																																																	
Leadwire Sets																																																	
700-0378-00	3 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>M</td> <td>Yellow</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C5R</td> <td>Green</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	205	8	510	20	865	34	M	Yellow	255	10	560	22	915	36	C5R	Green	205	8	510	20	865	34	Each								
Snap	Color				Short		Medium		Long																																								
			mm	Inch	mm	Inch	mm	Inch																																									
C5	Red	205	8	510	20	865	34																																										
M	Yellow	255	10	560	22	915	36																																										
C5R	Green	205	8	510	20	865	34																																										
700-0391-00	3 Lead snap on with yoke Medium		Each																																														
700-0390-00	3 Lead snap on with yoke Long		Each																																														
700-0379-00	4 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>C3</td> <td>Orange</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>M</td> <td>White</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C1</td> <td>Brown</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	230	9	510	20	865	34	C3	Orange	255	10	560	22	915	36	M	White	255	10	560	22	915	36	C1	Brown	230	9	510	20	865	34	Each
Snap	Color				Short		Medium		Long																																								
			mm	Inch	mm	Inch	mm	Inch																																									
C5	Red	230	9	510	20	865	34																																										
C3	Orange	255	10	560	22	915	36																																										
M	White	255	10	560	22	915	36																																										
C1	Brown	230	9	510	20	865	34																																										
700-0394-00	4 Lead snap on with yoke Medium		Each																																														
700-0393-00	4 Lead snap on with yoke Long		Each																																														
700-0396-00	10 Lead snap on with yoke <i>IEC color coding generally used in European Union</i>	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th colspan="2">Length</th> </tr> <tr> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>L</td> <td>610</td> <td>24</td> </tr> <tr> <td>F</td> <td>585</td> <td>23</td> </tr> <tr> <td>C6</td> <td>635</td> <td>25</td> </tr> <tr> <td>C5</td> <td>535</td> <td>21</td> </tr> </tbody> </table>	Snap	Length		mm	Inch	L	610	24	F	585	23	C6	635	25	C5	535	21	Each																													
Snap	Length																																																
	mm	Inch																																															
L	610	24																																															
F	585	23																																															
C6	635	25																																															
C5	535	21																																															
700-0377-00	10 Lead snap on with yoke <i>AHA color coding generally used in USA, outside of EU</i>	<table border="1"> <tbody> <tr> <td>C4</td> <td>460</td> <td>18</td> </tr> <tr> <td>C3</td> <td>355</td> <td>14</td> </tr> <tr> <td>C2</td> <td>330</td> <td>13</td> </tr> <tr> <td>C1</td> <td>330</td> <td>13</td> </tr> <tr> <td>N</td> <td>380</td> <td>15</td> </tr> <tr> <td>R</td> <td>610</td> <td>24</td> </tr> </tbody> </table>	C4	460	18	C3	355	14	C2	330	13	C1	330	13	N	380	15	R	610	24	Each																												
C4	460	18																																															
C3	355	14																																															
C2	330	13																																															
C1	330	13																																															
N	380	15																																															
R	610	24																																															
712-1052-01	ECLIPSE PRO ECG SENSORS <i>Not available for sale in USA</i>		Box of 10																																														
Accessories																																																	
706-0207-00	Eclipse Pro belt clip		Each																																														
073-0267-00	7 Patient Diary, English CE0123		Each																																														
392196-001	Adhesive tape remover pads, 100/box		100/Box																																														
717-0071-00	ECG electrode, Ambu Blue Sensor L		5/Pack																																														

Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.


Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

	Order Code	Description	Length/Style	Qty
Eclipse Mini				
	Eclipse Mini Supplies			
	712-1053-01	ECLIPSE MINI ECG SENSORS <i>Not available for sale in USA</i>		Box of 10
	712-1053-02	ECLIPSE MINI ECG SENSORS, DEFIB SAFE <i>Available for USA September 2023</i>		Box of 10