


SUPPLIES & ACCESSORIES

INVASIVE BLOOD PRESSURE



Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Integrated Adapter Cables (connects directly to transducer)

Edward Lifesciences



700-0293-01 Single pressure integrated cable

396 cm/13 ft

1 Each

015-0727-00 Single pressure integrated cable

400 cm/13 ft

1 Each



700-0292-01 Dual pressure integrated cable


396 cm/13 ft

1 Each

015-0728-00 Dual pressure integrated cable

400 cm/13 ft

1 Each



Abbott Laboratories/Hospira/ICU Medical



700-0295-01 Single pressure integrated cable

396 cm/13 ft

1 Each

015-0729-00 Single pressure integrated cable

400 cm/13 ft

1 Each



700-0294-01 Dual pressure integrated cable


396 cm/13 ft

1 Each

015-0730-00 Dual pressure integrated cable

400 cm/13 ft

1 Each



Bold font indicates a new part number

Indicates discontinued part, available while supplies last

Interface Cable & Adapters



Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)

690-0021-00 Requires adapter cable 700-0091-00 or 700-0028-00

366 cm/12 ft

1 Each



Interface Adapter- Transducer interface cable required from transducer company

700-0791-00 Single pressure cable

122 cm/4 ft

1 Each

700-0028-00 Dual pressure cable

122 cm/4 ft

1 Each

Transducer Kits & Accessories

Pressure Transducer Kits, Single use



690-0008-00 Transducer Kit, Patient Mount, DTX/Plus

152.5 cm/60 in

5/Box



690-0009-00 Transducer Kit, Pole Mount, DTX/Plus

152.5 cm/60 in

5/Box