

SUPPLIES & ACCESSORIES

PRESSURE INFUSORS



Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



| Order Code | Description | Size | Qty |
|------------|-------------|------|-----|
|------------|-------------|------|-----|

Unifusor Disposable Pressure Infusor

Single Patient Use Infusor with aneroid gauge - sold as one piece/complete unit



| | | | |
|-----------|--|--------|--------|
| 803FGA-24 | Unifusor with thumbwheel valve & aneroid gauge | 500mL | 24/Box |
| 804FGA-24 | Unifusor with thumbwheel valve & aneroid gauge | 1000mL | 24/Box |
| 805FGA-10 | Unifusor with thumbwheel valve & aneroid gauge | 3000mL | 10/Box |



| | | | |
|-----------|--|--------|--------|
| 803SGA-24 | Unifusor with stopcock valve & aneroid gauge | 500mL | 24/Box |
| 804SGA-24 | Unifusor with stopcock valve & aneroid gauge | 1000mL | 24/Box |
| 805SGA-10 | Unifusor with stopcock valve & aneroid gauge | 3000mL | 10/Box |

Disposable Infusor with piston gauge - sold as one piece/complete unit



| | | | |
|-----------|---|--------|--------|
| 903TGA-24 | Unifusor with thumbwheel valve & piston gauge | 500mL | 24/Box |
| 904TGA-24 | Unifusor with thumbwheel valve & piston gauge | 1000mL | 24/Box |
| 905TGA-10 | Unifusor with thumbwheel valve & piston gauge | 3000mL | 10/Box |



| | | | |
|-----------|---|--------|--------|
| 903SGA-24 | Unifusor with stopcock valve & piston gauge | 500mL | 24/Box |
| 904SGA-24 | Unifusor with stopcock valve & piston gauge | 1000mL | 24/Box |
| 905SGA-10 | Unifusor with stopcock valve & piston gauge | 3000mL | 10/Box |

Indicates discontinued part, available while supplies last

Unifusor II Disposable Pressure Infusor with Reusable Inflation System

Disposable Infusor - white Disposable bag with reusable inflation system



| | | | |
|------------|---|--------|--------|
| 2803SGA-24 | Unifusor II with aneroid gauge & stopcock valve | 500mL | 24/Box |
| 2804SGA-24 | Unifusor II with aneroid gauge & stopcock valve | 1000mL | 24/Box |
| 2805SGA-10 | Unifusor II with aneroid gauge & stopcock valve | 3000mL | 10/Box |



| | | | |
|------------|--|--------|--------|
| 2903SGA-24 | Unifusor II with piston gauge & stopcock valve | 500mL | 24/Box |
| 2904SGA-24 | Unifusor II with piston gauge & stopcock valve | 1000mL | 24/Box |
| 2905SGA-10 | Unifusor II with piston gauge & stopcock valve | 3000mL | 10/Box |



Disposable Replacement Parts

| | | | |
|----------|---|-------|--------|
| 2803X-24 | Unifusor II Disposable white infusion bag | 500mL | 24/Box |
|----------|---|-------|--------|

Unifusor® Pressure Infusors


Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



| Order Code | Description | Size | Qty |
|------------|-------------|------|-----|
|------------|-------------|------|-----|

Unifusor Plus Reusable Pressure Infusor

Reusable Infusor - blue bag with white hose



| | | | |
|------------|---|--------|--------|
| 1103SGA-24 | Unifusor Plus with aneroid gauge & stopcock valve | 500mL | 24/Box |
| 1104SGA-24 | Unifusor Plus with aneroid gauge & stopcock valve | 1000mL | 24/Box |
| 1105SGA-10 | Unifusor Plus with aneroid gauge & stopcock valve | 3000mL | 10/Box |
| 1113SGA-24 | Unifusor Plus with piston gauge & stopcock valve | 500mL | 24/Box |

Reusable replacement parts



| | | | |
|-----------|---|--------|-------|
| 1103XB-05 | Unifusor Plus blue reusable bag with white hose | 500mL | 5/Box |
| 1104XB-05 | Unifusor Plus blue reusable bag with white hose | 1000mL | 5/Box |
| 1105XB-05 | Unifusor Plus blue reusable bag with white hose | 3000mL | 5/Box |

Indicates discontinued part, available while supplies last

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with black hose

| | | | |
|------|---|--------|--------|
| 1103 | Unifusor Plus with trigger / gauge assembly | 500mL | 1 Each |
| 1104 | Unifusor Plus with trigger / gauge assembly | 1000mL | 1 Each |
| 1105 | Unifusor Plus with trigger / gauge assembly | 3000mL | 1 Each |

Reusable replacement parts



| | | | |
|-------|---|--------|--------|
| 1103X | Unifusor Plus blue reusable bag with black hose | 500mL | 1 Each |
| 1104X | Unifusor Plus blue reusable bag with black hose | 1000mL | 1 Each |
| 1105X | Unifusor Plus blue reusable bag with black hose | 3000mL | 1 Each |

| | | | |
|------|------------------------------|--|--------|
| 1199 | Trigger/ gauge assembly only | | 1 Each |
|------|------------------------------|--|--------|