

# SUPPLIES & ACCESSORIES PRODUCT CATALOG

September 2023

# SUPPLIES & ACCESSORIES

## NON-INVASIVE BLOOD PRESSURE CUFFS



# BLOOD PRESSURE CUFFS

## TruLink® Nylon

### TruLink® Nylon Reusable

TruLink reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. TruLink cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, TruLink cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. TruLink cuffs are naturally latex and PVC-free.



Range/  
Length Qty

Order Code Description Size

### TruLink Nylon Reusable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1111-10	Single tube, HP	Infant	8 - 14 cm	1 Each
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each



Bayonet

#### Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



Part Numbers have been changed

## TruLink® Soft Disposable

TruLink fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### TruLink Soft Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*



**Bayonet**

715-1113-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

#### Single Tube Hose



714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each

**Part Numbers have been changed**

# BLOOD PRESSURE CUFFS

## TruLink® Vinyl

### TruLink® Vinyl Disposable

TruLink vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### TruLink Vinyl Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*



Bayonet

715-1112-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1122-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



#### Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each

Part Numbers have been changed

## Blood Pressure Cuff Accessories

Spacelabs non-invasive blood pressure cuff accessories are compatible with all Spacelabs and TruLink branded non-invasive blood pressure cuffs. Choose from a variety of hoses, adapters and conversion accessories which offer ease-of-use and flexibility.

Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

### Non-Invasive Blood Pressure Hose

#### Single Tube Hose



*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

714-0018-02	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-02	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0019-02	90496/91496	Neonatal	275 cm/9 ft	1 Each

Part Numbers have been changed

#### Double Tube to Single Tube Hose



*For use with Command Modules 90496/91496*

714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each

### Hose Connectors & Adapters

#### Slip Lock Fitting



*Fitting for Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.*





714-0038-00	3mm ID Adult NIBP hoses	1 Each
714-0047-00	4mm ID Adult NIBP hoses	1 Each
714-0046-00	5mm ID Adult NIBP hoses	1 Each

Indicates discontinued part, available while supplies last

**UltraCheck® Nylon Reusable**

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.






	Order Code	Description	Size	Length	Qty
<b>UltraCheck Nylon Reusable BP Cuffs</b>					
<b>Single Tube</b>					
 <b>Plastic Bayonet</b>	US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
 <b>Metal Bayonet</b>	US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
	US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
	US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
	US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
	US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
	US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
	US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
	US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag
 <b>Threaded Screw</b>	US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
	US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
	US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
	US2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
	US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
	US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
 <b>Female Subminiature Twist Lock</b>	US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
	US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
	US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
	US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
	US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag

Indicates discontinued part, available while supplies last

**UltraCheck® Nylon Reusable**

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.








	Order Code	Description	Size	Length	Qty
<b>UltraCheck Nylon Reusable BP Cuffs</b>					
<b>Single Tube</b>					
 <b>Subminiature Twist Lock</b>	US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
	US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
	US2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
	US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
	US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
 <b>No Fitting</b>	US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
 <b>Manual Bulb</b>	US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
	US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
	US2938MB	Single tube, MB w/gauge	Adult Long	29 - 38 cm	1 Each
	US3544MB	Single tube, MB w/gauge	Large Adult Long	35 - 44 cm	1 Each
	US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each
Indicates discontinued part, available while supplies last					



**UltraCheck® Nylon Reusable**

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



	Order Code	Description	Size	Length	Qty
<b>Double Tube</b>					
 <b>Threaded Screw</b>	UD1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
	UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
	UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
 <b>Subminiature Twist Lock</b>	UD0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag
	UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
	UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
	UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag
 <b>Male &amp; Female Mated</b>	UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
	UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
	UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
	UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag
 <b>No Fitting</b>	UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
	UD1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
	UD2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
	UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
	UD4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag
 <b>Manual Bulb</b>	UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
	UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
	UD3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag

Indicates discontinued part, available while supplies last

**SoftCheck® Fabric Disposable**

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



**Range/  
Length**

**Qty**

**Order Code**

**Description**

**Size**

**SoftCheck Fabric Disposable BP Cuffs**

**Single Tube**

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
<b>ST0814HP-20</b>	<b>Single tube, HP</b>	<b>Infant</b>	<b>8 - 14 cm</b>	<b>20/Box</b>
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
<b>ST4250HP-20</b>	<b>Single tube, HP</b>	<b>Thigh</b>	<b>42 - 50 cm</b>	<b>20/Box</b>
<b>ST0814MQ-05</b>	<b>Single tube, MQ</b>	<b>Infant</b>	<b>8 - 14 cm</b>	<b>5/Bag</b>
<b>ST1320MQ-05</b>	<b>Single tube, MQ</b>	<b>Child</b>	<b>13 - 20 cm</b>	<b>5/Bag</b>
<b>ST1826MQ-05</b>	<b>Single tube, MQ</b>	<b>Small Adult</b>	<b>18 - 26 cm</b>	<b>5/Bag</b>
<b>ST2635MQ-05</b>	<b>Single tube, MQ</b>	<b>Adult</b>	<b>26 - 35 cm</b>	<b>5/Bag</b>
<b>ST2635MQ-20</b>	<b>Single tube, MQ</b>	<b>Adult</b>	<b>26 - 35 cm</b>	<b>20/Box</b>
<b>ST2938MQ-05</b>	<b>Single tube, MQ</b>	<b>Adult Long</b>	<b>29 - 38 cm</b>	<b>5/Bag</b>
<b>ST3242MQ-05</b>	<b>Single tube, MQ</b>	<b>Large Adult</b>	<b>32 - 42 cm</b>	<b>5/Bag</b>



**Plastic Bayonet**



**Female  
Subminiature**

**Indicates discontinued part, available while supplies last**

**SoftCheck® Fabric Disposable**

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



**Range/  
Length**

**Qty**

**Order Code**

**Description**

**Size**

**SoftCheck Fabric Single Patient Multiple Use BP Cuffs**



**No Fitting**

**Single Tube**

ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
ST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag



**Threaded Screw**

**Double Tube**

DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box



**Twist Lock**

DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box
DT2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag



**Manual Bulb**

DT2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Box
-------------	-----------------	-------	------------	-------

**No fitting**






DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
-----------	-------------------------	------------	------------	-------

Indicates discontinued part, available while supplies last

## SoftCheck® Vinyl Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



	Order Code	Description	Size	Range/ Length	Qty
<b>Single Tube</b>					
	VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
	VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
	VST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
	VST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
	VST0814MQML5-20	Single tube, MQ	Infant	8 - 14 cm	5/Bag
	VST1320MQML5-20	Single tube, MQ	Child	13 - 20 cm	5/Bag
	VST1826MQML5-20	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
	VST2635MQML5-20	Single tube, MQ	Adult	26 - 35 cm	5/Bag
	VST2938MQML5-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box

Indicates discontinued part, available while supplies last

# BLOOD PRESSURE CUFFS

## UltraCheck® & SoftCheck® Curve

### UltraCheck® Curve Nylon Reusable SoftCheck® Curve Fabric Disposable

Your patients come in all shapes and sizes. That's why we offer the Curve cuff for extra-large adults. Its conical-shaped design provides a better fit for the contours of your patient's arm. This can increase patient comfort while providing a more accurate reading.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

#### UltraCheck Curve Nylon Reusable BP Cuff

Single Tube				
BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MB	Single tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each
Double Tube				
BRUD3854	Double tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MB	Double tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

#### SoftCheck Curve Fabric Disposable BP Cuff

Single Tube				
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Box
BRST3854HP-05	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	5/Box

Indicates discontinued part, available while supplies last



# BLOOD PRESSURE CUFFS SoftCheck® Fabric Neonatal

## SoftCheck® Fabric Neonatal Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Fabric Neonatal Disposable BP Cuffs



Micro-taper  
slip-fit

Single Tube					
NN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box	
NN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box	
NN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box	
NN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box	
NN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box	
NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box	
NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box	
NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box	
NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box	
NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box	
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box	
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box	



Male Luer



Plastic Bayonet

Indicates discontinued part, available while supplies last



Male Luer

Double Tube					
NN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box	
NN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box	
NN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box	
NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box	








# BLOOD PRESSURE CUFFS

## SoftCheck® Vinyl Neonatal

### SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



	Order Code	Description	Size	Range/ Length	Qty
<b>SoftCheck Vinyl Neonatal Disposable BP Cuffs</b>					
<b>Single Tube</b>					
 <b>Micro-taper slip-fit</b>	VNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
 <b>SnapQuik</b>	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 <b>Male Luer</b>	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 <b>Female Twist Lock</b>	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box
<b>Double Tube</b>					
 <b>Male Luer</b>	VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
	VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
	VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
	VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
	VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

# SUPPLIES & ACCESSORIES

## PRESSURE INFUSORS





**Unifusor® Pressure Infusors**

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

**Unifusor Disposable Pressure Infusor**

**Single Patient Use Infusor with aneroid gauge - sold as one piece/complete unit**



803FGA-24	Unifusor with thumbwheel valve & aneroid gauge	500mL	24/Box
804FGA-24	Unifusor with thumbwheel valve & aneroid gauge	1000mL	24/Box
805FGA-10	Unifusor with thumbwheel valve & aneroid gauge	3000mL	10/Box



803SGA-24	Unifusor with stopcock valve & aneroid gauge	500mL	24/Box
804SGA-24	Unifusor with stopcock valve & aneroid gauge	1000mL	24/Box
805SGA-10	Unifusor with stopcock valve & aneroid gauge	3000mL	10/Box

**Disposable Infusor with piston gauge - sold as one piece/complete unit**



903TGA-24	Unifusor with thumbwheel valve & piston gauge	500mL	24/Box
904TGA-24	Unifusor with thumbwheel valve & piston gauge	1000mL	24/Box
905TGA-10	Unifusor with thumbwheel valve & piston gauge	3000mL	10/Box



903SGA-24	Unifusor with stopcock valve & piston gauge	500mL	24/Box
904SGA-24	Unifusor with stopcock valve & piston gauge	1000mL	24/Box
905SGA-10	Unifusor with stopcock valve & piston gauge	3000mL	10/Box

Indicates discontinued part, available while supplies last

**Unifusor II Disposable Pressure Infusor with Reusable Inflation System**

**Disposable Infusor - white Disposable bag with reusable inflation system**



2803SGA-24	Unifusor II with aneroid gauge & stopcock valve	500mL	24/Box
2804SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000mL	24/Box
2805SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000mL	10/Box



2903SGA-24	Unifusor II with piston gauge & stopcock valve	500mL	24/Box
2904SGA-24	Unifusor II with piston gauge & stopcock valve	1000mL	24/Box
2905SGA-10	Unifusor II with piston gauge & stopcock valve	3000mL	10/Box



**Disposable Replacement Parts**

2803X-24	Unifusor II Disposable white infusion bag	500mL	24/Box
----------	---	-------	--------

**Unifusor® Pressure Infusors**

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

**Unifusor Plus Reusable Pressure Infusor**

**Reusable Infusor - blue bag with white hose**



1103SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	500mL	24/Box
1104SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	1000mL	24/Box
1105SGA-10	Unifusor Plus with aneroid gauge & stopcock valve	3000mL	10/Box
1113SGA-24	Unifusor Plus with piston gauge & stopcock valve	500mL	24/Box

**Reusable replacement parts**



1103XB-05	Unifusor Plus blue reusable bag with white hose	500mL	5/Box
1104XB-05	Unifusor Plus blue reusable bag with white hose	1000mL	5/Box
1105XB-05	Unifusor Plus blue reusable bag with white hose	3000mL	5/Box

Indicates discontinued part, available while supplies last

**Unifusor Plus Reusable Pressure Infusor**



**Reusable Infusor - blue bag with black hose**

1103	Unifusor Plus with trigger / gauge assembly	500mL	1 Each
1104	Unifusor Plus with trigger / gauge assembly	1000mL	1 Each
1105	Unifusor Plus with trigger / gauge assembly	3000mL	1 Each

**Reusable replacement parts**



1103X	Unifusor Plus blue reusable bag with black hose	500mL	1 Each
1104X	Unifusor Plus blue reusable bag with black hose	1000mL	1 Each
1105X	Unifusor Plus blue reusable bag with black hose	3000mL	1 Each



1199	Trigger/ gauge assembly only		1 Each
------	------------------------------	--	--------

# SUPPLIES & ACCESSORIES

## ELECTROCARDIOGRAPHY


AAMI Color Coding



# ECG Combiner Lead Wires & Cables

# ECG COMBINER LEAD WIRE SYSTEM AAMI Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Disposable, DIN Shielded, Pinch, AAMI					
700-0006-50	3-Lead	Pinch	127cm /50 in	25 Sets	
* 700-0006-32	5 Lead	Pinch	127cm /50 in	25 Sets	
* 700-0006-33	5 Lead (V-Lead)	Pinch	127cm /50 in	25 Sets	



Disposable, DIN Shielded, Snap, AAMI					
700-0006-61	5 Lead	Snap	127cm /50 in	25 Sets	



Reusable, DIN Shielded, Pinch, AAMI					
* 700-0006-35	3 Lead	Pinch	74cm /29 in	1 Set	
* 700-0006-37	5 Lead	Pinch	74cm /29 in	1 Set	
* 700-0006-38	5 Lead (V-Lead)	Pinch	74cm /29 in	1 Set	
700-0010-00	3 Lead	Pinch	130cm/51 in	1 Set	
* 700-0010-01	5 Lead	Pinch	130cm/51 in	1 Set	
700-0010-02	5 Lead (V-Lead)	Pinch	130cm/51 in	1 Set	



Reusable, DIN Shielded, Snap, AAMI					
* 700-0007-35	3 Lead	Snap	74cm /29 in	1 Set	
* 700-0007-37	5 Lead	Snap	74cm /29 in	1 Set	
700-0010-19	5-Lead	Snap	130cm /51 in	1 Set	

\* Product not available for sale into CE Mark Countries (EU)

### ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter




Radiolucent, Reusable, DIN Unshielded, Pinch, AAMI					
700-0010-06	5 Lead, Radiolucent	Pinch	122cm/48 in	1 Set	
700-0010-08	5 Lead, Radiolucent	Pinch	183cm/72 in	1 Set	

# ECG Combiner Lead Wires & Cables

# ECG COMBINER LEAD WIRE SYSTEM AAMI Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Combiner Cables



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496					
* 700-0008-55	3 Lead	Spacelabs	305 cm/10 ft	1 Cable	
* 700-0008-57	5 Lead	Spacelabs	305 cm/10 ft	1 Cable	
700-0008-32	5 Lead	Spacelabs	457 cm/15 ft	1 Cable	
* 700-0008-56	10 Lead (5+5)	Spacelabs	305 cm/10 ft	1 Cable	

\* Product not available for sale into CE Mark Countries (EU)

### ECG Combiner Adapters



Spacelabs Telemetry, AAMI			
700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)	1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)	1 Adapter

### ECG Combiner Cables for Standardization



Other manufacturers					
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable	
700-0008-48	10 Lead (5+5)	Philips	305 cm/10 ft	1 Cable	
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable	

Indicates discontinued part, available while supplies last



# ECG LEAD WIRES & CABLES AAMI Color Coded - Shielded

## ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

#### Snap Sets

700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set



Pinch

#### Pinch Sets

700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set

### ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI



#### 24" Snap Individual Color

012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
<b>012-0498-34</b>	<b>Black</b>	<b>61 cm/24 in</b>	<b>5/Pack</b>
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

Indicates discontinued part, available while supplies last

# ECG LEAD WIRES & CABLES AAMI Color Coded - Shielded

## ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Cables, Shielded, AAMI



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496			
700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable


### ECG Accessories



Electrodes			
717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 pouches	Adult	50/Box

Indicates discontinued part, available while supplies last

### Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable			
700-0008-12	3 Lead	Neonatal	1 Cable
DIN Unshielded, Neonatal Electrodes.			
717-0080-00	Pre-wired neonatal electrodes, Pre-gelled, silver/silver chloride, single use	3-Lead	100 Sets/Box

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

## ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap

Snap Sets			
012-0283-11	3 Lead	61 cm/24 in	1 Set
015-0732-11	3 Lead	90 cm/35 in	1 Set
012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0605-10	5 Lead	64 cm/25.2 in	1 Set
015-0736-00	5 Lead	90 cm/35 in	1 Set



Pinch

Pinch Sets			
012-0284-11	3 Lead	61 cm/24 in	1 Set
015-0731-00	3 Lead	90 cm/35 in	1 Set
012-0286-11	5 Lead	61 cm/24 in	1 Set
015-0734-00	5 Lead	90 cm/35 in	1 Set

**Bold font indicates a new part number**

*Indicates discontinued part, available while supplies last*

### ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Snap

Snap Individual Color			
012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack



## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

## ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch

#### Pinch Individual Color

012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack

### Telemetry Transmitter Adapters



#### For use with ECG Combiner Lead Wires

700-0008-88	Aria Telemetry Combiner Adapter (96281)	1 Adapter
-------------	---	-----------



#### For use with DIN Unshielded Lead Wires

706-0191-00	Aria Telemetry Grouper (96281)	5/Bag
-------------	--------------------------------	-------

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

## ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Telemetry Accessories

#### Protective Solutions and Electrodes

010-2350-00	Aria Tele Bumper		5/Box
-------------	------------------	--	-------



015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
-------------	-----------------------------------	--	---------



717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches	Adult	50/Box
-------------	---	-------	--------

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

## ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### SpO<sub>2</sub> Cable & Sensors

Note: Only the below genuine Spacelabs and Nellcor branded SpO<sub>2</sub> sensors are compatible with Spacelabs Telemetry systems

#### SpO<sub>2</sub> Cable



700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

#### Reusable SpO<sub>2</sub> Sensors



015-0660-00	TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

#### Disposable SpO<sub>2</sub> Sensors



015-0662-00	TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO <sub>2</sub> vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last

# SUPPLIES & ACCESSORIES

## ELECTROCARDIOGRAPHY


IEC Color Coding



## ECG Combiner Lead Wires & Cables

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.

## ECG COMBINER LEAD WIRE SYSTEM IEC Color Coded



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Multi-pinch

#### Disposable, DIN Shielded, Multi-Pinch, IEC

700-0006-44	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-53	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets



#### Reusable, DIN Shielded, Multi-Pinch, IEC

700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set
-------------	-----------------	-------------	-------------	-------



Snap

#### Reusable, DIN Shielded, Snap, IEC

700-0007-41	3 Lead	Snap	74cm /29 in	1 Set
700-0007-44	5 Lead	Snap	74cm /29 in	1 Set

Indicates discontinued part, available while supplies last

### ECG Combiner Cables



#### Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)

700-0008-81	3 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------



700-0008-83	5 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------



# ECG LEAD WIRES & CABLES IEC Color Coded - Shielded

## ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Lead Wire Sets, Shielded, Reusable, IEC



Snap

#### Snap Sets

700-0007-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0007-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0007-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0007-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0007-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0007-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set



Pinch

#### Pinch Sets

700-0006-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0006-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0006-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0006-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0006-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0006-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set

### ECG Cables, Shielded, IEC



#### Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496

700-0008-09	3 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-07	5 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-01	10 Lead, shielded	305 cm/10 ft	1 Cable

Indicates discontinued part, available while supplies last

### ECG Accessories



#### Electrodes

717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches	Adult	50/Box
-------------	---	-------	--------

# ECG LEAD WIRES & CABLES IEC Color Coded - Shielded

## ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### Neonatal ECG Supplies & Accessories



#### DIN Unshielded, Neonatal ECG Cable

700-0008-13	3 Lead		1 Cable
-------------	--------	--	---------



#### DIN Unshielded, Neonatal Electrodes

717-0080-00	Pre-wired neonatal electrodes, pre-gelled, silver/silver chloride, CE	3-Lead	100 Sets/Box
-------------	---	--------	--------------

Indicates discontinued part, available while supplies last

## DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

## ECG LEAD WIRES & TELEMETRY IEC Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC



Snap

Snap Sets			
<b>015-0733-00</b>	<b>3 Lead</b>	<b>90 cm/35 in</b>	<b>1 Set</b>
<b>012-0328-10</b>	<b>5 Lead</b>	<b>61 cm/24 in</b>	<b>1 Set</b>
<b>015-0737-00</b>	<b>5 Lead</b>	<b>90 cm/35 in</b>	<b>1 Set</b>



Pinch

Pinch Sets			
<b>012-0327-10</b>	<b>5 Lead</b>	<b>61 cm/24 in</b>	<b>1 Set</b>
<b>015-0735-00</b>	<b>5 Lead</b>	<b>90 cm/35 in</b>	<b>1 Set</b>



For use with DIN Unshielded Lead Wires			
<b>706-0190-00</b>	<b>Aria Telemetry Grouper (96281), IEC</b>		<b>5/Bag</b>

### Telemetry Accessories



Disposable Pouches and Electrodes			
<b>010-2350-00</b>	<b>Aria Tele Bumper</b>		<b>5/box</b>
<b>015-0500-02</b>	<b>Ultraview Telemetry pouch (90347)</b>		<b>100/Box</b>
<b>717-0083-00</b>	<b>Electrode, Transparent Tape, 5 per pouch X 10 Pouches</b>	<b>Adult</b>	<b>50/Box</b>

**Bold font indicates a new part number**

**Indicates discontinued part, available while supplies last**



## DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

## ECG LEAD WIRES & TELEMETRY IEC Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### SpO<sub>2</sub> Cable & Sensors

#### SpO<sub>2</sub> Cable



700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

**Note: Only the below genuine Spacelabs and Nellcor branded SpO<sub>2</sub> sensors are compatible with Spacelabs Telemetry systems**

#### Reusable SpO<sub>2</sub> Sensors



015-0660-00	TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

#### Disposable SpO<sub>2</sub> Sensors



015-0662-00	TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO <sub>2</sub> vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last

# SUPPLIES & ACCESSORIES

## CAPNOGRAPHY & MULTIGAS



## NomLine Sampling Lines

NomLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



### NomLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

#### NomLine LH Nasal CO2 Cannula



704-3817-00	Adult	2 m/7 ft	25/Box
704-3819-00	Infant	2 m/7 ft	25/Box
704-3820-00	Adult w/O2	2 m/7 ft	25/Box
704-4624-00	Adult w/O2	4 m/13 ft	25/Box
704-3821-00	Pediatric w/O2	2 m/7 ft	25/Box



#### NomLine LH Nasal/Oral CO2 Cannula

704-3824-00	Adult w/O2	2 m/7 ft	25/Box
704-4628-00	Adult w/O2	4 m/13 ft	25/Box
704-3825-00	Pediatric w/O2	2 m/7 ft	25/Box



#### NomLine LH Single Nasal Prong CO2 Cannula

704-3826-00	Adult	2 m/7 ft	25/Box
-------------	-------	----------	--------

### NomLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



#### NomLine LH Airway Adapter Set

704-3814-00	Adult/Pediatric	2 m/7 ft	25/Box
704-4620-00	Adult/Pediatric	4 m/13 ft	25/Box

### NomLine LH (Low Humidity) Procedural Oxygen Masks

#### NomLine LH Nasal CO2 Cannula



704-5213-00 <sup>USA</sup>	Adult, Medium O2 Flow Rate	3 m/10 ft	30/Box
704-5214-00 <sup>USA</sup>	Pediatric, Medium O2 Flow Rate	3 m/10 ft	30/Box
704-5215-00 <sup>USA</sup>	Adult, High O2 Flow Rate	3 m/10 ft	30/Box
704-5216-00 <sup>USA</sup>	Pediatric, High O2 Flow Rate	3 m/10 ft	30/Box

*USA - Parts only available for sale in the USA*

**NomOLine Sampling Lines**

NomOLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomOLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



**NomOLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients**

**NomOLine HH Nasal CO2 Cannula**



704-3830-00	Adult	2 m/7 ft	25/Box
704-3832-00	Infant	2 m/7 ft	25/Box
704-3833-00	Adult w/O2	2 m/7 ft	25/Box
704-3834-00	Pediatric w/O2	2 m/7 ft	25/Box
704-4968-00 <sup>USA</sup>	Adult w/O2	4 m/13 ft	25/Box



**NomOLine HH Nasal/Oral CO2 Cannula**

704-3837-00	Adult w/O2	2 m/7 ft	25/Box
704-4933-00 <sup>USA</sup>	Infant w/O2	4 m/13 ft	25/Box
704-3838-00	Pediatric w/O2	2 m/7 ft	25/Box




**NomOLine HH Single Nasal Prong CO2 Cannula**

704-3839-00	Adult	2 m/7 ft	25/Box
-------------	-------	----------	--------

**NomOLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients**

**NomOLine HH Airway Adapter Set**



704-3827-00	Adult/Pediatric	2 m/7 ft	25/Box
704-3828-00	Adult/Pediatric	3 m/10 ft	25/Box
704-4989-00 <sup>USA</sup>	Adult/Pediatric	4 m/13 ft	25/Box
704-4367-00	Infant/Neonatal	2 m/7 ft	25/Box

*USA - Parts only available for sale in the USA*

## Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Sidestream Monitoring Supplies

**Please note:** the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.

#### CO2 Nasal Cannula with male luer connector



*	704-0176-00	Adult	2 m/7 ft	25/Box
CE	704-0176-01	Adult	2 m/7 ft	25/Box
**	704-0177-00	Pediatric	2 m/7 ft	25/Box
**	704-0178-00	Infant	2 m/7 ft	25/Box

#### CO2 Nasal Cannula w/O2 Delivery and male luer connector



*	704-0181-00	Adult	2 m/7 ft	25/Box
CE	704-0181-01	Adult	2 m/7 ft	25/Box
*	704-0180-00	Pediatric	2 m/7 ft	25/Box
**	704-0179-00	Infant	2 m/7 ft	25/Box
**	704-0186-00	Adult	3 m/10 ft	25/Box

#### CO2 Oral/Nasal Cannula w/O2 Delivery and male luer connector



*	704-0183-00	Adult	2 m/7 ft	25/Box
CE	704-0183-01	Adult	2 m/7 ft	25/Box
**	704-0187-00	Adult	3 m/10 ft	25/Box

Indicates discontinued part, available while supplies last

### Sidestream Monitoring Supplies - Intubated patients

#### Sampling Lines & Airway Adapters, single use



**	704-0191-00	Adult CO2 Sample Line w/ T-Adapter	2 m/7 ft	25/Box
**	704-0192-00	Adult CO2 Sample Line w/ T-Adapter	3 m/10 ft	25/Box



*	166-7121-00	Sidestream sampling line with male luer lock on both ends	3 m/10 ft	50/Box
CE	166-7121-02	Sidestream sampling line with male luer lock on both ends	3 m/10 ft	50/Box



*	704-0175-00	Sidestream airway adapter/T-Adapter.		25/Box
---	-------------	--------------------------------------	--	--------

\* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

\*\* Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients



**Sensor & Airway Adapters - required for mainstream monitoring**

010-1980-00	Mainstream sensor, reusable		1 Each
704-0173-00	Mainstream airway adapter, single use	Adult/Pediatric	25/Box
704-0174-00	Mainstream airway adapter, single use	Infant	10/Box







## Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)  
Dual Capnography Module

	Order Code	Description	Size	Qty
<b>Sidestream, Non-intubated, single patient use</b>				
	<b>CO<sub>2</sub> Nasal Cannula</b>			
	* 704-0010-00	Adult		10/Box
	CE 704-0010-10	Adult		10/Box
	* 704-0012-00	Pediatric		10/Box
	<b>CO<sub>2</sub> Nasal Cannula w/O<sub>2</sub> Delivery</b>			
	* 704-0014-00	Adult		10/Box
	CE 704-0014-10	Adult, (for sale into EU/UK only)		10/Box
	* 704-0016-00	Pediatric		10/Box
	CE 704-0016-10	Pediatric, (for sale into EU/UK only)		10/Box
	* 704-0015-00	Infant/Neonatal		10/Box
	<b>CO<sub>2</sub> Oral/Nasal Cannula</b>			
	* 704-0017-00	Adult		10/Box
	<b>CO<sub>2</sub> Oral/Nasal Cannula w/O<sub>2</sub> Delivery</b>			
	* 704-0020-00	Adult		10/Box
	CE 704-0020-10	Adult, (for sale into EU/UK only)		10/Box
	* 704-0022-00	Pediatric		10/Box
<b>Sidestream, Intubated, single patient use</b>				
	<b>Sidestream - short term up to 12 hours (patient non-humidified)</b>			
	* 704-0023-00	Airway adapter + sampling line, ET tube > 4.0 mm	Adult/Ped	10/Box
	CE 704-0023-10	Airway adapter + sampling line, ET tube > 4.0 mm	Adult/Ped	10/Box
	* 704-0027-00	CO <sub>2</sub> sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>		10/Box
	CE 704-0027-10	CO <sub>2</sub> sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>		10/Box

## Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



**(91517)**  
**Dual Capnography Module**

Order Code	Description	Size	Qty
<b>Sidestream - long term up to 120 hours (patient humidified)</b>			
* 704-0025-00	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	Adult/Ped	10/Box
CE 704-0025-10	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	Adult/Ped	10/Box
704-0026-00	Airway adapter + sampling line with Nafion	Ped/Infant	10/Box
* 704-0028-00	CO <sub>2</sub> sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>		10/Box
CE 704-0028-10	CO <sub>2</sub> sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>		10/Box

### Mainstream, for Intubated Patients Only

<b>Mainstream Airway Adapters, for use with CO<sub>2</sub> Sensor 704-0001-00</b>			
* 704-0005-00	CO <sub>2</sub> airway adapter, single use	Adult	24/Box
CE 704-0005-10	CO <sub>2</sub> airway adapter, single use	Adult	10/Box
* 704-0004-00	CO <sub>2</sub> airway adapter, single use	Neonate	24/Box
CE 704-0004-10	CO <sub>2</sub> airway adapter, single use	Neonate	10/Box
* 704-0002-00	CO <sub>2</sub> airway adapter, reusable	Adult	1 Each
CE 704-0002-10	CO <sub>2</sub> airway adapter, reusable	Adult	1 Each
CE 704-0003-10	CO <sub>2</sub> airway adapter, reusable	Neonate	1 Each

\* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

\*\* Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

# MULTIGAS ANALYZER

## Spacelabs Modules 90518 & 91518

### Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)




(91518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Multigas Supplies & Accessories

#### Anesthetic Agent Samples Lines



015-0128-00	SAMPLE LINE,MALE TO MALE,LUER,SPU,7'X.06",25/BX		25/Box
-------------	---	--	--------

015-0312-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518)	213 cm/7 ft.	50/Box
-------------	--	--------------	--------



015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518)	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518)	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



162-0047-01	Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518)	213 cm/7 ft.	50/Box
-------------	--	--------------	--------

Indicates discontinued part, available while supplies last

#### Water Traps, Reusable



717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)		10/Pack
-------------	--	--	---------

717-0448-00	D-FEND water trap cartridge, Green. Replace at least once every 24 hours or when occluded (91518)		10/Pack
-------------	--	--	---------

# SUPPLIES & ACCESSORIES











## PULSE OXIMETRY



**TruLink SpO2 Supplies & Accessories**  
**Spacelabs TruLink Technology**

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO2 supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can read from finger, ear, foot or toe.






Order Code	Description	Length/Size	Qty
<b>Spacelabs TruLink Technology</b>			
<b>Cables</b>			
 700-0030-00	TruLink SpO <sub>2</sub> cable	305 cm/10 ft.	1 Each
<b>Reusable sensors</b>			
 015-0660-00	TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
 015-0661-00	TruLink SpO <sub>2</sub> Y multi-site sensor <u>Includes:</u> (1) Sensor (1) ear clip (10) adult finger/neonate feet/hand wraps (1) pediatric finger or infant toe wraps	Universal Adult>30 kg Pediatric 10-50kg Infant 3-20kg Neonatal <3kg	1 Each
<b>Disposable Sensors</b>			
 015-0662-00	TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
 015-0663-00	TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
 015-0664-00	TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
 015-0665-00	TruLink SpO <sub>2</sub> vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
 015-0666-00	TruLink SpO <sub>2</sub> cloth sensor	Infant 3-20 kg/7-44 lbs.	24/Box
 015-0667-00	TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
<b>Dixtal Cables</b>			
 700-0029-00	TruLink SpO <sub>2</sub> cable compatible with Novametrix® sensors	30 cm/12 in.	1 Each

## Masimo® SpO2 Supplies & Accessories

### LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M).  
Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.



	Order Code	Description	Length/Size	Qty
<b>LNCS Masimo SET Technology</b>				
	<b>LNCS Adapter Cable</b>			
	700-0906-01	TruLink SpO <sub>2</sub> cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft	1 Each
	<b>LNCS Reusable Sensors</b>			
	690-0247-00	Masimo LNCS DC-I, SpO <sub>2</sub> digit sensor	Adult > 30 kg	1 Each
	690-0248-00	Masimo LNCS DC-IP, SpO <sub>2</sub> digit sensor	Pediatric 10-50 kg	1 Each
	690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO <sub>2</sub> ear sensor	Adult > 30 kg	1 Each
	690-0291-00	Masimo LNCS Y-I multisite sensor	> 1 kg	1 Each
<b>LNCS Masimo SET Technology</b>				
	<b>LNCS Disposable Sensors</b>			
	690-0251-00	Masimo LNCS Amtx, adhesive SpO <sub>2</sub> sensor, 18 in.	Adult > 30 kg	20/Box
	690-0252-00	Masimo LNCS Pmtx, adhesive SpO <sub>2</sub> sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
	690-0258-00	Masimo LNCS Inf-L, adhesive SpO <sub>2</sub> sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
	690-0259-00	Masimo LNCS Neo-L, adhesive SpO <sub>2</sub> sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box



## Masimo® SpO<sub>2</sub> Supplies & Accessories RD SET Technology

For use with Ultraview SL™ Command Module (91496-M). Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.



Order Code	Description	Length/Size	Qty
<b>RD SET Masimo Technology</b>			
<b>RD SET Adapter Cable</b>			
010-1649-00	Masimo RD SET SpO <sub>2</sub> Cable	12 ft.	1 Each
<b>RD SET Reusable Sensors</b>			
690-0267-00	Masimo RD SET DCI reusable SpO <sub>2</sub> sensor	Adult > 30 kg	1 Each
690-0270-00	Masimo RD SET DCI-P reusable SpO <sub>2</sub> sensor	Pediatric 10-50 kg	1 Each
690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO <sub>2</sub> sensor	Adult > 30 kg	1 Each
690-0281-00	Masimo RD SET Y-I multi-site reusable SpO <sub>2</sub> sensor	> 1 kg	1 Each
<b>RD SET Masimo Technology</b>			
<b>RD SET Adhesive Disposable Sensors</b>			
690-0283-00	Masimo RD SET Adt, adhesive SpO <sub>2</sub> sensor	Adult, > 30 kg finger or toe	20/Box
690-0285-00	Masimo RD SET Inf, adhesive SpO <sub>2</sub> sensor	Infant, 3—20 kg, thumb or great toe	20/Box
690-0286-00	Masimo RD SET Neo, adhesive SpO <sub>2</sub> sensor	Neonatal <3 kg, hand or foot Adult > 40 kg finger or toe	20/Box

**Nellcor® SpO<sub>2</sub> Supplies & Accessories**  
**Nellcor™ OxiMax® Technology**

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor. Cables and connectors are color-coded purple. **Note, use only genuine Nellcor sensors.**



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

**Nellcor OxiMax Technology**

**Cable**



700-0792-00	TruLink SpO <sub>2</sub> cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
-------------	---	--	--------

**Reusable Sensors**

**Note:** Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien), now part of Medtronic.



690-0003-01	Nellcor DS-100A, finger SpO <sub>2</sub> sensor	Adult	1 Each
-------------	---	-------	--------



690-0004-01	Nellcor OXI-A/N, SpO <sub>2</sub> sensor	Adult/Neonatal	1 Each
-------------	--	----------------	--------



690-0039-01	Nellcor OXI-P/I, SpO <sub>2</sub> sensor	Pediatric/Infant	1 Each
-------------	--	------------------	--------



690-0215-01	Nellcor DS-Y, multi-site SpO <sub>2</sub> sensor	Universal	1 Each
-------------	--	-----------	--------

# SUPPLIES & ACCESSORIES








## TEMPERATURE



## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.











Order Code	Description	Length/Size	Qty
<b>Monitor Adapter Cables</b>			
<b>Adapters</b>			
 700-0031-00	TruLink Dual Temperature Cable. *Required for all reusable probes and as directed for single-patient-use disposable system	30 cm/12 in.	1 Each
<b>Reusable Probes/Sensors</b>			
<b>Skin Surface</b>			
 690-0295-00	YSI 400 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0296-00	YSI 700 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0297-00	YSI 700 series compatible skin surface sensor temperature probe, pediatric size	305 cm/10 ft	1 Each
<b>Esophageal / Rectal</b>			
 690-0294-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
 690-0298-00	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
 690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each

Indicates discontinued part, available while supplies last

## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



	Order Code	Description	Length/Size	Qty
<b>Single-Patient use temperature system</b>				
<b>Extension cables for single use probes</b>				
	700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors, adapter required.	305 cm/10 ft	1 Each
	700-0401-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors. Direct connection to Command module, no adapter required.	305 cm/10 ft	1 Each
<b>Skin</b>				
	701-4000-00	Skin Temperature Sensor		20/Box
<b>Esophageal/ Stethoscope</b>				
	701-4009-00	9 French, YSI 400 series compatible		20/Box
	701-4012-00	12 French, YSI 400 series compatible		20/Box
	701-4018-00	18 French, YSI 400 series compatible		20/Box
<b>Esophageal/ Rectal</b>				
	701-4109-00	9 French, YSI 400 series compatible		20/Box
	701-4112-00	12 French, YSI 400 series compatible		20/Box
<b>Foley catheter with temperature probe</b>				
	701-4214-00	14 French, YSI 400 series compatible		20/Box
	701-4216-00	16 French, YSI 400 series compatible		20/Box
<b>Tympanic temperature sensor</b>				
	701-4300-00	Adult, YSI 400 series compatible		20/Box
<b>Exergen</b>				
<i>For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) &amp; DM4 (91331) Monitors</i>				
	010-2157-00	Temporal Artery Thermometer		1 Each

# SUPPLIES & ACCESSORIES

## INVASIVE BLOOD PRESSURE





## Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Integrated Adapter Cables (connects directly to transducer)

#### Edward Lifesciences



**700-0293-01** Single pressure integrated cable

396 cm/13 ft

1 Each

**015-0727-00** Single pressure integrated cable

400 cm/13 ft

1 Each



**700-0292-01** Dual pressure integrated cable

396 cm/13 ft

1 Each

**015-0728-00** Dual pressure integrated cable

400 cm/13 ft

1 Each



#### Abbott Laboratories/Hospira/ICU Medical



**700-0295-01** Single pressure integrated cable

396 cm/13 ft

1 Each

**015-0729-00** Single pressure integrated cable

400 cm/13 ft

1 Each



**700-0294-01** Dual pressure integrated cable


396 cm/13 ft

1 Each

**015-0730-00** Dual pressure integrated cable

400 cm/13 ft

1 Each



**Bold font indicates a new part number**

**Indicates discontinued part, available while supplies last**

### Interface Cable & Adapters



#### Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)

690-0021-00 Requires adapter cable 700-0091-00 or 700-0028-00

366 cm/12 ft

1 Each



#### Interface Adapter- Transducer interface cable required from transducer company

700-0791-00 Single pressure cable

122 cm/4 ft

1 Each

700-0028-00 Dual pressure cable

122 cm/4 ft

1 Each

### Transducer Kits & Accessories

#### Pressure Transducer Kits, Single use



690-0008-00 Transducer Kit, Patient Mount, DTX/Plus

152.5 cm/60 in

5/Box



690-0009-00 Transducer Kit, Pole Mount, DTX/Plus

152.5 cm/60 in

5/Box

# SUPPLIES & ACCESSORIES

## CARDIAC OUTPUT









# CARDIAC OUTPUT


## Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



	Order Code	Description	Length/Size	Qty
<b>Cardiac Output Cables &amp; Probes</b>				
	<b>Cables</b>			
	700-0027-00	TruLink Cardiac Output Adapter Cable	30 cm/12 in.	1 Each
	306655-003	Cardiac Output cable compatible with Baxter Edwards/Abbott injectate temperature probes and most models of thermodilution catheters		1 Each
	306655-004	Cardiac Output cable compatible with Merit Medical, Becton Dickinson/Datex-Ohmeda and Argon Medical injectate temperature probes and most models of thermodilution catheters		1 Each
	<b>Injectate Probes</b>			
	690-0023-00	In-line injectate temperature probe Compatible with Datex-Ohmeda SP4045		1 Each
	690-0024-00	Reference solution in-line injectate temperature probe. Compatible with Datex-Ohmeda SP5030		1 Each
	690-0033-00	In-line injectate temperature probe. Compatible with Datex-Ohmeda SP4042. Not compatible with Abbott or Baxter in-line injectate probes or delivery systems.		1 Each

## Cardiac Output Injectate Kits

<b>Closed Loop Injectate Kits</b>				
	690-0293-00		SP5045 Critikit In-Line Sensor, Control Syringe closed loop injectate kit, sterile	5/Box

# SUPPLIES & ACCESSORIES

## RECORDING & PRINTING



**Paper**

High quality Z-fold paper for use with the Bedside Printer Module (90449) and Central Printer (90469).

Roll paper for use with the Standalone printer (91449) and Spacelabs Portable Monitors.

Paper for Heartwave and Cardioline devices.



**(91449)  
Printer**

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

**Thermal Paper**



**Z-fold Paper for use with 90449 and 90469**

006-0196-00	TruLink Thermal Paper, Z-Fold Pads, 2-channel printing, blank paper, 230 Sheets/pad	27.4 m/90 ft	10 Pads/Box
-------------	---	--------------	-------------



**Roll Paper for use with 91449, 91367, 91369, 91370 & 91390**

307438-001	TruLink Thermal Paper, Rolls, single- or dual-channel printing, 50 mm/1.97 inch wide blank paper	30.1 m/100 ft	10 Rolls/Box
------------	--	---------------	--------------

**Recording Paper**

20257-004	Thermal recorder paper, black	HearTwave II	each
010-2246-00	PAPR,THERM,Z-FOLD, 210X295MM X 180 SHTS/10RM,IEC	Cardioline	each

## SUPPLIES & ACCESSORIES BISPECTRAL INDEX





Bispectral Index Supplies

For use with Spacelabs BISx™ Module (91482).



Order Code	Description	Qty
<b>BISx Cables &amp; Sensors</b>		
<b>BISx Cables</b>		
719-0004-01	BISx Patient Interface cable	1 Each
719-0006-20	BISx Pod	1 Each
719-0007-01	BISx host monitor cable	1 Each

**BISx Sensors**

**Note:** Sales of BISx sensors are limited to one box with the purchase of capital equipment. All future purchases should be made directly through Covidien, now part of Medtronic.

719-0003-10	BISx Extend Sensor, single patient use	25/Box
719-0002-10	BISx Quatro Sensor, single patient use	25/Box

# SUPPLIES & ACCESSORIES

## ULTRAVIEW DM3 DUAL MODE VITAL SIGNS MONITOR



# VITAL SIGNS MONITOR Ultraview DM3

## Vital Sign Monitoring

Ultraview DM3 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

See previous sections for NIBP and SpO2 supplies:

- The DM3 uses the TruLink NIBP cuffs (single tube, bayonet fitting) and the standard air hose 714-0018-00.
- The DM3 used the same cables and sensors, by brand, for pulse oximetry.



DM3 91330

Order Code	Description	Size	Length/Size	Qty
<b>Supplies unique to the DM3 monitor:</b>				
019-0539-00	Replacement Temperature			1 Each
Indicates discontinued part, available while supplies last				

# SUPPLIES & ACCESSORIES

## DM4 VITAL SIGNS MONITOR



# VITAL SIGNS MONITOR

## Ultraview DM4









### Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO<sub>2</sub>, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

The DM4 uses TruLink branded NIBP cuffs (single tube, bayonet fitting) and air hose (714-0018-00).



91331

Order Code	Description	Length/Size	Qty
<b>Supplies unique to the DM4 monitor</b>			
<b>LNCS Adapter Cable</b>			
 010-2227-00	DM4 Masimo LNCS SpO <sub>2</sub> cable compatible with Masimo LNCS sensors, see section on Masimo LNCS sensors	305 cm/10 ft	1 Each
<b>RD SET Adapter Cable</b>			
 010-2193-00	DM4 Masimo RD SET SpO <sub>2</sub> cable compatible with Masimo RD SET sensors, see section on RD SET sensors	305 cm/10 ft	1 Each
<b>Nellcor OxiMax Cable</b>			
 010-2192-00	DM4 Nellcor OxiMax SpO <sub>2</sub> cable compatible with Nellcor OxiMax sensors, see section on Nellcor sensors	305 cm/10 ft	1 Each
<b>Filac</b>			
 010-2174-00	Filac dock and insulation chamber		1 Each
 010-2175-00	Filac temperature probe		1 Each
 010-2228-00	Filac temperature oral probe covers, 20/bx		20/box 25/boxes
<b>Exergen</b>			
 010-2157-00	Temporal Artery Thermometer		1 Each
<b>NeoNatal NIBP supplies</b>			
 <b>714-0087-00</b>	<b>NIBP,HOSE ASSY,NEONATE,DM4</b> Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	<b>2.5 meters</b>	<b>1 Each</b>

**Bold font indicates a new part number**

# SUPPLIES & ACCESSORIES

## ELANCE





## Elace Vital Signs Monitor







The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add Microstream end tidal CO2.

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below



Elace Vital Signs Monitor  
93300

Order Code	Description		Qty	
<b>Supplies unique to the Elance monitor</b>				
<b>ECG Lead Wires &amp; Cables (can also use standard TruLink ECG cables and leadwire systems)</b>				
<b>ECG Integrated Cables and Lead Wires</b>				
	700-0008-22	5-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
	700-0008-23	5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
	700-0008-24	3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
<b>Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)</b>				
	<b>SpO2 Cable</b>			
700-0287-00	SpO2 adapter cable	305 cm/10 ft	1 Each	
<b>LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)</b>				
	<b>LNCS Adapter Cable</b>			
690-0256-00	SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft	1 Each	
<b>Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)</b>				
	<b>Cable</b>			
700-0290-00	SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft	1 Each	
<b>Temperature</b>				
	<b>Temperature Probe</b>			
690-0261-01	Skin/Surface Temp probe, reusable, YSI-400	305cm/10ft	1 Each	
<b>Invasive Blood Pressure (See IBP section for compatible sensors and transducers)</b>				
	<b>Interface Cable</b>			
690-0021-00	Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	366 cm/12 ft	1 Each	
	Indicates discontinued part, available while supplies last			

# SUPPLIES & ACCESSORIES

## C50 VITAL SIGNS MONITOR













## Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
<b>ECG leadwire systems</b>				
 040-000908-00	CABLE,INTEGRATED 3-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000909-00	CABLE,INTEGRATED 5-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000910-00	CABLE,INTEGRATED 10-WIRE CLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000911-00	CABLE,INTEGRATED 3-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000912-00	CABLE,INTEGRATED 5-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000913-00	CABLE,INTEGRATED 10-WIRE CLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
<b>Pulse Oximetry Supplies</b>				
 040-000869-00	CABLE,INTEGRATED ADT FINGER CLIP SENSOR,C50	3 meters	Reusable	1 Each
 040-000769-00	CABLE,SPO2,12PIN TO DB9F,C50 <i>Sensors sold seperately</i>	3 meters	Reusable	1 Each
 040-000726-00	SENSOR,SPO2,ADT,FINGER CLIP,C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each
 040-000730-00	SENSOR,SPO2,ADT/NEO/PED,WRAP, 3',C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each









## Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
<b>Non-Invasive Blood Pressure supplies</b>				
 040-000808-00	ADULT NIBP TUBE (WHITE)2M	2 meters	Reusable	1 Each
	Reusable or Disposable Cuffs See NIBP section under TruLink Reusable or Disposable cuffs. These are compatible with the C50. <i>Requires use of the 040-00808-00 air hose (above).</i>			
 <b>714-0086-00</b>	<b>NIBP,HOSE ASSY,NEONATE,C50</b> Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	<b>2.5 meters</b>	<b>Reusable</b>	<b>1 Each</b>
	<b>Bold font indicates a new part number</b>			
<b>Capnography Supplies: requires external Phasein ISA module, 099-000007-00</b>				
 040-000400-00	<b>CABLE,CO2/AG,C50</b> <b>Connects the Phasein ISA module to the C50</b>	<b>1 meter</b>	<b>Reusable</b>	<b>1 Each</b>
	For airway adapter, cannula and sampling line sets, please see the capnography supply section under Nomoline SPU supplies. These are compatible with the C50.			
<b>Temperature sensors</b>				
 040-000246-00	<b>PROBE,TEMP,ADT,SKIN,10K,C50</b>	<b>3 meters</b>	<b>Reusable</b>	<b>1 Each</b>
 040-000387-00	PROBE,TEMP,NEO SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000385-00	PROBE,TEMP,ADT,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each
 040-000386-00	<b>PROBE,TEMP,NEO,RECTAL/ESOPHAGEAL,10K,C50</b>	<b>3 meters</b>	<b>Reusable</b>	<b>1 Each</b>









## Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
<b>Invasive Blood Pressure Cables and Transducer kits</b>				
 040-000454-00	INTERFACE CABLE,IBP,EDWARD CONNCTR,C50 <i>Requires 040-000211-00 transducers</i>	 3 meters	Reusable	1 Each
 040-000360-00	INTERFACE CABLE,IBP,ABBOTT CONNCTR,C50 <i>Requires 040-000870-00 transducers</i>	 3 meters	Reusable	1 Each
 040-000870-00	DISPOSABLE PRESSURE TRANSDUCER KIT(COPPER),QTY2 <i>Abbott compatible, requires 040-000360-00 cable</i>	sterile kit	Disposable	2 Each
 040-000456-00	INTERFACE CABLE,IBP,BD CONNCTR,C50 <i>Requires 040-000403-00 transducers</i>	 3 meters	Reusable	1 Each
<b>Recording Paper</b>				
 040-000116-00	RECORDER PAPER,THERMAL,C50	50mm X 20 meters		1 roll Each

# SUPPLIES & ACCESSORIES

## ANESTHESIA





## Anesthesia

Anesthesia System For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseSirius Anesthesia Machines.



BleasFocus



BleasSirius



ARKON

Order Code	Description	Qty
------------	-------------	-----

### Anesthesia Supplies

#### Patient Flow Sensors

*For use with BleasSirius and BleasFocus Anesthesia Machines only*



* 10110090	Adult Patient flow sensor	1 Each
717-0428-00	D-Lite reusable air flow sensor, may require luer caps	1 Each
050-0842-00	Luer caps for D-Lite air flow sensor	Pack of 5
<i>See image for assembly, available for sale to EU</i>		
10110218	2M Patient flow sensor tubing assembly	1 Each
10110102	3M Patient flow sensor tubing assembly	1 Each
<b>10110103</b>	<b>6M Patient flow sensor tubing assembly</b>	<b>1 Each</b>
70300001	O2 probe with tee	1 Each
376-0561-01	Pediatric flow sensor kit *EU/UK ONLY*	1 Each
FR101029	Pediatric flow sensor and tubing kit, 2.0 meters	1 Each
FR101030	Pediatric flow sensor and tubing kit, 3.0 meters	1 Each

**Indicates discontinued part, available while supplies last**

\* Product not available for sale into CE Mark Countries (EU)

# SUPPLIES & ACCESSORIES

## AMBULATORY BLOOD PRESSURE MONITOR



## Ambulatory Blood Pressure

Ambulatory blood pressure is a leading diagnostic tool. Rather than relying on one isolated reading, it compiles a series of measurements resulting in increased accuracy. Spacelabs offers a line of reusable and disposable cuffs, covers and pouches for use with the Spacelabs 90227, 90207 and 90217A APB Monitors.



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

### Reusable ABP Cuffs




#### UltraCheck Reusable Nylon ABP Cuffs

016-1320-02	Statcorp reusable ABP cuff	child	13 - 18 cm	1 Each
016-1826-02	Statcorp reusable ABP cuff	small adult	18 - 26 cm	1 Each
016-2635-02	Statcorp reusable ABP cuff	adult	26 - 35 cm	1 Each
016-3242-02	Statcorp reusable ABP cuff	large adult	32 - 42 cm	1 Each



#### Spacelabs Reusable Cloth ABP Cuffs (includes bladder, wrap, hose, fittings)

015-0118-05Q	Spacelabs reusable ABP cuff kit	child	12 - 20 cm	1 Each
015-0067-05Q	Spacelabs reusable ABP cuff kit	small adult	17 - 26 cm	1 Each
015-0068-06Q	Spacelabs reusable ABP cuff kit	adult	24 - 32 cm	1 Each
016-0077-06Q	Spacelabs reusable ABP cuff kit	large adult	32 - 42 cm	1 Each
016-0109-05Q	Spacelabs reusable ABP cuff kit	XL adult	38 - 50 cm	1 Each



### Disposable Covers



#### Disposable Cuff Covers (for use with Spacelabs reusable ABP cuffs above)

706-0182-00	Universal cuff cover	Standard	50/Box
706-0183-00	Universal cuff cover	Large	50/Box

### Disposable ABP Cuffs



#### Spacelabs Soft Fabric Disposable ABP Cuffs (no fitting)

015-1826-00	Spacelabs disposable ABP cuff	small adult	18 - 26 cm	20/Box
015-2635-00	Spacelabs disposable ABP cuff	adult	26 - 35 cm	20/Box
015-3242-00	Spacelabs disposable ABP cuff	large adult	32 - 42 cm	20/Box
015-4250-00	Spacelabs disposable ABP cuff	XL adult	42 - 50 cm	20/Box
ZFT-QC8220	Quick release metal fittings, female bayonet (HP), 1/8 in. barb, For use with above single tube, disposable cuffs with no fitting.			each

## Ambulatory Blood Pressure

Ambulatory blood pressure is a leading diagnostic tool. Rather than relying on one isolated reading, it compiles a series of measurements resulting in increased accuracy. Spacelabs offers a line of reusable and disposable cuffs, covers and pouches for use with the Spacelabs 90227, 90207 and 90217A APB Monitors.



	Order Code	Description	Size	Range	Qty
<b>ABP Accessories</b> 	<b>Pouches - For use with ABP Monitor 90207</b>				
	016-0410-01	Cloth carrying pouch with Velcro closure, black			1 Each
	<b>Pouches - For use with ABP Monitor 90217</b>				
	016-0340-01	Cloth carrying pouch with Velcro closure, black			1 Each
	<b>Pouches - For use with ABP Monitor OnTrak 90227</b>				
	706-0200-00	Disposable pouch			50/Pack
	015-0715-01	Reusable pouch			1 Each
	015-0716-00	Reusable shoulder strap (for use with OnTrak ABP Recorder pouch Item No. 015-0715-00)			1 Each
	<b>Support</b>				
	015-0070-00	Cuff support harness			1 Each
016-0080-00	Nylon web adjustable waist belt, d-ring adjustable fastener			1 Each	
PK0001	Nylon webbing adjustable 66 inch x 1 inch with slide release buckle. For use as waist belt or shoulder strap.			1 Each	
<b>Miscellaneous ABP Accessories</b>					
712-0773-00	Adapter, converts ABP Luer Monitors (90207 and 90217) to quick disconnect			1 Each	
712-0794-00	Quick Disconnect metal hose connector for ABP Monitors 90207 & 90217			1 Each	
712-0060-00	QUICK CONNECTOR,1/8 TUBE,NIBP,CUFF			1 Each	

Indicates discontinued part, available while supplies last

# SUPPLIES & ACCESSORIES

## HOLTER AND EVENT MONITORING





# Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.











Lifecard CF

evo

Aria

CardioCall

Order Code	Description	For Use With	Qty
<b>Aria/evo Hook-Up Kits</b>			
 23189-106-401	4 electrode hook-up kit <b>Includes:</b> 4 clear tape electrodes, 1 AA battery, 4 prep pads, 1 prep scrub, 1 Aria disposable pouch, 6 tape strips, 1 prep razor and 1 patient diary	Aria/evo	20/Case
 50-0049	4 electrode hook-up kit with disposable evo chest pouch <b>Includes:</b> 4 clear tape electrodes, 4 prep pads, 1 prep scrub, 1 evo disposable pouch, 6 tape strips, 1 prep razor, 1 patient diary and pencil	Aria/evo	Each
23-0255	Lanyard clip & assembly	CardioCall	Each
<b>Lead Wires</b>			
 25452-525	4-Lead wire set, snap, 3 - 14" leads (orange/brown/red) and 1 - 7" lead (white)	Aria/evo	4LW/set
 25452-535	4-Lead wire set, snap, 14" leads (red/white/orange/brown)	Aria/evo	4LW/set
 25452-529	4-Lead wire set, snap, 40" leads (red/white/orange/brown)	Aria/evo	4LW/set
 25452-519	7-Lead wire set, snap, 18" leads (white/black/red/green/orange/blue/brown)	423, 423A, 459, 461, 461A, 464, 483, 485	7 LW/set
700-0366-00	White, snap, 14" lead wires	Aria/evo	10 /pack
19-7509	6-Lead wire set, snap, 24 in, AHA (requires 700-0316-00)	Lifecard CF	6LW/Set
23-0059	Patient cable, 1 channel, 2 electrodes (black/white), US	CardioCall	Each
23-0051	Patient cable, 2 channel, 3 electrodes (red/yellow/green)	CardioCall	Each
<b>Electrodes and Prep supplies</b>			
 392196-001	Adhesive tape remover pads, 100/box		100/Box
 717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches		50/Box
717-0071-00	ECG electrode, Ambu Blue Sensor L		5/Pack



# Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF



evo



Aria



CardioCall

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

### Lifecard CF Patient Cables



700-0316-00	3-channel, 6-lead, 22 in, AHA (requires 19-7509)	Lifecard CF	Each
700-0317-00	3-channel, 3-lead, snap, long	Lifecard CF	Each
700-0318-00	3-channel, 3-lead, snap, extra-long	Lifecard CF	Each
700-0319-00	3-channel, 3-lead, snap, short	Lifecard CF	Each
700-0320-00	CABLE,PX,LCF,2SNAP,CLIP,RD/WH.ROHS	Lifecard CF	Each
700-0321-00	CABLE,PX,LCF,2SNAP,SHORT,WO/CLIP,RD/WH.ROHS	Lifecard CF	Each
<b>700-0322-00</b>	<b>CABLE,PX,LCF,7RECEPTACLE,AHA,WO/LEADWIRES,ROHS</b>	<b>Lifecard CF</b>	<b>Each</b>
700-0323-00	2-channel, 4-lead, snap, long	Lifecard CF	Each
700-0324-00	2-channel, 4-lead, snap, short	Lifecard CF	Each

### Lifecard 12 Patient Cables



<b>*700-0325-00</b>	<b>Varios Yoke</b>	<b>Lifecard 12</b>	<b>Each</b>
<b>*700-0326-00</b>	<b>cable with varios yoke, 12-lead patient cable, EU</b>	<b>Lifecard 12</b>	<b>Each</b>
<b>*700-0327-00</b>	<b>cable with varios yoke, 12-lead patient cable, US</b>	<b>Lifecard 12</b>	<b>Each</b>
700-0328-00	12-lead patient cable, EU	Lifecard 12	Each
700-0329-00	12-lead patient cable, US	Lifecard 12	Each
700-0330-00	12-lead patient cable, long, EU	Lifecard 12	Each
700-0331-00	12-lead patient cable, long, US	Lifecard 12	Each

**\*Available for sale until April 30th, 2024**

### Pouches, Shoulder Straps & Belts



27249-7-401	Disposable pouch	Aria	50/Pack
50-0191	HOOK UP KIT WITH SINGLEUSE POUCH, CASE OF 20	evo	20/Case
50-0042	Belt holster	evo	Each
PK0001	Fabric waist belt	evo	Each
706-0189-00	Disposable pouch	evo	100/Pack
46-0005	Disposable cord (use with 706-0189-00)	evo	50m roll
50-0050	Disposable pouch	evo	50/Pack
46-0408	Fabric pouch set, 2 pouches with laces, Cordlock	Lifecard CF	Each
19-7130	Laces for pouch set	Lifecard CF	Each
19-7131	Cordlock for pouch set	Lifecard CF	Each
46-0092	Fabric pouch	Lifecard CF	Each
46-0512	Disposable pouch	Lifecard CF	100/Pack

**Indicates discontinued part, available while supplies last**

# Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF



evo



Aria



CardioCall

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

### Recorder Memory and Download

CFC90MB	Compact flash card memory, 90MByte	Lifecard CF	Each
CFC256MB	Compact flash card memory, 256MB	LC12	Each
46-0063	Flash card plastic protective case	CFCxxMB	Each
P60-7056	Compact USB Flashcard Reader	CFCxxMB	Each

### Patient Diaries

073-0267-00	7 Patient Diary, English CE0123	Aria, evo, Lifecard CF	Each
-------------	---------------------------------	---------------------------	------



# SUPPLIES & ACCESSORIES ECLIPSE MINI AND ECLIPSE PRO



Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

*Note: Patches may not be available in all regions at this time.*

## Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code	Description	Length/Style	Qty																																														
<b>Eclipse Pro</b>																																																	
<b>Leadwire Sets</b>																																																	
700-0378-00	3 Lead snap on with yoke <b>Short</b>	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>M</td> <td>Yellow</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C5R</td> <td>Green</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	205	8	510	20	865	34	M	Yellow	255	10	560	22	915	36	C5R	Green	205	8	510	20	865	34	Each								
Snap	Color				Short		Medium		Long																																								
			mm	Inch	mm	Inch	mm	Inch																																									
C5	Red	205	8	510	20	865	34																																										
M	Yellow	255	10	560	22	915	36																																										
C5R	Green	205	8	510	20	865	34																																										
700-0391-00	3 Lead snap on with yoke <b>Medium</b>		Each																																														
700-0390-00	3 Lead snap on with yoke <b>Long</b>		Each																																														
700-0379-00	4 Lead snap on with yoke <b>Short</b>	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>C3</td> <td>Orange</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>M</td> <td>White</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C1</td> <td>Brown</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	230	9	510	20	865	34	C3	Orange	255	10	560	22	915	36	M	White	255	10	560	22	915	36	C1	Brown	230	9	510	20	865	34	Each
Snap	Color				Short		Medium		Long																																								
			mm	Inch	mm	Inch	mm	Inch																																									
C5	Red	230	9	510	20	865	34																																										
C3	Orange	255	10	560	22	915	36																																										
M	White	255	10	560	22	915	36																																										
C1	Brown	230	9	510	20	865	34																																										
700-0394-00	4 Lead snap on with yoke <b>Medium</b>		Each																																														
700-0393-00	4 Lead snap on with yoke <b>Long</b>		Each																																														
700-0396-00	10 Lead snap on with yoke <i>IEC color coding generally used in European Union</i>	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th colspan="2">Length</th> </tr> <tr> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>L</td> <td>610</td> <td>24</td> </tr> <tr> <td>F</td> <td>585</td> <td>23</td> </tr> <tr> <td>C6</td> <td>635</td> <td>25</td> </tr> <tr> <td>C5</td> <td>535</td> <td>21</td> </tr> </tbody> </table>	Snap	Length		mm	Inch	L	610	24	F	585	23	C6	635	25	C5	535	21	Each																													
Snap	Length																																																
	mm	Inch																																															
L	610	24																																															
F	585	23																																															
C6	635	25																																															
C5	535	21																																															
700-0377-00	10 Lead snap on with yoke <i>AHA color coding generally used in USA, outside of EU</i>	<table border="1"> <tbody> <tr> <td>C4</td> <td>460</td> <td>18</td> </tr> <tr> <td>C3</td> <td>355</td> <td>14</td> </tr> <tr> <td>C2</td> <td>330</td> <td>13</td> </tr> <tr> <td>C1</td> <td>330</td> <td>13</td> </tr> <tr> <td>N</td> <td>380</td> <td>15</td> </tr> <tr> <td>R</td> <td>610</td> <td>24</td> </tr> </tbody> </table>	C4	460	18	C3	355	14	C2	330	13	C1	330	13	N	380	15	R	610	24	Each																												
C4	460	18																																															
C3	355	14																																															
C2	330	13																																															
C1	330	13																																															
N	380	15																																															
R	610	24																																															
712-1052-01	ECLIPSE PRO ECG SENSORS <i>Not available for sale in USA</i>		Box of 10																																														
<b>Accessories</b>																																																	
706-0207-00	Eclipse Pro belt clip		Each																																														
073-0267-00	7 Patient Diary, English CE0123		Each																																														
392196-001	Adhesive tape remover pads, 100/box		100/Box																																														
717-0071-00	ECG electrode, Ambu Blue Sensor L		5/Pack																																														

Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

*Note: Patches may not be available in all regions at this time.*


## Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

	Order Code	Description	Length/Style	Qty
<b>Eclipse Mini</b>				
	<b>Eclipse Mini Supplies</b>			
	712-1053-01	ECLIPSE MINI ECG SENSORS <i>Not available for sale in USA</i>		Box of 10
	712-1053-02	ECLIPSE MINI ECG SENSORS, DEFIB SAFE <i>Available for USA September 2023</i>		Box of 10

# SUPPLIES & ACCESSORIES

## ECG STRESS TESTING





## Brand Name Stress Test Devices

As part of the Spacelabs Sentinel system Spacelabs offers supplies for the Suntech Tango M2 ABP system. Below you can find the related supplies and accessories for these devices.



Suntech Tango M2

Order Code	Description	Dimensions	Qty
<b>Suntech Tango NIBP Supplies for Stress Test</b>			
010-2238-00	CUFF, ORBIT-K, SMALL ADULT		Each
98-0061-02	Orbit-K cuff, adult, 23-35 cm		Each
98-0061-05	Orbit-K cuff, adult plus, 27-40 cm		Each
98-0061-03	Orbit-K cuff, large adult, 32-44 cm		Each
010-2237-00	K-SOUND MICROPHONE,14" FOR SUNTECH TANGO		Each
010-2241-00	K-SOUND MICROPHONE,12" FOR SUNTECH TANGO		Each
98-0002-50	Microphone pads		5/Pack



Indicates discontinued part, available while supplies last

# SUPPLIES & ACCESSORIES

## CardioPulse ECG














Compatible supplies and accessories for the portfolio of CardioPulse ECG 12-lead systems. Note: CardioPulse is only available for sale in the USA.



CardioPulse Prime



CardioPulse Go

	Order Code	Description	For Use With	Qty
<b>Paper</b>				
	010-2280-00	PAPER,THERM,Z-FOLD,210 MM,100SHX10PKS	Prime	1000 sheets
	010-2281-00	PAPER,THERMAL,Z-FOLD,114 MM,60SHX10PKS	Go	600 sheets
<b>Lead Wires, Cables, Adapters</b>				
	700-0386-00	CABLE,PATIENT,AHA,SNAP TYPE,CARDIOPULSE GO	Go	Each
	700-0387-00	CABLE,PATIENT,10-WIRE,AHA,BANANA PLUG,CP GO	Go	Each
	700-0388-00	CABLE,PATIENT,10-WIRE,AHA,CLIP,CP PRIME,1/BG	Prime	Each
	700-0389-00	CABLE,PATIENT,10-WIRE,AHA,BANANAPLUG,CP PRIME	Prime	Each
	700-0399-00	ADAPTER,D15,PATIENT CABLE COMPATIBLE,CP GO	Go	Each
<b>Tab Electrodes, for use with electrode clips</b>				
	700-0400-00	ELECTRODES,WHITE,DISPOSABLE,500 PCS,1/BOX		Box of 500
	SU00C	Disposable tab electrodes for resting ECG set		Box of 100
	58100	Easiclip, Skintact, tab electrode clips <i>For use with banana 4mm sized plugs</i>		Pack of 10
<b>Barcode Scanner</b>				
	010-2282-00	BARCODE SCANNER,LASER,INSTRUCTIONS,CARDIOPULSE		Each

# SUPPLIES & ACCESSORIES

## RESTING ECG



## CardioExpress resting ECG

## CardioDirect resting ECG

CardioExpress SL6A, SL12A and SL18A provide fast and accurate 12-lead ECG acquisition. Advanced design and ergonomics enable the operator to quickly and easily control the unit to meet individual requirements. They provide storage direct to the units and connectivity via Sentinel, Cardiology Information Management System (CIMS).



CardioExpress

## Resting ECG systems 6, 12 and 18 leads



CardioDirect 12

Order Code	Description	For Use With	Qty
<b>ECG Patient Cables &amp; Hook-up Sets</b>			
MS1R-107584-AO	ECG cable, grabber connector, AHA	CardioExpress	Each
MS1R-107048-AO	ECG cable, 4MM, banana connector, AHA	CardioExpress	Each
010-2106-00	ECG cable, banana plugs, AHA	CardioExpress SL18A	Each
010-2107-00	ECG CABLE, SL18A, IEC, BANANA PLUGS	CardioExpress SL18A	Each
MS1R-106902-AO	ECG cable, 4MM, banana connector, IEC	CardioExpress	Each
25-0004	10 electrode patient cable, 4mm, EU	CD12	Each
25-0005	10 electrode patient cable, pinch, EU	CD12	Each
M15R-040172	CLIP/SNAP, BANANA SOCKET ADAPTERS, 10 PCS	CardioExpress	Each
MS1R-107449-AO	SNAP/BANANA SOCKET ADAPTERS (10 PCS/SET)	CardioExpress	Each
M15R-040163	ADULT PRECORDIAL SUCTION ELECTRODES (6 PCS/SET)	CardioExpress	Each
M15R-040168	CHILD PRECORDIAL SUCTION ELECTRODES (6PCS/SET)	CardioExpress	Each
012-0907-00	ELECTRODE,ADULT LIMB CLAMP, QTY 4	CardioExpress	Each
58100	Easiclip, Skintact, 10/pk		10/Pack
<b>AHA CardioExpress SL18A lead wires</b>			
012-0963-00	RA ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0964-00	LA ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0965-00	V5R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0966-00	V4R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0967-00	V3R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0968-00	V1+V2+V3 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0969-00	V4+V5+V6 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0970-00	V7 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0972-00	V9 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0973-00	LL ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0974-00	RL ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
<b>ECG Patient Cables &amp; Hook-up Sets</b>			
<b>IEC CardioExpress SL18A lead wires</b>			
012-0975-00	R (RA) ECG Lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0976-00	L (LA) ECG LEADWIRE FORSL-18A SAMPLING BOX IEC	CardioExpress SL18A	Each
012-0977-00	C5R ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0978-00	C4R ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0979-00	C8 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0980-00	C9 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0981-00	C7 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0982-00	C1+C2+C3 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0983-00	C4+C5+C6 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0984-00	C3R ECG lead wire for SL-18A sampling box,, IEC	CardioExpress SL18A	Each
012-0985-00	F (LL) ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0986-00	N (RL) ECG LEADWIRE FORSL-18A SAMPLING BOX IEC	CardioExpress SL18A	Each

Indicates discontinued part, available while supplies last

## CardioExpress resting ECG

## CardioDirect resting ECG

CardioExpress SL6A, SL12A and SL18A provide fast and accurate 12-lead ECG acquisition. Advanced design and ergonomics enable the operator to quickly and easily control the unit to meet individual requirements. They provide storage direct to the units and connectivity via Sentinel, Cardiology Information Management System (CIMS).



CardioExpress

## Resting ECG systems 6, 12 and 18 leads



CardioDirect 12

Order Code	Description	For Use With	Qty
<b>Miscellaneous ECG Accessories</b>			
<b>Recording Paper</b>			
MS1R-107371-A0	Z-fold, 210 MM x 295 MM, 100 sheets	CardioExpress	Each
MS1R-32462	Z-fold, 210 MM x 295 MM, 200 sheets	CardioExpress	Each
MS1R-107560-A0	Z-fold, 110 MM x 140 MM, 144 sheets	CardioExpress	Each
MS1R-107451-A0	Z-fold, 215 MM x 280MM, 100 sheets	CardioExpress	Each
MS1R-32461	Roll, 210 MM x 30M	CardioExpress	Each
010-1897-00	RECORDING PAPER,ROLL,80MM X 20 M	CardioExpress	Each
MS1R-32463	Roll, 215 MM x 30M	CardioExpress	Each
<b>Others</b>			
010-2127-00	Carry case, SL6A	CardioExpress SL6A	Each
010-2128-00	Carry case, SL12A	CardioExpress SL12A	Each
690-0275-00	Barcode scanner, one dimensional	CardioExpress	Each
690-0276-00	Barcode scanner, two dimensional	CardioExpress	Each
SU00C	Disposable tab electrodes for resting ECG set	CD12USB	Each

Indicates discontinued part, available while supplies last



# SUPPLIES & ACCESSORIES

## Non-Medical



**Non-Medical Accessories**

Spacelabs carries a convenient assortment of Non-Medical supplies to help support your use of our monitoring systems.

Order Code	Description	For Use With	Qty
			

Note: Bumper is black. Yellow color used in image to highlight design and fit of the bumper.

**Bumper Kits**

010-2350-00	Aria Tele Bumpers	\$395.00	5/Box
050-0816-00	Qube Mini Bumper Kit (black)	\$80.00	Each



**Cable Management Systems**

016-1027-00	CORD CADDY, 4 PAIRS OF HOOKS	\$224.89	Each
016-1028-00	CORD CADDY, 6 PAIRS OF HOOKS	\$265.78	Each
016-1029-00	MOUNT,DOWNPOST,GCX,12 INCH (required to mount caddy)	\$65.52	Each

**AV and Networking Cables**

2287169	TRIPP LITE 10' HDMI M/M CABLE	\$19.14	Each
2195215	DP mini to DP cable, male/male, 10ft	\$26.14	Each
1973897	DP mini to DP cable, male/male, 6ft	\$20.22	Each
384343	USB 2.0 A to B Printer Cable, 10ft	\$11.84	Each
1561876	DP to DP cable, male/male, 10ft	\$28.82	Each
54403	DP to DP cable, male/male, 15ft	\$51.82	Each

**Computer Supplies**



5491141	NETGEAR GS308V3 8PRT SWITCH	\$43.10	Each
5615760	NETGEAR GS305 5PRT SWITCH	\$35.26	Each
3006282	LOGITECH MK270 WRLS KBD/MS SET	\$47.50	Each
010-2327-00	Poly Blackwire 5210 Monaural headset for use with SafeNSound	\$120.66	Each

**Barcode Scanners**



690-0275-00	Zebra LS2208 Barcode Scanner	\$250	Each
690-0276-00	Honeywell 1900GSR-2USB-C Barcode Scanner	\$400	Each

### Non-Medical Accessories

Spacelabs carries a convenient assortment of Non-Medical supplies to help support your use of our monitoring systems.

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----



#### Displays and Mounts

210326	ERGOTRON DS100, DUAL DISP MOUNT	\$668.66	Each
6358609	NEC 43" LED-BACKLIT LCD DISPLAY	\$1,575.74	Each
6358610	NEC 50" LED-BACKLIT LCD DISPLAY	\$2,064.32	Each



#### Exergen TAT and Mounts: for use with Xprezzon, Qube, Qube Mini, DM4

010-2157-00	Exergen TAT 5000S, with coiled USB cable	Call for price	Each
134201	Exergen TAT 5000S Wall mount, black	\$67.20	Each
134305	Exergen TAT 5000S locking wall mount, white	\$54.60	Each
134306	Exergen TAT 5000S locking wall mount with cap holder, white	\$63.00	Each

#### Universal Power Supplies and Network Switches



317698	APC,CS500,UPS,300WATT,500VA	\$273.40	Each
6064874	APC,BE850G2,UPS,450WATT,850VA	\$253.08	Each
5491141	NETGEAR GS308V3 8PRT SWITCH	\$43.10	Each
5615760	NETGEAR GS305 5PRT SWITCH	\$35.26	Each