

## Cardiac Output



# CARDIAC OUTPUT

## Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Cardiac Output Cables & Probes




Cables			
700-0027-00	TruLink Cardiac Output Adapter Cable	30 cm/12 in.	1 Each
306655-003	Cardiac Output cable compatible with Baxter Edwards/Abbott injectate temperature probes and most models of thermodilution catheters		1 Each
306655-004	Cardiac Output cable compatible with Merit Medical, Becton Dickinson/Datex-Ohmeda and Argon Medical injectate temperature probes and most models of thermodilution catheters		1 Each

*Note: Spacelabs is working to release new cables to replace the above which have been discontinued by the manufacturer. The above remain available while supplies last.*



Injectate Probes			
690-0023-00	In-line injectate temperature probe Compatible with Datex-Ohmeda SP4045		1 Each
690-0024-00	Reference solution in-line injectate temperature probe. Compatible with Datex-Ohmeda SP5030		1 Each
690-0033-00	In-line injectate temperature probe. Compatible with Datex-Ohmeda SP4042. Not compatible with Abbott or Baxter in-line injectate probes or delivery systems.		1 Each

### Cardiac Output Injectate Kits

Closed Loop Injectate Kits			
690-0293-00		SP5045 Critikit In-Line Sensor, Control Syringe closed loop injectate kit, sterile	5/Box