

Eclipse Mini and Eclipse Pro



Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.

Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code	Description	Length/Style	Qty
------------	-------------	--------------	-----

Eclipse Pro

Leadwire Sets

700-0378-00	3 Lead snap on with yoke Short	<table><tr><th>Snap</th><th>Color</th><th colspan="2">Short</th><th colspan="2">Medium</th><th colspan="2">Long</th></tr><tr><th></th><th></th><th>mm</th><th>Inch</th><th>mm</th><th>Inch</th><th>mm</th><th>Inch</th></tr><tr><td>C5</td><td>Red</td><td>205</td><td>8</td><td>510</td><td>20</td><td>865</td><td>34</td></tr><tr><td>M</td><td>Yellow</td><td>255</td><td>10</td><td>560</td><td>22</td><td>915</td><td>36</td></tr><tr><td>C5R</td><td>Green</td><td>205</td><td>8</td><td>510</td><td>20</td><td>865</td><td>34</td></tr></table>	Snap	Color	Short		Medium		Long				mm	Inch	mm	Inch	mm	Inch	C5	Red	205	8	510	20	865	34	M	Yellow	255	10	560	22	915	36	C5R	Green	205	8	510	20	865	34	Each
Snap	Color	Short		Medium		Long																																					
		mm	Inch	mm	Inch	mm	Inch																																				
C5	Red	205	8	510	20	865	34																																				
M	Yellow	255	10	560	22	915	36																																				
C5R	Green	205	8	510	20	865	34																																				
700-0391-00	3 Lead snap on with yoke Medium		Each																																								
700-0390-00	3 Lead snap on with yoke Long		Each																																								
700-0379-00	4 Lead snap on with yoke Short	<table><tr><th>Snap</th><th>Color</th><th colspan="2">Short</th><th colspan="2">Medium</th><th colspan="2">Long</th></tr><tr><th></th><th></th><th>mm</th><th>Inch</th><th>mm</th><th>Inch</th><th>mm</th><th>Inch</th></tr><tr><td>C5</td><td>Red</td><td>230</td><td>9</td><td>510</td><td>20</td><td>865</td><td>34</td></tr><tr><td>C3</td><td>Orange</td><td>255</td><td>10</td><td>560</td><td>22</td><td>915</td><td>36</td></tr></table>	Snap	Color	Short		Medium		Long				mm	Inch	mm	Inch	mm	Inch	C5	Red	230	9	510	20	865	34	C3	Orange	255	10	560	22	915	36	Each								
Snap	Color	Short		Medium		Long																																					
		mm	Inch	mm	Inch	mm	Inch																																				
C5	Red	230	9	510	20	865	34																																				
C3	Orange	255	10	560	22	915	36																																				
700-0394-00	4 Lead snap on with yoke Medium		Each																																								
700-0393-00	4 Lead snap on with yoke Long	<table><tr><td>M</td><td>White</td><td>255</td><td>10</td><td>560</td><td>22</td><td>915</td><td>36</td></tr><tr><td>C1</td><td>Brown</td><td>230</td><td>9</td><td>510</td><td>20</td><td>865</td><td>34</td></tr></table>	M	White	255	10	560	22	915	36	C1	Brown	230	9	510	20	865	34	Each																								
M	White	255	10	560	22	915	36																																				
C1	Brown	230	9	510	20	865	34																																				
700-0396-00	10 Lead snap on with yoke <i>IEC color coding generally used in European Union</i>	<table><tr><th>Snap</th><th colspan="2">Length</th></tr><tr><th></th><th>mm</th><th>Inch</th></tr><tr><td>L</td><td>610</td><td>24</td></tr><tr><td>F</td><td>585</td><td>23</td></tr><tr><td>C6</td><td>635</td><td>25</td></tr><tr><td>C5</td><td>535</td><td>21</td></tr></table>	Snap	Length			mm	Inch	L	610	24	F	585	23	C6	635	25	C5	535	21	Each																						
Snap	Length																																										
	mm	Inch																																									
L	610	24																																									
F	585	23																																									
C6	635	25																																									
C5	535	21																																									
700-0377-00	10 Lead snap on with yoke <i>AHA color coding generally used in USA, outside of EU</i>	<table><tr><td>C4</td><td>460</td><td>18</td></tr><tr><td>C3</td><td>355</td><td>14</td></tr><tr><td>C2</td><td>330</td><td>13</td></tr><tr><td>C1</td><td>330</td><td>13</td></tr><tr><td>N</td><td>380</td><td>15</td></tr><tr><td>R</td><td>610</td><td>24</td></tr></table>	C4	460	18	C3	355	14	C2	330	13	C1	330	13	N	380	15	R	610	24	Each																						
C4	460	18																																									
C3	355	14																																									
C2	330	13																																									
C1	330	13																																									
N	380	15																																									
R	610	24																																									

712-1052-01	ECLIPSE PRO ECG SENSORS <i>Not available for sale in USA</i>	Box of 10
-------------	---	-----------

Accessories

706-0207-00	Eclipse Pro belt clip	Each
073-0267-00	7 Patient Diary, English CE0123	Each
392196-001	Adhesive tape remover pads, 100/box	100/Box

Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.


Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code		Description	Length/Style	Qty
Eclipse Mini				
	Eclipse Mini Supplies			
	712-1053-01	ECLIPSE MINI ECG SENSORS <i>Not available for sale in USA</i>		Box of 10
	712-1053-02	ECLIPSE MINI ECG SENSORS, DEFIB SAFE Now available for sale in USA		Box of 10