


## Non-Invasive Blood Pressure Cuffs



TruLink® Nylon Reusable

TruLink reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. TruLink cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, TruLink cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. TruLink cuffs are naturally latex and PVC-free.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

TruLink Nylon Reusable BP Cuffs

Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1111-10	Single tube, HP	Infant	8 - 14 cm	1 Each
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each



Bayonet

Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



### TruLink® Soft Disposable

TruLink fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### TruLink Soft Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1113-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



#### Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



### TruLink® Vinyl Disposable

TruLink vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### TruLink Vinyl Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors*

715-1112-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1122-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Bayonet

#### Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



Blood Pressure Cuff Accessories

Spacelabs non-invasive blood pressure cuff accessories are compatible with all Spacelabs and TruLink branded non-invasive blood pressure cuffs. Choose from a variety of hoses, adapters and conversion accessories which offer ease-of-use and flexibility.

Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

Non-Invasive Blood Pressure Hose



Single Tube Hose

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

714-0018-02	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-02	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0019-02	90496/91496	Neonatal	275 cm/9 ft	1 Each

Double Tube to Single Tube Hose

For use with Command Modules 90496/91496

714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each



Hose Connectors & Adapters



Slip Lock Fitting

Fitting for Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.

714-0038-00	3mm ID Adult NIBP hoses	1 Each
714-0047-00	4mm ID Adult NIBP hoses	1 Each
714-0046-00	5mm ID Adult NIBP hoses	1 Each

Indicates discontinued part, available while supplies last

### UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

#### UltraCheck Nylon Reusable BP Cuffs

##### Single Tube



Plastic Bayonet

US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag



Metal Bayonet

US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag



Threaded Screw

US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag



Female  
Subminiature  
Twist Lock

US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag

Indicates discontinued part, available while supplies last

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

UltraCheck Nylon Reusable BP Cuffs



Twist Lock

Single Tube

US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag

No Fitting



US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag



Manual Bulb

US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
US2938MB	Single tube, MB w/gauge	Adult Long	29 - 38 cm	1 Each
US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each

Indicates discontinued part, available while supplies last








### UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

#### UltraCheck Nylon Reusable BP Cuffs

Double Tube					
 Threaded Screw	UD1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
	UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
	UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
 Subminiature Twist Lock	UD0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag
	UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
	UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
	UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag
 Male & Female Mated	UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
	UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
	UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
	UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag
 No Fitting	UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
	UD1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
	UD2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
	UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
	UD4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag
 Manual Bulb	UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
	UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
	UD3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
	Indicates discontinued part, available while supplies last				



### SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

#### SoftCheck Fabric Disposable BP Cuffs

##### Single Tube

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box



ST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
ST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
ST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag



Indicates discontinued part, available while supplies last

### SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Fabric Single Patient Multiple Use BP Cuffs



No Fitting

Single Tube					
ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag	
ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag	
ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag	
ST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag	



Threaded Screw

Double Tube					
DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag	
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box	



Twist Lock

DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box	
DT2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag	



Manual Bulb

DT2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Box	
-------------	-----------------	-------	------------	-------	--

No fitting



DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag	
-----------	-------------------------	------------	------------	-------	--

Indicates discontinued part, available while supplies last

### SoftCheck® Vinyl Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Vinyl Disposables BP Cuffs



Plastic Bayonet



Male  
Twist Lock

Single Tube					
VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag	
VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag	
VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag	
VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag	
VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box	
VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag	
VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag	
VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag	
VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag	
VST1320MQML5-20	Single tube, MQ	Child	13 - 20 cm	5/Bag	
VST1826MQML5-20	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag	
VST2635MQML5-20	Single tube, MQ	Adult	26 - 35 cm	5/Bag	

Indicates discontinued part, available while supplies last

# BLOOD PRESSURE CUFFS

## UltraCheck® & SoftCheck® Curve

### UltraCheck® Curve Nylon Reusable SoftCheck® Curve Fabric Disposable

Your patients come in all shapes and sizes. That's why we offer the Curve cuff for extra-large adults. Its conical-shaped design provides a better fit for the contours of your patient's arm. This can increase patient comfort while providing a more accurate reading.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

#### UltraCheck Curve Nylon Reusable BP Cuff

Single Tube				
BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MB	Single tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each
Double Tube				
BRUD3854	Double tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MB	Double tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

#### SoftCheck Curve Fabric Disposable BP Cuff

Single Tube				
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Box
BRST3854HP-05	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	5/Box

Indicates discontinued part, available while supplies last



Plastic Bayonet  
(HP)

Metal Bayonet  
(HP)

Female Screw  
Threaded (DM)

Female Locking  
(MQ)

Male Locking  
(MQML)

Mated Locking  
(MQMF)

Manual Bulb  
w/Gauge  
(MB)

# BLOOD PRESSURE CUFFS SoftCheck® Fabric Neonatal

## SoftCheck® Fabric Neonatal Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

### SoftCheck Fabric Neonatal Disposable BP Cuffs



Micro-taper  
slip-fit

#### Single Tube

NN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
NN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
NN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
NN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
NN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box



Male Luer

NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



Plastic Bayonet

NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Indicates discontinued part, available while supplies last



#### Double Tube

NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box
----------	-----------------	--------	---------	--------







# BLOOD PRESSURE CUFFS

## SoftCheck® Vinyl Neonatal

### SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



	Order Code	Description	Size	Range/ Length	Qty
<b>SoftCheck Vinyl Neonatal Disposable BP Cuffs</b>					
 <b>Micro-taper slip-fit</b>	<b>Single Tube</b>				
	VNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
 <b>SnapQuik</b>	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 <b>Male Luer</b>	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 <b>Plastic Bayonet</b>	VNN1ST-HP-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN3ST-HP-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-HP-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-HP-10	Single tube, MQ	Size 4	8-15 cm	10/Box
 <b>Female Twist Lock</b>	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box
 <b>Male Luer</b>	<b>Double Tube</b>				
	VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
	VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
	VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
	VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
	VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box