


# Cardiac Output



# CARDIAC OUTPUT




## Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.






Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Cardiac Output Cables & Probes

Cables			
	700-0027-00	TruLink Cardiac Output Adapter Cable. Connects CO cable 30 cm/12 in. to 91496 command module.	1 Each
	306655-003	Cardiac Output cable compatible connects with 690-0033-00 temperature probe. 700-0027-00 adapter required.	1 Each
	306655-004	Cardiac Output cable compatible connects with 690-0023-00 or 690-0024-00 temperature probe. 700-0027-00 adapter required.	1 Each


*Indicates discontinued part, available while supplies last. Spacelabs is releasing new parts to replace the 306655-003 and 306655-004. See July 2024 catalog when available.*

### Injectate Probes

	690-0023-00	In-line injectate temperature probe. Connects to 306655-004 cable and 690-0293-00 injectate kit integrated thermistor.	1 Each
	690-0024-00	Ice bath temperature probe. Connects to 306655-004 cable and needle end is submerged in an ice bath.	1 Each
	690-0033-00	In-line injectate temperature probe. Connects to 306655-003 cable and 690-0293-00 injectate kit integrated thermistor.	1 Each

### Cardiac Output Injectate Kits

#### Closed Loop Injectate Kits

690-0293-00		SP5045 Critikit In-Line Sensor, Control Syringe closed loop injectate kit, sterile	5/Box
-------------	---	--	-------

# CARDIAC OUTPUT

## Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### 91496-C Module, Cardiac Output Connection guide

