


DM4

VITAL SIGNS MONITOR



VITAL SIGNS MONITOR

Ultraview DM4









Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

The DM4 uses TruLink branded NIBP cuffs (single tube, bayonet fitting) and air hose (714-0018-00)



91331

Order Code	Description	Length/Size	Qty
Supplies unique to the DM4 monitor			
LNCS Adapter Cable			
 010-2227-00	DM4 Masimo LNCS SpO ₂ cable compatible with Masimo LNCS sensors, see section on Masimo LNCS sensors	305 cm/10 ft	1 Each
RD SET Adapter Cable			
 010-2193-00	DM4 Masimo RD SET SpO ₂ cable compatible with Masimo RD SET sensors, see section on RD SET sensors	305 cm/10 ft	1 Each
Nellcor OxiMax Cable			
 010-2192-00	DM4 Nellcor OxiMax SpO ₂ cable compatible with Nellcor OxiMax sensors, see section on Nellcor sensors	305 cm/10 ft	1 Each
Filac			
 010-2174-00	Filac dock and insulation chamber		1 Each
 010-2175-00	Filac temperature probe		1 Each
 010-2228-00	Filac temperature oral probe covers, 20/bx		20/box 25/boxes
Exergen			
 010-2157-00	Temporal Artery Thermometer		1 Each
NeoNatal NIBP supplies			
 714-0087-00	NIBP,HOSE ASSY,NEONATE,DM4 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	1 Each