

Pulse Oximetry



TruLink SpO₂ Supplies & Accessories
Spacelabs TruLink Technology

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO₂ supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can



	Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology				
	Cables			
	700-0030-00	TruLink SpO ₂ cable	305 cm/10 ft.	1 Each
	Reusable sensors			
	015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
	015-0661-00	TruLink SpO ₂ Y multi-site sensor <u>Includes:</u> (1) Sensor (1) ear clip (10) adult finger/neonate feet/hand wraps (1) pediatric finger or infant toe wraps	Universal Adult>30 kg Pediatric 10-50kg Infant 3-20kg Neonatal <3kg	1 Each
Disposable Sensors				
	015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
	015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
	015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0666-00	TruLink SpO ₂ cloth sensor	Infant 3-20 kg/7-44 lbs.	24/Box
	015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
	Dixtal Cables			
	700-0029-00	TruLink SpO ₂ cable compatible with Novametrix® sensors	30 cm/12 in.	1 Each

Masimo® SpO2 Supplies & Accessories

LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.



Order Code	Description	Length/Size	Qty
LNCS Masimo SET Technology			
LNCS Adapter Cable			
 700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft	1 Each
LNCS Reusable Sensors			
 690-0247-00	Masimo LNCS DC-I, SpO ₂ digit sensor	Adult > 30 kg	1 Each
 690-0248-00	Masimo LNCS DC-IP, SpO ₂ digit sensor	Pediatric 10-50 kg	1 Each
 690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO ₂ ear sensor	Adult > 30 kg	1 Each
 690-0291-00	Masimo LNCS Y-I multisite sensor	> 1 kg	1 Each
LNCS Masimo SET Technology			
LNCS Disposable Sensors			
 690-0251-00	Masimo LNCS Amtx, adhesive SpO ₂ sensor, 18 in.	Adult > 30 kg	20/Box
 690-0252-00	Masimo LNCS Pmtx, adhesive SpO ₂ sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
 690-0258-00	Masimo LNCS Inf-L, adhesive SpO ₂ sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
 690-0259-00	Masimo LNCS Neo-L, adhesive SpO ₂ sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box

Masimo® SpO₂ Supplies & Accessories RD SET Technology


For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.



Order Code	Description	Length/Size	Qty
RD SET Masimo Technology			
RD SET Adapter Cable			
 010-1649-00	Masimo RD SET SpO ₂ Cable	12 ft.	1 Each
RD SET Reusable Sensors			
 690-0267-00	Masimo RD SET DCI reusable SpO ₂ sensor	Adult > 30 kg	1 Each
 690-0270-00	Masimo RD SET DCI-P reusable SpO ₂ sensor	Pediatric 10-50 kg	1 Each
 690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO ₂ sensor	Adult > 30 kg	1 Each
 690-0281-00	Masimo RD SET Y-I multi-site reusable SpO ₂ sensor	> 1 kg	1 Each
RD SET Masimo Technology			
RD SET Adhesive Disposable Sensors			
 690-0283-00	Masimo RD SET Adt, adhesive SpO ₂ sensor	Adult, > 30 kg finger or toe	20/Box
 690-0285-00	Masimo RD SET Inf, adhesive SpO ₂ sensor	Infant, 3—20 kg, thumb or great toe	20/Box
 690-0286-00	Masimo RD SET Neo, adhesive SpO ₂ sensor	Neonatal <3 kg, hand or foot Adult > 40 kg finger or toe	20/Box

Nellcor® SpO₂ Supplies & Accessories Nellcor™ OxiMax® Technology

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor. Cables and connectors are color-coded purple. **Note, use only genuine Nellcor sensors.**



Order Code	Description	Length/Size	Qty
Nellcor OxiMax Technology			
Cable			
	700-0792-00 TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
Reusable Sensors			
Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien), now part of Medtronic.			
	690-0003-01 Nellcor DS-100A, finger SpO ₂ sensor	Adult	1 Each
	690-0004-01 Nellcor OXI-A/N, SpO ₂ sensor	Adult/Neonatal	1 Each
	690-0039-01 Nellcor OXI-P/I, SpO ₂ sensor	Pediatric/Infant	1 Each
	690-0215-01 Nellcor DS-Y, multi-site SpO ₂ sensor	Universal	1 Each