


# Temperature










# TEMPERATURE

## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
<b>Monitor Adapter Cables</b>			
<b>Adapters</b>			
 700-0031-00	TruLink Dual Temperature Cable. *Required for all reusable probes and as directed for single-patient-use disposable system	30 cm/12 in.	1 Each
<b>Reusable Probes/Sensors</b>			
<b>Skin Surface</b>			
 690-0295-00	YSI 400 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0296-00	YSI 700 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0297-00	YSI 700 series compatible skin surface sensor temperature probe, pediatric size	305 cm/10 ft	1 Each
<b>Esophageal / Rectal</b>			
 690-0294-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
 690-0298-00	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
 690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each













Indicates discontinued part, available while supplies last

# TEMPERATURE

## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
<b>Single-Patient use temperature system</b>			
<b>Extension cables for single use probes</b>			
 700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors, adapter required.	305 cm/10 ft	1 Each
 700-0401-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors. Direct connection to Command module, no adapter required.	305 cm/10 ft	1 Each
<b>Skin</b>			
 701-4000-00	Skin Temperature Sensor		20/Box
<b>Esophageal/ Stethoscope</b>			
 701-4009-00	9 French, YSI 400 series compatible		20/Box
 701-4012-00	12 French, YSI 400 series compatible		20/Box
 701-4018-00	18 French, YSI 400 series compatible		20/Box
<b>Esophageal/ Rectal</b>			
 701-4109-00	9 French, YSI 400 series compatible		20/Box
 701-4112-00	12 French, YSI 400 series compatible		20/Box
<b>Foley catheter with temperature probe</b>			
 701-4214-00	14 French, YSI 400 series compatible		20/Box
 701-4216-00	16 French, YSI 400 series compatible		20/Box
<b>Tympanic temperature sensor</b>			
 701-4300-00	Adult, YSI 400 series compatible		20/Box
<b>Exergen</b>			
<i>For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) &amp; DM4 (91331) Monitors</i>			
 010-2157-00	Temporal Artery Thermometer		1 Each

# TEMPERATURE

## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
<b>Reusable Temperature Probes</b>			
690-0294-00	YSI400 Adult Esophageal/Rectal		
690-0298-00	YSI700 Adult Esophageal/Rectal		
690-0295-00	YSI400 Adult Skin		
690-0296-00	YSI700 Adult Skin		
690-0297-00	YSI700 Pediatric Skin		
<b>Single-Patient Use Temperature Probes</b>			
701-4009-00	w/ Stethoscope 9Fr		
701-4012-00	12Fr		
701-4018-00	18Fr		
701-4109-00	Esophageal/Rectal 9Fr		
701-4112-00	12Fr		
701-4214-00	w/ Foley Catheter 14Fr		
701-4216-00	16Fr		
701-4000-00	Skin		
701-4300-00	Tympanic		

91496 w/ T1-2	700-0031-00 Dual temp cable	→	Reusable Temperature Probes
91496-L w/ T3-4	700-0401-00 Reusable extension cable	→	Single-Patient Use Temperature Probes
91496-L w/ T3-4	700-0031-00 Dual temp cable	→	Single-Patient Use Temperature Probes
91496-L w/ T3-4	700-4000-00 Reusable extension cable	→	Single-Patient Use Temperature Probes