



MAY 2024

Supplies & Accessories Product Catalog

An OSI Systems Company


© 2024 Spacelabs Healthcare Specifications subject to change without notice.

Spacelabs Healthcare reserves the right to make changes to the product at any time without notice. All brands and product names and trademarks of their respective owners. Product may not be available for sale in all locations. Check with your Spacelabs representative.

SPACELABS
HEALTHCARE

www.spacelabshealthcare.com

Non-Invasive Blood Pressure Cuffs



TruLink® Nylon Reusable

TruLink reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. TruLink cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, TruLink cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. TruLink cuffs are naturally latex and PVC-free.



**Range/
Length** **Qty**

Order Code **Description**

Size

TruLink Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

715-1111-10	Single tube, HP	Infant	8 - 14 cm	1 Each
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each



Bayonet

Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



BLOOD PRESSURE CUFFS
TruLink® Soft

TruLink® Soft Disposable

TruLink fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

TruLink Soft Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors



715-1113-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each

TruLink® Vinyl Disposable

TruLink vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

TruLink Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

715-1112-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1122-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Bayonet

Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



Blood Pressure Cuff Accessories

Spacelabs non-invasive blood pressure cuff accessories are compatible with all Spacelabs and TruLink branded non-invasive blood pressure cuffs. Choose from a variety of hoses, adapters and conversion accessories which offer ease-of-use and flexibility.

Order Code	Description	Size	Length	Qty
------------	-------------	------	--------	-----

Non-Invasive Blood Pressure Hose

Single Tube Hose



For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

714-0018-02	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-02	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0019-02	90496/91496	Neonatal	275 cm/9 ft	1 Each

Double Tube to Single Tube Hose



For use with Command Modules 90496/91496

714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each
-------------	-------------	----------	-------------	--------

Hose Connectors & Adapters

Slip Lock Fitting



Fitting for Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.





714-0038-00	3mm ID Adult NIBP hoses	1 Each
714-0047-00	4mm ID Adult NIBP hoses	1 Each
714-0046-00	5mm ID Adult NIBP hoses	1 Each

Indicates discontinued part, available while supplies last

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.






	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
Single Tube					
 Plastic Bayonet	US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
 Metal Bayonet	US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
	US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
	US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
	US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
	US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
	US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
	US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
	US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag
 Threaded Screw	US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
	US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
	US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
	US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
	US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
 Female Subminiature Twist Lock	US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
	US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
	US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
	US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
	US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
<i>Indicates discontinued part, available while supplies last</i>					

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.








	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
Single Tube					
 Twist Lock	US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
	US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
	US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
	US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
 No Fitting	US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
 Manual Bulb	US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
	US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
	US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each
	Indicates discontinued part, available while supplies last				

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
Double Tube					
	UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
Threaded Screw	UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
	UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
	UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
	UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
Subminiature Twist Lock	UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
	UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
	UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
	UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag
	UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
	UD2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
	UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
	No Fitting	UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm
	UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag

Indicates discontinued part, available while supplies last

SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



**Range/
Length Qty**

Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box



Plastic Bayonet



**Female
Subminiature**





Indicates discontinued part, available while supplies last

BLOOD PRESSURE CUFFS SoftCheck® Fabric

SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.





	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Fabric Single Patient Multiple Use BP Cuffs					
 No Fitting	Single Tube				
	ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
	ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
 Threaded Screw	Double Tube				
	DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box
 Twist Lock	DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box
	 No fitting	DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm

Indicates discontinued part, available while supplies last

SoftCheck® Vinyl Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Vinyl Disposables BP Cuffs					
	Single Tube				
	VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
	VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
					
Plastic Bayonet					
					
Male Twist Lock					
	VST1826MQML5-20	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
	VST2635MQML5-20	Single tube, MQ	Adult	26 - 35 cm	5/Bag

Indicates discontinued part, available while supplies last

BLOOD PRESSURE CUFFS UltraCheck® & SoftCheck® Curve

UltraCheck® Curve Nylon Reusable SoftCheck® Curve Fabric Disposable

Your patients come in all shapes and sizes. That's why we offer the Curve cuff for extra-large adults. Its conical-shaped design provides a better fit for the contours of your patient's arm. This can increase patient comfort while providing a more accurate reading.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

UltraCheck Curve Nylon Reusable BP Cuff

Single Tube				
BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MB	Single tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each
Double Tube				
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MB	Double tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

SoftCheck Curve Fabric Disposable BP Cuff

Single Tube				
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Box
BRST3854HP-05	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	5/Box

Indicates discontinued part, available while supplies last







BLOOD PRESSURE CUFFS SoftCheck® Fabric Neonatal

SoftCheck® Fabric Neonatal Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Fabric Neonatal Disposable BP Cuffs					
 Micro-taper slip-fit	Single Tube				
	NN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	NN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	NN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	NN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
NN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box	
 Male Luer	NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Plastic Bayonet	NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
	NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
	NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box
<i>Indicates discontinued part, available while supplies last</i>					
 Double Tube	NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box







BLOOD PRESSURE CUFFS

SoftCheck® Vinyl Neonatal

SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Vinyl Neonatal Disposable BP Cuffs					
Single Tube					
 Micro-taper slip-fit	VNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
 SnapQuik	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 Male Luer	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Plastic Bayonet	VNN1ST-HP-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN3ST-HP-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-HP-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-HP-10	Single tube, MQ	Size 4	8-15 cm	10/Box
 Female Twist Lock	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box
Double Tube					
 Male Luer	VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
	VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
	VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
	VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
	VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

Pressure Infusors



Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Disposable Pressure Infusor

Single Patient Use Infusor with aneroid gauge - sold as one piece/complete unit



803FGA-24	Unifusor with thumbwheel valve & aneroid gauge	500mL	24/Box
804FGA-24	Unifusor with thumbwheel valve & aneroid gauge	1000mL	24/Box
805FGA-10	Unifusor with thumbwheel valve & aneroid gauge	3000mL	10/Box



803SGA-24	Unifusor with stopcock valve & aneroid gauge	500mL	24/Box
804SGA-24	Unifusor with stopcock valve & aneroid gauge	1000mL	24/Box
805SGA-10	Unifusor with stopcock valve & aneroid gauge	3000mL	10/Box

Disposable Infusor with piston gauge - sold as one piece/complete unit



903TGA-24	Unifusor with thumbwheel valve & piston gauge	500mL	24/Box
904TGA-24	Unifusor with thumbwheel valve & piston gauge	1000mL	24/Box
905TGA-10	Unifusor with thumbwheel valve & piston gauge	3000mL	10/Box



903SGA-24	Unifusor with stopcock valve & piston gauge	500mL	24/Box
904SGA-24	Unifusor with stopcock valve & piston gauge	1000mL	24/Box
905SGA-10	Unifusor with stopcock valve & piston gauge	3000mL	10/Box

Indicates discontinued part, available while supplies last

Unifusor II Disposable Pressure Infusor with Reusable Inflation System

Disposable Infusor - white Disposable bag with reusable inflation system



2803SGA-24	Unifusor II with aneroid gauge & stopcock valve	500mL	24/Box
2804SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000mL	24/Box
2805SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000mL	10/Box



2903SGA-24	Unifusor II with piston gauge & stopcock valve	500mL	24/Box
2904SGA-24	Unifusor II with piston gauge & stopcock valve	1000mL	24/Box
2905SGA-10	Unifusor II with piston gauge & stopcock valve	3000mL	10/Box



Disposable Replacement Parts

2803X-24	Unifusor II Disposable white infusion bag	500mL	24/Box
----------	---	-------	--------

Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Plus Reusable Pressure Infusor

Reusable Infusor - blue bag with white hose



1103SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	500mL	24/Box
1104SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	1000mL	24/Box
1105SGA-10	Unifusor Plus with aneroid gauge & stopcock valve	3000mL	10/Box
1113SGA-24	Unifusor Plus with piston gauge & stopcock valve	500mL	24/Box

Reusable replacement parts



1103XB-05	Unifusor Plus blue reusable bag with white hose	500mL	5/Box
1104XB-05	Unifusor Plus blue reusable bag with white hose	1000mL	5/Box
1105XB-05	Unifusor Plus blue reusable bag with white hose	3000mL	5/Box

Indicates discontinued part, available while supplies last

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with black hose

1105	Unifusor Plus with trigger / gauge assembly	3000mL	1 Each
------	---	--------	--------

Reusable replacement parts

1103X	Unifusor Plus blue reusable bag with black hose	500mL	1 Each
1104X	Unifusor Plus blue reusable bag with black hose	1000mL	1 Each
1105X	Unifusor Plus blue reusable bag with black hose	3000mL	1 Each

Electrocardiography


AAMI COLOR CODING



ECG Combiner Lead Wires & Cables

ECG COMBINER LEAD WIRE SYSTEM AAMI Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Disposable, DIN Shielded, Pinch, AAMI				
700-0006-50	3-Lead	Pinch	127cm /50 in	25 Sets
* 700-0006-32	5 Lead	Pinch	127cm /50 in	25 Sets
* 700-0006-33	5 Lead (V-Lead)	Pinch	127cm /50 in	25 Sets



Disposable, DIN Shielded, Snap, AAMI				
700-0006-61	5 Lead	Snap	127cm /50 in	25 Sets



Reusable, DIN Shielded, Pinch, AAMI				
* 700-0006-35	3 Lead	Pinch	74cm /29 in	1 Set
* 700-0006-37	5 Lead	Pinch	74cm /29 in	1 Set
* 700-0006-38	5 Lead (V-Lead)	Pinch	74cm /29 in	1 Set
700-0010-00	3 Lead	Pinch	130cm/51 in	1 Set
* 700-0010-01	5 Lead	Pinch	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Pinch	130cm/51 in	1 Set



Reusable, DIN Shielded, Snap, AAMI				
* 700-0007-35	3 Lead	Snap	74cm /29 in	1 Set
* 700-0007-37	5 Lead	Snap	74cm /29 in	1 Set
700-0010-19	5-Lead	Snap	130cm /51 in	1 Set

* Product not available for sale into CE Mark Countries (EU)

ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter




Radiolucent, Reusable, DIN Unshielded, Pinch, AAMI				
700-0010-06	5 Lead, Radiolucent	Pinch	122cm/48 in	1 Set
700-0010-08	5 Lead, Radiolucent	Pinch	183cm/72 in	1 Set

ECG Combiner Lead Wires & Cables

ECG COMBINER LEAD WIRE SYSTEM AAMI Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496					
* 700-0008-55	3 Lead	Spacelabs	305 cm/10 ft	1 Cable	
* 700-0008-57	5 Lead	Spacelabs	305 cm/10 ft	1 Cable	
700-0008-32	5 Lead	Spacelabs	457 cm/15 ft	1 Cable	
* 700-0008-56	10 Lead (5+5)	Spacelabs	305 cm/10 ft	1 Cable	

* Product not available for sale into CE Mark Countries (EU)

ECG Combiner Adapters



Spacelabs Telemetry, AAMI					
700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)		1 Adapter	
700-0008-88	5 Lead	Aria Telemetry (96281)		1 Adapter	

ECG Combiner Cables for Standardization




Other manufacturers					
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable	
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable	

Indicates discontinued part, available while supplies last

ECG LEAD WIRES & CABLES AAMI Color Coded - Shielded

ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Order Code Description Length/Size Qty

ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

Snap Sets

Order Code	Description	Length/Size	Qty
700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set



Pinch

Pinch Sets

Order Code	Description	Length/Size	Qty
700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set

ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI



24" Snap Individual Color


Order Code	Description	Length/Size	Qty
012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-34	Black	61 cm/24 in	5/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

Indicates discontinued part, available while supplies last

ECG LEAD WIRES & CABLES AAMI Color Coded - Shielded

ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Cables, Shielded, AAMI



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496

700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

ECG Accessories



Electrodes

717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 pouches	Adult	50/Box
-------------	---	-------	--------

Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable

700-0008-12	3 Lead	Neonatal	1 Cable
-------------	--------	----------	---------

DIN Unshielded, Neonatal Electrodes.

717-0080-00	Pre-wired neonatal electrodes, Pre-gelled, silver/silver chloride, single use	3-Lead	100 Sets/Box
-------------	---	--------	--------------

DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap Sets

015-0732-00	3 Lead	90 cm/35 in	1 Set
015-0736-00	5 Lead	90 cm/35 in	1 Set



Pinch Sets

015-0731-00	3 Lead	90 cm/35 in	1 Set
015-0734-00	5 Lead	90 cm/35 in	1 Set

Bold font indicates a new part number

Telemetry Transmitter Adapters



For use with DIN Unshielded Lead Wires




706-0191-00	Aria Telemetry Grouper (96281)	5/Bag
-------------	--------------------------------	-------

DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded




Order Code	Description	Length/Size	Qty
Telemetry Accessories			
Protective Solutions and Electrodes			
	010-2350-00 Aria Tele Bumper		5/Box
	015-0500-02 Ultraview Telemetry pouch (90347)		100/Box
	717-0083-00 Electrode, Transparent Tape, 5 per pouch X 10 Pouches	Adult	50/Box

DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

SpO₂ Cable & Sensors

Note: Only the below genuine Spacelabs and Nellcor branded SpO₂ sensors are compatible with Spacelabs Telemetry systems

SpO₂ Cable



700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

Reusable SpO₂ Sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

Disposable SpO₂ Sensors



015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last

Electrocardiography


IEC COLOR CODING



ECG Combiner Lead Wires & Cables

ECG COMBINER LEAD WIRE SYSTEM IEC Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Multi-pin

Disposable, DIN Shielded, Multi-Pinch, IEC

700-0006-44	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-53	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets



Reusable, DIN Shielded, Multi-Pinch, IEC

700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set
-------------	-----------------	-------------	-------------	-------



Snap

Reusable, DIN Shielded, Snap, IEC

700-0007-41	3 Lead	Snap	74cm /29 in	1 Set
-------------	--------	------	-------------	-------

Indicates discontinued part, available while supplies last

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)

700-0008-81	3 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------







700-0008-83	5 Lead	IEC	305 cm/10 ft	1 Cable
-------------	--------	-----	--------------	---------

ECG LEAD WIRES & CABLES IEC Color Coded - Shielded

ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.






	Order Code	Description	Length	Qty
ECG Lead Wire Sets, Shielded, Reusable, IEC				
 <p>Snap</p>	Snap Sets			
	700-0007-01	3 Lead, shielded	61 cm/24 in	1 Set
	700-0007-03	3 Lead, shielded	102 cm/40 in	1 Set
	700-0007-09	5 Lead, shielded	61 cm/24 in	1 Set
	700-0007-11	5 Lead, shielded	102 cm/40 in	1 Set
	700-0007-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
	700-0007-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set
 <p>Pinch</p>	Pinch Sets			
	700-0006-01	3 Lead, shielded	61 cm/24 in	1 Set
	700-0006-03	3 Lead, shielded	102 cm/40 in	1 Set
	700-0006-09	5 Lead, shielded	61 cm/24 in	1 Set
	700-0006-11	5 Lead, shielded	102 cm/40 in	1 Set
	700-0006-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
	700-0006-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set
ECG Cables, Shielded, IEC				
	Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496			
	700-0008-09	3 Lead, shielded	305 cm/10 ft	1 Cable
	700-0008-07	5 Lead, shielded	305 cm/10 ft	1 Cable
	700-0008-01	10 Lead, shielded	305 cm/10 ft	1 Cable
<i>Indicates discontinued part, available while supplies last</i>				
ECG Accessories				
	Electrodes			
	717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches	Adult	50/Box

ECG LEAD WIRES & CABLES IEC Color Coded - Shielded

ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Order Code	Description	Length	Qty
Neonatal ECG Supplies & Accessories			
	DIN Unshielded, Neonatal ECG Cable		
700-0008-13	3 Lead		1 Cable
	DIN Unshielded, Neonatal Electrodes		
717-0080-00	Pre-wired neonatal electrodes, pre-gelled, silver/silver chloride, CE	3-Lead	100 Sets/Box

DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY IEC Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC



Snap

Snap Sets			
015-0733-00	3 Lead	90 cm/35 in	1 Set
012-0328-10	5 Lead	61 cm/24 in	1 Set
015-0737-00	5 Lead	90 cm/35 in	1 Set



Pinch Sets			
015-0735-00	5 Lead	90 cm/35 in	1 Set



For use with DIN Unshielded Lead Wires			
706-0190-00	Aria Telemetry Grouper (96281), IEC		5/Bag

Telemetry Accessories



Disposable Pouches and Electrodes			
010-2350-00	Aria Tele Bumper		5/box
015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches	Adult	50/Box

Bold font indicates a new part number

Indicates discontinued part, available while supplies last

DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY IEC Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

SpO₂ Cable & Sensors

SpO₂ Cable



700-0014-00	TruLink SpO ₂ Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO ₂ sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

Note: Only the below genuine Spacelabs and Nellcor branded SpO₂ sensors are compatible with Spacelabs Telemetry systems

Reusable SpO₂ Sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO ₂ sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

Disposable SpO₂ Sensors



015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last

Capnography & Multigas



NomoLine Sampling Lines

NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



Capno Pod (92516) Capnography Module (92517) Multi-Gas Module (92518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

NomoLine LH Nasal CO2 Cannula



704-3817-00	Adult	2 m/7 ft	25/Box
704-3819-00	Infant	2 m/7 ft	25/Box
704-3820-00	Adult w/O2	2 m/7 ft	25/Box
704-4624-00	Adult w/O2	4 m/13 ft	25/Box
704-3821-00	Pediatric w/O2	2 m/7 ft	25/Box



NomoLine LH Nasal/Oral CO2 Cannula

704-3824-00	Adult w/O2	2 m/7 ft	25/Box
704-4628-00	Adult w/O2	4 m/13 ft	25/Box
704-3825-00	Pediatric w/O2	2 m/7 ft	25/Box



NomoLine LH Single Nasal Prong CO2 Cannula

704-3826-00	Adult	2 m/7 ft	25/Box
-------------	-------	----------	--------

NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine LH Airway Adapter Set

704-3814-00	Adult/Pediatric	2 m/7 ft	25/Box
704-4620-00	Adult/Pediatric	4 m/13 ft	25/Box

NomoLine LH (Low Humidity) Procedural Oxygen Masks



NomoLine LH Nasal CO2 Cannula

704-5213-00 ^{USA}	Adult, Medium O2 Flow Rate	3 m/10 ft	30/Box
704-5214-00 ^{USA}	Pediatric, Medium O2 Flow Rate	3 m/10 ft	30/Box
704-5215-00 ^{USA}	Adult, High O2 Flow Rate	3 m/10 ft	30/Box
704-5216-00 ^{USA}	Pediatric, High O2 Flow Rate	3 m/10 ft	30/Box

USA - Parts only available for sale in the USA

NomoLine Sampling Lines

NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



Capno Pod (92516) Capnography Module (92517) Multi-Gas Module (92518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

NomoLine HH Nasal CO2 Cannula



704-3830-00	Adult	2 m/7 ft	25/Box
704-3832-00	Infant	2 m/7 ft	25/Box
704-3833-00	Adult w/O2	2 m/7 ft	25/Box
704-3834-00	Pediatric w/O2	2 m/7 ft	25/Box
704-4933-00 ^{USA}	Infant w/O2	4 m/13 ft	25/Box



NomoLine HH Nasal/Oral CO2 Cannula


704-3837-00	Adult w/O2	2 m/7 ft	25/Box
704-4968-00 ^{USA}	Adult w/O2	4 m/13 ft	25/Box
704-3838-00	Pediatric w/O2	2 m/7 ft	25/Box



NomoLine HH Single Nasal Prong CO2 Cannula

704-3839-00	Adult	2 m/7 ft	25/Box
-------------	-------	----------	--------

NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine HH Airway Adapter Set

704-3827-00	Adult/Pediatric	2 m/7 ft	25/Box
704-3828-00	Adult/Pediatric	3 m/10 ft	25/Box
704-4989-00 ^{USA}	Adult/Pediatric	4 m/13 ft	25/Box
704-4367-00	Infant/Neonatal	2 m/7 ft	25/Box

USA - Parts only available for sale in the USA

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Capno Pod
(92516)

Capnography Module
(92517)

Multi-Gas Module
(92518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Sidestream Monitoring Supplies

Please note: the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.

CO2 Nasal Cannula with male luer connector



*	704-0176-00	Adult	2 m/7 ft	25/Box
**	704-0177-00	Pediatric	2 m/7 ft	25/Box
**	704-0178-00	Infant	2 m/7 ft	25/Box

CO2 Nasal Cannula w/O2 Delivery and male luer connector



*	704-0181-00	Adult	2 m/7 ft	25/Box
CE	704-0181-01	Adult	2 m/7 ft	25/Box
*	704-0180-00	Pediatric	2 m/7 ft	25/Box
**	704-0179-00	Infant	2 m/7 ft	25/Box
**	704-0186-00	Adult	3 m/10 ft	25/Box

CO2 Oral/Nasal Cannula w/O2 Delivery and male luer connector



*	704-0183-00	Adult	2 m/7 ft	25/Box
CE	704-0183-01	Adult	2 m/7 ft	25/Box
**	704-0187-00	Adult	3 m/10 ft	25/Box

Indicates discontinued part, available while supplies last

Sidestream Monitoring Supplies - Intubated patients

Sampling Lines & Airway Adapters, single use



**	704-0191-00	Adult CO2 Sample Line w/ T-Adapter	2 m/7 ft	25/Box
**	704-0192-00	Adult CO2 Sample Line w/ T-Adapter	3 m/10 ft	25/Box



*	166-7121-00	Sidestream sampling line with male luer lock on both ends	3 m/10 ft	50/Box
---	-------------	---	-----------	--------



*	704-0175-00	Sidestream airway adapter/T-Adapter.		25/Box
---	-------------	--------------------------------------	--	--------

* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Capno Pod
(92516)

Capnography Module
(92517)

Multi-Gas Module
(92518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients.



Sensor & Airway Adapters - required for mainstream monitoring

010-1980-00	Mainstream sensor, reusable		1 Each
704-0173-00	Mainstream airway adapter, single use	Adult/Pediatric	25/Box
704-0174-00	Mainstream airway adapter, single use	Infant	10/Box

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module


Order Code	Description	Size	Qty
Sidestream, Non-intubated, single patient use			
CO₂ Nasal Cannula			
* 704-0010-00	Adult		10/Box
* 704-0012-00	Pediatric		10/Box
CO₂ Nasal Cannula w/O₂ Delivery			
* 704-0014-00	Adult		10/Box
CE 704-0014-10	Adult, (for sale into EU/UK only)		10/Box
* 704-0016-00	Pediatric		10/Box
* 704-0015-00	Infant/Neonatal		10/Box
CO₂ Oral/Nasal Cannula			
* 704-0017-00	Adult		10/Box
CO₂ Oral/Nasal Cannula w/O₂ Delivery			
* 704-0020-00	Adult		10/Box
CE 704-0020-10	Adult, (for sale into EU/UK only)		10/Box
* 704-0022-00	Pediatric		10/Box
Sidestream, Intubated, single patient use			
Sidestream - short term up to 12 hours (patient non-humidified)			
* 704-0023-00	Airway adapter + sampling line, ET tube > 4.0 mm	Adult/Ped	10/Box
CE 704-0023-10	Airway adapter + sampling line, ET tube > 4.0 mm	Adult/Ped	10/Box
* 704-0027-00	CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>		10/Box
CE 704-0027-10	CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>		10/Box

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream - long term up to 120 hours (patient humidified)			
 * 704-0025-00	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	Adult/Ped	10/Box
CE 704-0025-10	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	Adult/Ped	10/Box
* 704-0026-00	Airway adapter + sampling line with Nafion	Ped/Infant	10/Box
* 704-0028-00	CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>		10/Box
CE 704-0028-10	CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>		10/Box

Mainstream, for Intubated Patients Only



Mainstream Airway Adapters, for use with CO₂ Sensor 704-0001-00			
* 704-0005-00	CO ₂ airway adapter, single use	Adult	24/Box
CE 704-0005-10	CO ₂ airway adapter, single use	Adult	10/Box
* 704-0004-00	CO ₂ airway adapter, single use	Neonate	24/Box
* 704-0002-00	CO ₂ airway adapter, reusable	Adult	1 Each

* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

MULTIGAS ANALYZER

Spacelabs Modules 90518 & 91518

Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).

These two legacy multi-gas modules are obsolete. However, Spacelabs continues to provide supplies that are available.



(90518)



(91518)

Order Code	Description	Length/Size	Qty
Multigas Supplies & Accessories			
Anesthetic Agent Samples Lines			
015-0128-00	SAMPLE LINE,MALE TO MALE,LUER,SPU,7'X.04",25/BX		25/Box
015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line:0.048 inches ID tubing (90518/91518)	305 cm/10 ft.	50/Box
162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 0.048 in ID tubing (90518/91518)	305 cm/10 ft.	50/Box
Water Traps, Reusable			
717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)		10/Pack
717-0448-00	D-FEND water trap cartridge, Green. Replace at least once every 24 hours or when occluded (91518)		10/Pack
<i>Note: manufacturer (GE) has discontinued these parts. Available while supplies last.</i>			







Pulse Oximetry



TruLink SpO₂ Supplies & Accessories
Spacelabs TruLink Technology

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO₂ supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can












	Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology				
	Cables			
	700-0030-00	TruLink SpO ₂ cable	305 cm/10 ft.	1 Each
	Reusable sensors			
	015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
	015-0661-00	TruLink SpO ₂ Y multi-site sensor <u>Includes:</u> (1) Sensor (1) ear clip (10) adult finger/neonate feet/hand wraps (1) pediatric finger or infant toe wraps	Universal Adult>30 kg Pediatric 10-50kg Infant 3-20kg Neonatal <3kg	1 Each
Disposable Sensors				
	015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
	015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
	015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
	015-0666-00	TruLink SpO ₂ cloth sensor	Infant 3-20 kg/7-44 lbs.	24/Box
	015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
	Dixtal Cables			
	700-0029-00	TruLink SpO ₂ cable compatible with Novametrix® sensors	30 cm/12 in.	1 Each

Masimo® SpO2 Supplies & Accessories

LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.



Order Code	Description	Length/Size	Qty
LNCS Masimo SET Technology			
LNCS Adapter Cable			
 700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft	1 Each
LNCS Reusable Sensors			
 690-0247-00	Masimo LNCS DC-I, SpO ₂ digit sensor	Adult > 30 kg	1 Each
 690-0248-00	Masimo LNCS DC-IP, SpO ₂ digit sensor	Pediatric 10-50 kg	1 Each
 690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO ₂ ear sensor	Adult > 30 kg	1 Each
 690-0291-00	Masimo LNCS Y-I multisite sensor	> 1 kg	1 Each
LNCS Masimo SET Technology			
LNCS Disposable Sensors			
 690-0251-00	Masimo LNCS Amtx, adhesive SpO ₂ sensor, 18 in.	Adult > 30 kg	20/Box
 690-0252-00	Masimo LNCS Pmtx, adhesive SpO ₂ sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
 690-0258-00	Masimo LNCS Inf-L, adhesive SpO ₂ sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
 690-0259-00	Masimo LNCS Neo-L, adhesive SpO ₂ sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box

Masimo® SpO₂ Supplies & Accessories RD SET Technology

For use with Ultraview SL™ Command Module (91496-M). Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.








Order Code	Description	Length/Size	Qty
RD SET Masimo Technology			
RD SET Adapter Cable			
010-1649-00	Masimo RD SET SpO ₂ Cable	12 ft.	1 Each
RD SET Reusable Sensors			
690-0267-00	Masimo RD SET DCI reusable SpO ₂ sensor	Adult > 30 kg	1 Each
690-0270-00	Masimo RD SET DCI-P reusable SpO ₂ sensor	Pediatric 10-50 kg	1 Each
690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO ₂ sensor	Adult > 30 kg	1 Each
690-0281-00	Masimo RD SET Y-I multi-site reusable SpO ₂ sensor	> 1 kg	1 Each
RD SET Masimo Technology			
RD SET Adhesive Disposable Sensors			
690-0283-00	Masimo RD SET Adt, adhesive SpO ₂ sensor	Adult, > 30 kg finger or toe	20/Box
690-0285-00	Masimo RD SET Inf, adhesive SpO ₂ sensor	Infant, 3—20 kg, thumb or great toe	20/Box
690-0286-00	Masimo RD SET Neo, adhesive SpO ₂ sensor	Neonatal <3 kg, hand or foot Adult > 40 kg finger or toe	20/Box


Nellcor® SpO₂ Supplies & Accessories Nellcor™ OxiMax® Technology

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor. Cables and connectors are color-coded purple. **Note, use only genuine Nellcor sensors.**



Order Code	Description	Length/Size	Qty
Nellcor OxiMax Technology			
Cable			
	700-0792-00 TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
Reusable Sensors			
Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien), now part of Medtronic.			
	690-0003-01 Nellcor DS-100A, finger SpO ₂ sensor	Adult	1 Each
	690-0004-01 Nellcor OXI-A/N, SpO ₂ sensor	Adult/Neonatal	1 Each
	690-0039-01 Nellcor OXI-P/I, SpO ₂ sensor	Pediatric/Infant	1 Each
	690-0215-01 Nellcor DS-Y, multi-site SpO ₂ sensor	Universal	1 Each

Temperature











TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Monitor Adapter Cables			
Adapters			
 700-0031-00	TruLink Dual Temperature Cable. *Required for all reusable probes and as directed for single-patient-use disposable system	30 cm/12 in.	1 Each
Reusable Probes/Sensors			
Skin Surface			
 690-0295-00	YSI 400 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0296-00	YSI 700 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0297-00	YSI 700 series compatible skin surface sensor temperature probe, pediatric size	305 cm/10 ft	1 Each
Esophageal / Rectal			
 690-0294-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
 690-0298-00	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
 690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
 118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each













Indicates discontinued part, available while supplies last

TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Single-Patient use temperature system			
Extension cables for single use probes			
 700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors, adapter required.	305 cm/10 ft	1 Each
 700-0401-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors. Direct connection to Command module, no adapter required.	305 cm/10 ft	1 Each
Skin			
 701-4000-00	Skin Temperature Sensor		20/Box
Esophageal/ Stethoscope			
 701-4009-00	9 French, YSI 400 series compatible		20/Box
 701-4012-00	12 French, YSI 400 series compatible		20/Box
 701-4018-00	18 French, YSI 400 series compatible		20/Box
Esophageal/ Rectal			
 701-4109-00	9 French, YSI 400 series compatible		20/Box
 701-4112-00	12 French, YSI 400 series compatible		20/Box
Foley catheter with temperature probe			
 701-4214-00	14 French, YSI 400 series compatible		20/Box
 701-4216-00	16 French, YSI 400 series compatible		20/Box
Tympanic temperature sensor			
 701-4300-00	Adult, YSI 400 series compatible		20/Box
Exergen			
<i>For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) & DM4 (91331) Monitors</i>			
 010-2157-00	Temporal Artery Thermometer		1 Each

TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Reusable Temperature Probes			
700-0031-00	Dual temp cable		
690-0294-00	YSI400 Adult Esophageal/Rectal		
690-0298-00	YSI700 Adult Esophageal/Rectal		
690-0295-00	YSI400 Adult Skin		
690-0296-00	YSI700 Adult Skin		
690-0297-00	YSI700 Pediatric Skin		
Single-Patient Use Temperature Probes			
700-0401-00	Reusable extension cable		
700-0031-00	Dual temp cable		
700-4000-00	Reusable extension cable		
701-4009-00	w/ Stethoscope 9Fr		
701-4012-00	12Fr		
701-4018-00	18Fr		
701-4109-00	Esophageal/Rectal 9Fr		
701-4112-00	12Fr		
701-4214-00	w/ Foley Catheter 14Fr		
701-4216-00	16Fr		
701-4000-00	Skin		
701-4300-00	Tympanic		







Invasive Blood Pressure



Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.




Order Code	Description	Length/Size	Qty
Integrated Adapter Cables (connects directly to transducer)			
Edward Lifesciences			
 015-0727-00	Single pressure integrated cable	400 cm/13 ft	1 Each
015-0728-00	Dual pressure integrated cable	400 cm/13 ft	1 Each
Abbott Laboratories/Hospira/ICU Medical			
 015-0729-00	Single pressure integrated cable	400 cm/13 ft	1 Each
700-0294-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
015-0730-00	Dual pressure integrated cable	400 cm/13 ft	1 Each
Bold font indicates a new part number			
Indicates discontinued part, available while supplies last			
Interface Cable & Adapters			
Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)			
 690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
Interface Adapter- Transducer interface cable required from transducer company			
 700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
Transducer Kits & Accessories			
Pressure Transducer Kits, Single use			
 690-0008-00	Transducer Kit, Patient Mount, DTX/Plus	152.5 cm/60 in	5/Box
 690-0009-00	Transducer Kit, Pole Mount, DTX/Plus	152.5 cm/60 in	5/Box

Invasive Pressure Supplies


Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



Cardiac Output



CARDIAC OUTPUT




Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.






Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Cardiac Output Cables & Probes

Cables			
	700-0027-00	TruLink Cardiac Output Adapter Cable. Connects CO cable 30 cm/12 in. to 91496 command module.	1 Each
	306655-003	Cardiac Output cable compatible connects with 690-0033-00 temperature probe. 700-0027-00 adapter required.	1 Each
	306655-004	Cardiac Output cable compatible connects with 690-0023-00 or 690-0024-00 temperature probe. 700-0027-00 adapter required.	1 Each


Indicates discontinued part, available while supplies last. Spacelabs is releasing new parts to replace the 306655-003 and 306655-004. See July 2024 catalog when available.

Injectate Probes

	690-0023-00	In-line injectate temperature probe. Connects to 306655-004 cable and 690-0293-00 injectate kit integrated thermistor.	1 Each
	690-0024-00	Ice bath temperature probe. Connects to 306655-004 cable and needle end is submerged in an ice bath.	1 Each
	690-0033-00	In-line injectate temperature probe. Connects to 306655-003 cable and 690-0293-00 injectate kit integrated thermistor.	1 Each

Cardiac Output Injectate Kits

Closed Loop Injectate Kits

690-0293-00		SP5045 Critikit In-Line Sensor, Control Syringe closed loop injectate kit, sterile	5/Box
-------------	---	--	-------

CARDIAC OUTPUT


Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

91496-C Module, Cardiac Output Connection guide



Recording & Printing



Paper

High quality Z-fold paper for use with the Bedside Printer Module (90449) and Central Printer (90469).

Roll paper for use with the Standalone printer (91449) and Spacelabs Portable Monitors.



(91449)
Printer

Order Code	Description	Length/Size	Qty
Thermal Paper			
Z-fold Paper for use with 90449 and 90469			
006-0196-00	TruLink Thermal Paper, Z-Fold Pads, 2-channel printing, blank paper, 230 Sheets/pad	27.4 m/90 ft	10 Pads/Box
Roll Paper for use with 91449, 91367, 91369, 91370 & 91390			
307438-001	TruLink Thermal Paper, Rolls, single- or dual-channel printing, 50 mm/1.97 inch wide blank paper	30.1 m/100 ft	10 Rolls/Box
Recording Paper			
20257-004	Thermal recorder paper, black	HearTwave II	each
010-2246-00	PAPR,THERM,Z-FOLD, 210X295MM X 180 SHTS/10RM,IEC	Cardioline	each

Bispectral Index



Bispectral Index Supplies

For use with Spacelabs BISx™ Module (91482).



Order Code	Description	Qty
BISx Cables & Sensors		
BISx Cables		
719-0004-01	BISx Patient Interface cable	1 Each
719-0006-20	BISx Pod	1 Each
719-0007-01	BISx host monitor cable	1 Each


BISx Sensors

Note: Sales of BISx sensors are limited to one box with the purchase of capital equipment. All future purchases should be made directly through Covidien, now part of Medtronic.

719-0003-10	BISx Extend Sensor, single patient use	25/Box
719-0002-10	BISx Quatro Sensor, single patient use	25/Box

Ultraview DM3

DUAL MODE VITAL SIGNS MONITOR



VITAL SIGNS MONITOR

Ultraview DM3

Vital Sign Monitoring

Ultraview DM3 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

See previous sections for NIBP and SpO2 supplies:

- The DM3 uses the TruLink NIBP cuffs (single tube, bayonet fitting) and the standard air hose 714-0018-00.
- The DM3 used the same cables and sensors, by brand, for pulse oximetry.



DM3 91330

Order Code	Description	Size	Length/Size	Qty
Supplies unique to the DM3 monitor:				
019-0539-00	Replacement Temperature			1 Each
Indicates discontinued part, available while supplies last				

DM4

VITAL SIGNS MONITOR



VITAL SIGNS MONITOR

Ultraview DM4









Vital Sign Monitoring

Ultraview DM4 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

The DM4 uses TruLink branded NIBP cuffs (single tube, bayonet fitting) and air hose (714-0018-00)




91331

Order Code	Description	Length/Size	Qty
Supplies unique to the DM4 monitor			
LNCS Adapter Cable			
 010-2227-00	DM4 Masimo LNCS SpO ₂ cable compatible with Masimo LNCS sensors, see section on Masimo LNCS sensors	305 cm/10 ft	1 Each
RD SET Adapter Cable			
 010-2193-00	DM4 Masimo RD SET SpO ₂ cable compatible with Masimo RD SET sensors, see section on RD SET sensors	305 cm/10 ft	1 Each
Nellcor OxiMax Cable			
 010-2192-00	DM4 Nellcor OxiMax SpO ₂ cable compatible with Nellcor OxiMax sensors, see section on Nellcor sensors	305 cm/10 ft	1 Each
Filac			
 010-2174-00	Filac dock and insulation chamber		1 Each
 010-2175-00	Filac temperature probe		1 Each
 010-2228-00	Filac temperature oral probe covers, 20/bx		20/box 25/boxes
Exergen			
 010-2157-00	Temporal Artery Thermometer		1 Each
NeoNatal NIBP supplies			
 714-0087-00	NIBP,HOSE ASSY,NEONATE,DM4 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	1 Each

Elance

VITAL SIGNS MONITOR



Elace Vital Signs Monitor







The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add Microstream end tidal CO2.

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below
- TruLink ECG cables and leadwire systems are compatible




Elace Vital Signs Monitor
93300

Order Code	Description	Qty
Supplies unique to the Elance monitor		
ECG Lead Wires & Cables (can also use standard TruLink ECG cables and leadwire systems)		
ECG Integrated Cables and Lead Wires		
	700-0008-22 5-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
	700-0008-23 5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
	700-0008-24 3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)		
SpO2 Cable		
	700-0287-00 SpO2 adapter cable	305 cm/10 ft 1 Each
LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)		
LNCS Adapter Cable		
	690-0256-00 SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft 1 Each
Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)		
Cable		
	700-0290-00 SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft 1 Each
Temperature		
Temperature Probe		
	690-0261-01 Skin/Surface Temp probe, reusable, YSI-400	305cm/10ft 1 Each
Invasive Blood Pressure (See IBP section for compatible sensors and transducers)		
Interface Cable		
	690-0021-00 Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	366 cm/12 ft 1 Each
Indicates discontinued part, available while supplies last		

C50

VITAL SIGNS MONITOR













Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
ECG leadwire systems				
 040-000908-00	CABLE,INTEGRATED 3-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000909-00	CABLE,INTEGRATED 5-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000910-00	CABLE,INTEGRATED 10-WIRE CLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000911-00	CABLE,INTEGRATED 3-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000912-00	CABLE,INTEGRATED 5-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000913-00	CABLE,INTEGRATED 10-WIRE CLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
Pulse Oximetry Supplies				
 040-000869-00	CABLE,INTEGRATED ADT FINGER CLIP SENSOR,C50	3 meters	Reusable	1 Each
 040-000769-00	CABLE,SPO2,12PIN TO DB9F,C50 <i>Sensors sold seperately</i>	3 meters	Reusable	1 Each
 040-000726-00	SENSOR,SPO2,ADT,FINGER CLIP,C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each
 040-000730-00	SENSOR,SPO2,ADT/NEO/PED,WRAP, 3',C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each









Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Non-Invasive Blood Pressure supplies				
 040-000808-00	ADULT NIBP TUBE (WHITE)2M	2 meters	Reusable	1 Each
	Reusable or Disposable Cuffs See NIBP section under TruLink Reusable or Disposable cuffs. These are compatible with the C50. <i>Requires use of the 040-00808-00 air hose (above).</i>			
 714-0086-00	NIBP,HOSE ASSY,NEONATE,C50 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	Reusable	1 Each
Capnography Supplies: requires external Phasein ISA module, 099-000007-00				
 040-000400-00	CABLE,CO2/AG,C50 Connects the Phasein ISA module to the C50	1 meter	Reusable	1 Each
	For airway adapter, cannula and sampling line sets, please see the capnography supply section under Nomoline SPU supplies. These are compatible with the C50.			
Temperature sensors				
 040-000246-00	PROBE,TEMP,ADT,SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000387-00	PROBE,TEMP,NEO SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000385-00	PROBE,TEMP,ADT,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each
 040-000386-00	PROBE,TEMP,NEO,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each








Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Invasive Blood Pressure Cables and Transducer kits				
	040-000454-00 INTERFACE CABLE,IBP,EDWARD CONNCTR,C50 <i>Requires 040-000211-00 transducers</i>	 3 meters	Reusable	1 Each
	040-000360-00 INTERFACE CABLE,IBP,ABBOTT CONNCTR,C50 <i>Requires 040-000870-00 transducers</i>	 3 meters	Reusable	1 Each
	040-000456-00 INTERFACE CABLE,IBP,BD CONNCTR,C50 <i>Requires 040-000403-00 transducers</i>	 3 meters	Reusable	1 Each
Recording Paper				
	040-000116-00 RECORDER PAPER,THERMAL,C50	50mm X 20 meters		1 roll Each

Anesthesia



Anesthesia

Anesthesia System For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseSirius Anesthesia Machines.



BleaseFocus



BleaseSirius



ARKON

Order Code	Description	Qty
------------	-------------	-----

Anesthesia Supplies

Patient Flow Sensors

For use with BleaseSirius and BleaseFocus Anesthesia Machines only



* 10110090	Adult Patient flow sensor	1 Each
717-0428-00	D-Lite reusable air flow sensor, may require luer caps	1 Each
050-0842-00	Luer caps for D-Lite air flow sensor	Pack of 5
<i>See image for assembly, available for sale to EU</i>		
10110218	2M Patient flow sensor tubing assembly	1 Each
10110102	3M Patient flow sensor tubing assembly	1 Each
10110103	6M Patient flow sensor tubing assembly	1 Each
70300001	O2 probe with tee	1 Each
376-0561-01	Pediatric flow sensor kit *EU/UK ONLY*	1 Each
FR101029	Pediatric flow sensor and tubing kit, 2.0 meters	1 Each
FR101030	Pediatric flow sensor and tubing kit, 3.0 meters	1 Each

Indicates discontinued part, available while supplies last

* Product not available for sale into CE Mark Countries (EU)

Ambulatory Blood Pressure Monitor



Ambulatory Blood Pressure

Ambulatory blood pressure is a leading diagnostic tool. Rather than relying on one isolated reading, it compiles a series of measurements resulting in increased accuracy. Spacelabs offers a line of reusable and disposable cuffs, covers and pouches for use with the Spacelabs 90227, 90207 and 90217A APB Monitors.



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Reusable ABP Cuffs



UltraCheck Reusable Nylon ABP Cuffs

016-1320-02	Statcorp reusable ABP cuff	child	13 - 18 cm	1 Each
016-1826-02	Statcorp reusable ABP cuff	small adult	18 - 26 cm	1 Each
016-2635-02	Statcorp reusable ABP cuff	adult	26 - 35 cm	1 Each
016-3242-02	Statcorp reusable ABP cuff	large adult	32 - 42 cm	1 Each



Spacelabs Reusable Cloth ABP Cuffs (includes bladder, wrap, hose, fittings)

015-0118-05Q	Spacelabs reusable ABP cuff kit	child	12 - 20 cm	1 Each
015-0067-05Q	Spacelabs reusable ABP cuff kit	small adult	17 - 26 cm	1 Each
015-0068-06Q	Spacelabs reusable ABP cuff kit	adult	24 - 32 cm	1 Each
016-0077-06Q	Spacelabs reusable ABP cuff kit	large adult	32 - 42 cm	1 Each
016-0109-05Q	Spacelabs reusable ABP cuff kit	XL adult	38 - 50 cm	1 Each



Disposable Covers



Disposable Cuff Covers (for use with Spacelabs reusable ABP cuffs above)

706-0182-00	Universal cuff cover	Standard	50/Box
706-0183-00	Universal cuff cover	Large	50/Box

Disposable ABP Cuffs



Spacelabs Soft Fabric Disposable ABP Cuffs (no fitting)






015-1826-00	Spacelabs disposable ABP cuff	small adult	18 - 26 cm	20/Box
015-2635-00	Spacelabs disposable ABP cuff	adult	26 - 35 cm	20/Box
015-3242-00	Spacelabs disposable ABP cuff	large adult	32 - 42 cm	20/Box
015-4250-00	Spacelabs disposable ABP cuff	XL adult	42 - 50 cm	20/Box
ZFT-QC8220	Quick release metal fittings, female bayonet (HP), 1/8 in. barb, For use with above single tube, disposable cuffs with no fitting.			each

AMBULATORY BLOOD PRESSURE

Ambulatory Blood Pressure

Ambulatory blood pressure is a leading diagnostic tool. Rather than relying on one isolated reading, it compiles a series of measurements resulting in increased accuracy. Spacelabs offers a line of reusable and disposable cuffs, covers and pouches for use with the Spacelabs 90227, 90207 and 90217A APB Monitors.



	Order Code	Description	Size	Range	Qty
ABP Accessories					
	Pouches - For use with ABP Monitor 90207				
	016-0410-01	Black polypropylene reusable pouch			1 Each
	Pouches - For use with ABP Monitor 90217				
	016-0340-01	Black polypropylene reusable pouch			1 Each
	Pouches - For use with ABP Monitor OnTrak 90227				
	015-0715-01	Black polypropylene reusable pouch			1 Each
	Support				
	015-0070-00	Cuff support harness			1 Each
	PK0001	1 inch with slide release buckle nylon webbing adjustable to max 59" inch length. For use as waist belt or shoulder strap.			1 Each
Miscellaneous ABP Accessories					
	712-0773-00	Adapter, converts ABP Luer Monitors (90207 and 90217) to quick disconnect			1 Each
	712-0794-00	Quick Disconnect metal hose connector for ABP Monitors 90207 & 90217			1 Each
	712-0060-00	QUICK CONNECTOR,1/8 TUBE,NIBP,CUFF			1 Each

Indicates discontinued part, available while supplies last

Holter and Event Monitoring



Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.











Lifecard CF

Evo

Aria

CardioCall

Order Code	Description	For Use With	Qty
Aria/evo Hook-Up Kits			
 23189-106-401	4 electrode hook-up kit Includes: 4 clear tape electrodes, 1 AA battery, 4 prep pads, 1 prep scrub, 1 Aria disposable pouch, 6 tape strips, 1 prep razor and 1 patient diary	Aria/evo	20/Case
 50-0049	4 electrode hook-up kit with disposable evo chest pouch Includes: 4 clear tape electrodes, 4 prep pads, 1 prep scrub, 1 evo disposable pouch, 6 tape strips, 1 prep razor, 1 patient diary and pencil	Aria/evo	Each
23-0255	Lanyard clip & assembly	CardioCall	Each
Lead Wires			
 25452-525	4-Lead wire set, snap, 3 - 14" leads (orange/brown/red) and 1 - 7" lead (white)	Aria/evo	4LW/set
 25452-535	4-Lead wire set, snap, 14" leads (red/white/orange/brown)	Aria/evo	4LW/set
 25452-529	4-Lead wire set, snap, 40" leads (red/white/orange/brown)	Aria/evo	4LW/set
 25452-519	7-Lead wire set, snap, 18" leads (white/black/red/green/orange/blue/brown)	423, 423A, 459, 461, 461A, 464, 483, 485	7 LW/set
700-0366-00	White, snap, 14" lead wires	Aria/evo	10 /pack
19-7509	6-Lead wire set, snap, 24 in, AHA (requires 700-0316-00)	Lifecard CF	6LW/Set
23-0059	Patient cable, 1 channel, 2 electrodes (black/white), US	CardioCall	Each
23-0051	Patient cable, 2 channel, 3 electrodes (red/yellow/green)	CardioCall	Each
Electrodes and Prep supplies			
 392196-001	Adhesive tape remover pads, 100/box		100/Box
 717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches		50/Box
717-0071-00	ECG electrode, Ambu Blue Sensor L		5/Pack

Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF

Evo

Aria

CardioCall

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Lifecard CF Patient Cables



700-0316-00	3-channel, 6-lead, 22 in, AHA (requires 19-7509)	Lifecard CF	Each
700-0317-00	3-channel, 3-lead, snap, long	Lifecard CF	Each
700-0318-00	3-channel, 3-lead, snap, extra-long	Lifecard CF	Each
700-0319-00	3-channel, 3-lead, snap, short	Lifecard CF	Each
700-0320-00	CABLE,PX,LCF,2SNAP,CLIP,RD/WH.ROHS	Lifecard CF	Each
700-0321-00	CABLE,PX,LCF,2SNAP,SHORT,WO/CLIP,RD/WH.ROHS	Lifecard CF	Each
700-0322-00	CABLE,PX,LCF,7RECEPTACLE,AHA,WO/LEADWIRES,ROHS	Lifecard CF	Each
700-0323-00	2-channel, 4-lead, snap, long	Lifecard CF	Each
700-0324-00	2-channel, 4-lead, snap, short	Lifecard CF	Each

Lifecard 12 Patient Cables



700-0328-00	12-lead patient cable, EU	Lifecard 12	Each
700-0329-00	12-lead patient cable, US	Lifecard 12	Each
700-0330-00	12-lead patient cable, long, EU	Lifecard 12	Each
700-0331-00	12-lead patient cable, long, US	Lifecard 12	Each

Pouches, Shoulder Straps & Belts



27249-7-401	Disposable pouch	Aria	50/Pack
50-0191	HOOK UP KIT WITH SINGLEUSE POUCH, CASE OF 20	Evo	20/Case
50-0042	Belt holster	Evo	Each
PK0001	Adjustable 1" nylon webbing w/ buckle, max 56" length	Evo	Each
46-0005	Disposable cord (use with 706-0189-00)	Evo	50m roll
50-0050	Disposable pouch	Evo	50/Pack
46-0408	Fabric pouch set, 2 pouches with laces, Cordlock	Lifecard CF	Each
19-7130	Replacement laces for 46-0408	Lifecard CF	Each
19-7131	Replacement cord lock for 46-0408	Lifecard CF	Each
46-0092	Replacement pouch for 46-0408	Lifecard CF	Each
46-0512	Disposable pouch	Lifecard CF	100/Pack



Indicates discontinued part, available while supplies last

Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF

Evo

Aria

CardioCall



Order Code	Description	For Use With	Qty
Recorder Memory and Download			
CFC90MB	Compact flash card memory, 90MByte	Lifecard CF	Each
CFC256MB	Compact flash card memory, 256MB	LC12	Each
46-0063	Flash card plastic protective case	CFCxxMB	Each
19-7118	Compact USB Flashcard Reader	CFCxxMB	Each
<i>Note: Spacelabs is working to release an alternate card reader due to manufacturer design change</i>			
Patient Diaries			
073-0267-00	7 Patient Diary, English CE0123	Aria, Evo, Lifecard CF	Each

Eclipse Mini and Eclipse Pro



Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.

Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code	Description	Length/Style	Qty																						
Eclipse Pro																									
Leadwire Sets																									
700-0378-00	3 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	205	8	510	20	865	34	Each
Snap	Color	Short			Medium		Long																		
		mm	Inch	mm	Inch	mm	Inch																		
C5	Red	205	8	510	20	865	34																		
700-0391-00	3 Lead snap on with yoke Medium	<table border="1"> <tbody> <tr> <td>M</td> <td>Yellow</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> </tbody> </table>	M	Yellow	255	10	560	22	915	36	Each														
M	Yellow	255	10	560	22	915	36																		
700-0390-00	3 Lead snap on with yoke Long	<table border="1"> <tbody> <tr> <td>C5R</td> <td>Green</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	C5R	Green	205	8	510	20	865	34	Each														
C5R	Green	205	8	510	20	865	34																		
700-0379-00	4 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	230	9	510	20	865	34	Each
Snap	Color	Short			Medium		Long																		
		mm	Inch	mm	Inch	mm	Inch																		
C5	Red	230	9	510	20	865	34																		
700-0394-00	4 Lead snap on with yoke Medium	<table border="1"> <tbody> <tr> <td>C3</td> <td>Orange</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> </tbody> </table>	C3	Orange	255	10	560	22	915	36	Each														
C3	Orange	255	10	560	22	915	36																		
700-0393-00	4 Lead snap on with yoke Long	<table border="1"> <tbody> <tr> <td>M</td> <td>White</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C1</td> <td>Brown</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	M	White	255	10	560	22	915	36	C1	Brown	230	9	510	20	865	34	Each						
M	White	255	10	560	22	915	36																		
C1	Brown	230	9	510	20	865	34																		
700-0396-00	10 Lead snap on with yoke <i>IEC color coding generally used in European Union</i>	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th colspan="2">Length</th> </tr> <tr> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>L</td> <td>610</td> <td>24</td> </tr> <tr> <td>F</td> <td>585</td> <td>23</td> </tr> <tr> <td>C6</td> <td>635</td> <td>25</td> </tr> <tr> <td>C5</td> <td>535</td> <td>21</td> </tr> </tbody> </table>	Snap	Length		mm	Inch	L	610	24	F	585	23	C6	635	25	C5	535	21	Each					
Snap	Length																								
	mm	Inch																							
L	610	24																							
F	585	23																							
C6	635	25																							
C5	535	21																							
700-0377-00	10 Lead snap on with yoke <i>AHA color coding generally used in USA, outside of EU</i>	<table border="1"> <tbody> <tr> <td>C4</td> <td>460</td> <td>18</td> </tr> <tr> <td>C3</td> <td>355</td> <td>14</td> </tr> <tr> <td>C2</td> <td>330</td> <td>13</td> </tr> <tr> <td>C1</td> <td>330</td> <td>13</td> </tr> <tr> <td>N</td> <td>380</td> <td>15</td> </tr> <tr> <td>R</td> <td>610</td> <td>24</td> </tr> </tbody> </table>	C4	460	18	C3	355	14	C2	330	13	C1	330	13	N	380	15	R	610	24	Each				
C4	460	18																							
C3	355	14																							
C2	330	13																							
C1	330	13																							
N	380	15																							
R	610	24																							
712-1052-01	ECLIPSE PRO ECG SENSORS <i>Not available for sale in USA</i>		Box of 10																						
Accessories																									
706-0207-00	Eclipse Pro belt clip		Each																						
073-0267-00	7 Patient Diary, English CE0123		Each																						
392196-001	Adhesive tape remover pads, 100/box		100/Box																						

Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.

Eclipse Pro and Mini Holter




Eclipse Pro



Eclipse Mini

	Order Code	Description	Length/Style	Qty
Eclipse Mini				
	Eclipse Mini Supplies			
	712-1053-01	ECLIPSE MINI ECG SENSORS <i>Not available for sale in USA</i>		Box of 10
	712-1053-02	ECLIPSE MINI ECG SENSORS, DEFIB SAFE Now available for sale in USA		Box of 10

ECG Stress Testing



Brand Name Stress Test Devices

ECG Stress Testing

Spacelabs offers supplies for the Suntech Tango M2 ABP system. Spacelabs has discontinued sale of these systems but continues to offer some accessories, while supplies last. Note: not available for sale into USA but is available for sale in the European Union and United Kingdom.



Suntech Tango M2

Order Code	Description	Dimensions	Qty
------------	-------------	------------	-----



Suntech Tango NIBP Supplies for Stress Test

010-2238-00	CUFF, ORBIT-K, SMALL ADULT		Each
98-0061-02	Orbit-K cuff, adult, 23-35 cm		Each
98-0061-05	Orbit-K cuff, adult plus, 27-40 cm		Each
98-0061-03	Orbit-K cuff, large adult, 32-44 cm		Each
010-2237-00	K-SOUND MICROPHONE,14" FOR SUNTECH TANGO		Each
010-2241-00	K-SOUND MICROPHONE,12" FOR SUNTECH TANGO ^{Tá}		Each
98-0002-50	Microphone pads		5/Pack


Indicates discontinued part, available while supplies last

CardioPulse ECG



CardioPulse Resting ECG












Compatible supplies and accessories for the portfolio of CardioPulse ECG 12-lead systems. Note: CardioPulse is only available for sale in the USA.



CardioPulse Prime



CardioPulse Go

	Order Code	Description	For Use With	Qty
Paper				
	010-2280-00	PAPER,THERM,Z-FOLD,210 MM,100SHX10PKS	Prime	1000 sheets
	010-2281-00	PAPER,THERMAL,Z-FOLD,114 MM,60SHX10PKS	Go	600 sheets
Lead Wires, Cables, Adapters				
	700-0386-00	CABLE,PATIENT,AHA,SNAP TYPE,CARDIOPULSE GO	Go	Each
	700-0387-00	CABLE,PATIENT,10-WIRE,AHA,BANANA PLUG,CP GO	Go	Each
	700-0388-00	CABLE,PATIENT,10-WIRE,AHA,CLIP,CP PRIME,1/BG	Prime	Each
	700-0389-00	CABLE,PATIENT,10-WIRE,AHA,BANANAPLUG,CP PRIME	Prime	Each
	700-0399-00	ADAPTER,D15,PATIENT CABLE COMPATIBLE,CP GO	Go	Each
Tab Electrodes, for use with electrode clips				
	700-0400-00	ELECTRODES,WHITE,DISPOSABLE,500 PCS,1/BOX		Box of 500
	SU00C	Disposable tab electrodes for resting ECG set		Box of 100
	58100	Easiclip, Skintact, tab electrode clips <i>For use with banana 4mm sized plugs</i>		Pack of 10
Barcode Scanner				
	010-2282-00	BARCODE SCANNER,LASER,INSTRUCTIONS,CARDIOPULSE		Each

Indicates discontinued part, available while supplies last

Resting ECG



CardioExpress resting ECG

CardioDirect resting ECG

CardioExpress SL6A, SL12A and SL18A provide fast and accurate 12-lead ECG acquisition. Advanced design and ergonomics enable the operator to quickly and easily control the unit to meet individual requirements. They provide storage direct to the units and connectivity via Sentinel, Cardiology Information Management System (CIMS).



CardioExpress



CardioDirect 12

Order Code	Description	For Use With	Qty
ECG Patient Cables & Hook-up Sets			
MS1R-107584-AO	ECG cable, grabber connector, AHA	CardioExpress	Each
MS1R-107048-AO	ECG cable, 4MM, banana connector, AHA	CardioExpress	Each
010-2106-00	ECG cable, banana plugs, AHA	CardioExpress SL18A	Each
010-2107-00	ECG CABLE, SL18A, IEC, BANANA PLUGS	CardioExpress SL18A	Each
MS1R-106902-AO	ECG cable, 4MM, banana connector, IEC	CardioExpress	Each
25-0004	10 electrode patient cable, 4mm, EU	CD12	Each
25-0005	10 electrode patient cable, pinch, EU	CD12	Each
M15R-040172	CLIP/SNAP, BANANA SOCKET ADAPTERS, 10 PCS	CardioExpress	Each
MS1R-107449-AO	SNAP/BANANA SOCKET ADAPTERS (10 PCS/SET)	CardioExpress	Each
M15R-040163	ADULT PRECORDIAL SUCTION ELECTRODES (6 PCS/SET)	CardioExpress	Each
M15R-040168	CHILD PRECORDIAL SUCTION ELECTRODES (6PCS/SET)	CardioExpress	Each
012-0907-00	ELECTRODE,ADULT LIMB CLAMP, QTY 4	CardioExpress	Each
58100	Easiclip, Skintact, 10/pk		10/Pack
AHA CardioExpress SL18A lead wires			
012-0963-00	RA ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0964-00	LA ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0965-00	V5R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0966-00	V4R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0967-00	V3R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0968-00	V1+V2+V3 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0969-00	V4+V5+V6 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0970-00	V7 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0972-00	V9 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0973-00	LL ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0974-00	RL ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
ECG Patient Cables & Hook-up Sets			
IEC CardioExpress SL18A lead wires			
012-0975-00	R (RA) ECG Lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0976-00	L (LA) ECG LEADWIRE FORSL-18A SAMPLING BOX IEC	CardioExpress SL18A	Each
012-0982-00	C1+C2+C3 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0983-00	C4+C5+C6 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0985-00	F (LL) ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0986-00	N (RL) ECG LEADWIRE FORSL-18A SAMPLING BOX IEC	CardioExpress SL18A	Each

Indicates discontinued part, available while supplies last

CardioExpress resting ECG

CardioDirect resting ECG

CardioExpress SL6A, SL12A and SL18A provide fast and accurate 12-lead ECG acquisition. Advanced design and ergonomics enable the operator to quickly and easily control the unit to meet individual requirements. They provide storage direct to the units and connectivity via Sentinel, Cardiology Information Management System (CIMS).

Resting ECG systems 6, 12 and 18 leads



CardioExpress



CardioDirect 12

Order Code	Description	For Use With	Qty
Miscellaneous ECG Accessories			
Recording Paper			
MS1R-107371-A0	Z-fold, 210 MM x 295 MM, 100 sheets	CardioExpress	Each
MS1R-32462	Z-fold, 210 MM x 295 MM, 200 sheets	CardioExpress	Each
MS1R-107560-A0	Z-fold, 110 MM x 140 MM, 144 sheets	CardioExpress	Each
MS1R-107451-A0	Z-fold, 215 MM x 280MM, 100 sheets	CardioExpress	Each
MS1R-32461	Roll, 210 MM x 30M	CardioExpress	Each
010-1897-00	RECORDING PAPER,ROLL,80MM X 20 M	CardioExpress	Each
Others			
010-2127-00	Carry case, SL6A	CardioExpress SL6A	Each
010-2128-00	Carry case, SL12A	CardioExpress SL12A	Each
690-0275-00	Barcode scanner, one dimensional	CardioExpress	Each
690-0276-00	Barcode scanner, two dimensional	CardioExpress	Each
SU00C	Disposable tab electrodes for resting ECG set	CD12USB	Each

Indicates discontinued part, available while supplies last

Non-Medical



Non-Medical Accessories

Spacelabs carries a convenient assortment of Non-Medical supplies to help support your use of our monitoring systems.

Order Code Description With Qty



Note: Bumper is black. Yellow color used in image to highlight design and fit of the bumper.

Bumper Kits

010-2350-00	Aria Tele Bumpers	Aria Tele	5/Box
050-0816-00	Qube Mini Bumper Kit (black)	Qube Mini	Each



Cable Management Systems

016-1027-00	CORD CADDY, 4 PAIRS OF HOOKS	\$224.89	Each
016-1028-00	CORD CADDY, 6 PAIRS OF HOOKS	\$265.78	Each
016-1029-00	MOUNT,DOWNPOST,GCX,12 INCH (required to mount caddy)	\$65.52	Each

AV and Networking Cables

2287169	TRIPP LITE 10' HDMI M/M CABLE	\$19.14	Each
2195215	DP mini to DP cable, male/male, 10ft	\$26.14	Each
1973897	DP mini to DP cable, male/male, 6ft	\$20.22	Each
384343	USB 2.0 A to B Printer Cable, 10ft	\$11.84	Each
1561876	DP to DP cable, male/male, 10ft	\$28.82	Each
54403	DP to DP cable, male/male, 15ft	\$51.82	Each

Computer Supplies



5491141	NETGEAR GS308V3 8PRT SWITCH	\$43.10	Each
5615760	NETGEAR GS305 5PRT SWITCH	\$35.26	Each
3006282	LOGITECH MK270 WRLS KBD/MS SET	\$47.50	Each
010-2327-00	Poly Blackwire 5210 Monaural headset for use with SafeNSound	\$120.66	Each
6113257	USB-A port blockers (4) with reusable key (1)	\$56.46	Pack
6113260	USB-A port blockers (10)	\$66.80	Pack

Barcode Scanners



690-0275-00	Zebra LS2208 Barcode Scanner	\$250	Each
690-0276-00	Honeywell 1900GSR-2USB-C Barcode Scanner	\$400	Each

Non-Medical Accessories

Spacelabs carries a convenient assortment of Non-Medical supplies to help support your use of our monitoring systems.

Order Code	Description	With	Qty
------------	-------------	------	-----



Displays and Mounts

210326	ERGOTRON DS100, DUAL DISP MOUNT	\$668.66	Each
6358609	NEC 43" LED-BACKLIT LCD DISPLAY	\$1,575.74	Each
6358610	NEC 50" LED-BACKLIT LCD DISPLAY	\$2,064.32	Each



Exergen TAT and Mounts: for use with Xprezzon, Qube, Qube Mini, DM4

010-2157-00	Exergen TAT 5000S, with coiled USB cable		Each
134201	Exergen TAT 5000S Wall mount, black		Each
134305	Exergen TAT 5000S locking wall mount, white		Each
134306	Exergen TAT 5000S locking wall mount with cap holder, white		Each

Universal Power Supplies and Network Switches



317698	APC,CS500,UPS,300WATT,500VA	\$273.40	Each
6064874	APC,BE850G2,UPS,450WATT,850VA	\$253.08	Each
5491141	NETGEAR GS308V3 8PRT SWITCH	\$43.10	Each
5615760	NETGEAR GS305 5PRT SWITCH	\$35.26	Each