


C50

VITAL SIGNS MONITOR













Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
ECG leadwire systems				
	040-000908-00 CABLE,INTEGRATED 3-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000909-00 CABLE,INTEGRATED 5-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000910-00 CABLE,INTEGRATED 10-WIRE CLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000911-00 CABLE,INTEGRATED 3-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000912-00 CABLE,INTEGRATED 5-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000913-00 CABLE,INTEGRATED 10-WIRE CLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
Pulse Oximetry Supplies				
	040-000869-00 CABLE,INTEGRATED ADT FINGER CLIP SENSOR,C50	3 meters	Reusable	1 Each
	040-000769-00 CABLE,SPO2,12PIN TO DB9F,C50 <i>Sensors sold seperately</i>	3 meters	Reusable	1 Each
	040-000726-00 SENSOR,SPO2,ADT,FINGER CLIP,C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each
	040-000730-00 SENSOR,SPO2,ADT/NEO/PED,WRAP, 3',C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each









Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Non-Invasive Blood Pressure supplies				
 040-000808-00	ADULT NIBP TUBE (WHITE)2M	2 meters	Reusable	1 Each
	Reusable or Disposable Cuffs See NIBP section under TruLink Reusable or Disposable cuffs. These are compatible with the C50. <i>Requires use of the 040-00808-00 air hose (above).</i>			
 714-0086-00	NIBP,HOSE ASSY,NEONATE,C50 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	Reusable	1 Each
Capnography Supplies: requires external Phasein ISA module, 099-000007-00				
 040-000400-00	CABLE,CO2/AG,C50 Connects the Phasein ISA module to the C50	1 meter	Reusable	1 Each
	For airway adapter, cannula and sampling line sets, please see the capnography supply section under Nomoline SPU supplies. These are compatible with the C50.			
Temperature sensors				
 040-000246-00	PROBE,TEMP,ADT,SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000387-00	PROBE,TEMP,NEO SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000385-00	PROBE,TEMP,ADT,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each
 040-000386-00	PROBE,TEMP,NEO,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each








Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Invasive Blood Pressure Cables and Transducer kits				
	040-000454-00 INTERFACE CABLE,IBP,EDWARD CONNCTR,C50 <i>Requires 040-000211-00 transducers</i>	 3 meters	Reusable	1 Each
	040-000360-00 INTERFACE CABLE,IBP,ABBOTT CONNCTR,C50 <i>Requires 040-000870-00 transducers</i>	 3 meters	Reusable	1 Each
	040-000456-00 INTERFACE CABLE,IBP,BD CONNCTR,C50 <i>Requires 040-000403-00 transducers</i>	 3 meters	Reusable	1 Each
Recording Paper				
	040-000116-00 RECORDER PAPER,THERMAL,C50	50mm X 20 meters		1 roll Each