

Eclipse Mini and Eclipse Pro



Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.








Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code	Description	Length/Style	Qty																																						
Eclipse Pro																																									
Leadwire Sets																																									
700-0378-00	3 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>M</td> <td>Yellow</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C5R</td> <td>Green</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	205	8	510	20	865	34	M	Yellow	255	10	560	22	915	36	C5R	Green	205	8	510	20	865	34	Each
Snap	Color	Short			Medium		Long																																		
		mm	Inch	mm	Inch	mm	Inch																																		
C5	Red	205	8	510	20	865	34																																		
M	Yellow	255	10	560	22	915	36																																		
C5R	Green	205	8	510	20	865	34																																		
	700-0391-00	3 Lead snap on with yoke Medium	Each																																						
	700-0390-00	3 Lead snap on with yoke Long	Each																																						
	700-0379-00	4 Lead snap on with yoke Short	Each																																						
	700-0394-00	4 Lead snap on with yoke Medium	Each																																						
	700-0393-00	4 Lead snap on with yoke Long	Each																																						
	700-0396-00	10 Lead snap on with yoke <i>IEC color coding generally used in European Union</i>	Each																																						
	700-0377-00	10 Lead snap on with yoke <i>AHA color coding generally used in USA, outside of EU</i>	Each																																						
	712-1052-01	ECLIPSE PRO ECG SENSORS <i>Not available for sale in USA</i>	Box of 10																																						
Accessories																																									
	706-0207-00	Eclipse Pro belt clip	Each																																						
	073-0267-00	7 Patient Diary, English CE0123	Each																																						
	392196-001	Adhesive tape remover pads, 100/box	100/Box																																						

Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.

Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

	Order Code	Description	Length/Style	Qty
Eclipse Mini				
	Eclipse Mini Supplies			
	712-1053-01	ECLIPSE MINI ECG SENSORS <i>Not available for sale in USA</i>		Box of 10
	712-1053-02	ECLIPSE MINI ECG SENSORS, DEFIB SAFE Now available for sale in USA		Box of 10