

Elance

VITAL SIGNS MONITOR



Elace Vital Signs Monitor







The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add Microstream end tidal CO2.

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below
- TruLink ECG cables and leadwire systems are compatible



Elace Vital Signs Monitor
93300

Order Code	Description	Qty
Supplies unique to the Elance monitor		
ECG Lead Wires & Cables (can also use standard TruLink ECG cables and leadwire systems)		
ECG Integrated Cables and Lead Wires		
	700-0008-22 5-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
	700-0008-23 5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
	700-0008-24 3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)		
SpO2 Cable		
	700-0287-00 SpO2 adapter cable	305 cm/10 ft 1 Each
LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)		
LNCS Adapter Cable		
	690-0256-00 SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft 1 Each
Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)		
Cable		
	700-0290-00 SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft 1 Each
Temperature		
Temperature Probe		
	690-0261-01 Skin/Surface Temp probe, reusable, YSI-400	305cm/10ft 1 Each
Invasive Blood Pressure (See IBP section for compatible sensors and transducers)		
Interface Cable		
	690-0021-00 Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	366 cm/12 ft 1 Each
Indicates discontinued part, available while supplies last		