

Capnography & Multigas



NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients



NomoLine LH Nasal CO2 Cannula		Shelf Life*		
704-3817-00	Adult	24 months	2 m/7 ft	25/Box
704-3819-00	Infant	24 months	2 m/7 ft	25/Box
704-3820-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4624-00	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3821-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box




NomoLine LH Nasal/Oral CO2 Cannula				
704-3824-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4628-00	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3825-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box



NomoLine LH Single Nasal Prong CO2 Cannula				
704-3826-00	Adult	24 months	2 m/7 ft	25/Box

NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine LH Airway Adapter Set		Shelf Life*		
704-3814-00	Adult/Pediatric	36 months	2 m/7 ft	25/Box
704-4620-00	Adult/Pediatric	36 months	4 m/13 ft	25/Box

NomoLine LH (Low Humidity) Procedural Oxygen Masks



NomoLine LH Nasal CO2 Cannula		Shelf Life*		
704-5213-00 ^{USA}	Adult, Medium O2 Flow Rate	24 months	3 m/10 ft	30/Box
704-5214-00 ^{USA}	Pediatric, Medium O2 Flow	24 months	3 m/10 ft	30/Box
704-5215-00 ^{USA}	Adult, High O2 Flow Rate	24 months	3 m/10 ft	30/Box
704-5216-00 ^{USA}	Pediatric, High O2 Flow Rate	24 months	3 m/10 ft	30/Box

* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA

NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients



NomoLine HH Nasal CO2 Cannula		Shelf Life*		
704-3830-00	Adult	24 months	2 m/7 ft	25/Box
704-3832-00	Infant	24 months	2 m/7 ft	25/Box
704-3833-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-3834-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box
704-4933-00 ^{USA}	Infant w/O2	24 months	4 m/13 ft	25/Box



NomoLine HH Nasal/Oral CO2 Cannula				
704-3837-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4968-00 ^{USA}	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3838-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box



NomoLine HH Single Nasal Prong CO2 Cannula				
704-3839-00	Adult	24 months	2 m/7 ft	25/Box

NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine HH Airway Adapter Set		Shelf Life*		
704-3827-00	Adult/Pediatric	36 months	2 m/7 ft	25/Box
704-3828-00	Adult/Pediatric	36 months	3 m/10 ft	25/Box
704-4989-00 ^{USA}	Adult/Pediatric	36 months	4 m/13 ft	25/Box
704-4367-00	Infant/Neonatal	36 months	2 m/7 ft	25/Box

* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Sidestream Monitoring Supplies

Please note: the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.



CO2 Nasal Cannula with luer		Shelf Life*		
* 704-0176-00	Adult	60 months	2 m/7 ft	25/Box
CO2 Nasal Cannula w/O2 and luer				
* 704-0181-00	Adult	60 months	2 m/7 ft	25/Box
CE 704-0181-01	Adult	36 months	2 m/7 ft	25/Box
* 704-0180-00	Pediatric	60 months	2 m/7 ft	25/Box
** 704-0179-00	Infant	60 months	2 m/7 ft	25/Box
** 704-0186-00	Adult	60 months	3 m/10 ft	25/Box
CO2 Oral/Nasal Cannula w/O2 and luer				
* 704-0183-00	Adult	60 months	2 m/7 ft	25/Box
CE 704-0183-01	Adult	36 months	2 m/7 ft	25/Box
** 704-0187-00	Adult	60 months	3 m/10 ft	25/Box

Indicates discontinued part, available while supplies last
 * - Shelf life is from date of manufacture, printed on product label as an expiration date

Sidestream Monitoring Supplies - Intubated patients



Airway adapters and lines		Shelf Life*		
** 704-0191-00	Adult Line w/ T-Adapter	60 months	2 m/7 ft	25/Box
** 704-0192-00	Adult Line w/ T-Adapter	60 months	3 m/10 ft	25/Box
* 166-7121-00	sampling line	60 months	3 m/10 ft	50/Box
* 704-0175-00	T-Adapter	60 months		25/Box

* Product not available for sale into CE Mark Countries (EU)
 CE product is available for sale into CE Mark Countries (EU)
 ** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients.
Featuring Masimo IRMA mainstream technology



Sensor & Airway Adapters - required for mainstream monitoring

010-1980-00	Masimo IRMA mainstream sensor, reusable		1 Each
704-0173-00	Masimo IRMA mainstream airway adapter, single use <i>36 month shelf life from date of manufacture</i>	Adult/Pediatric	25/Box
704-0174-00	Masimo IRMA mainstream airway adapter, single use <i>36 month shelf life from date of manufacture</i>	Infant	10/Box

CAPNOGRAPHY





Spacelabs Module 91517

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCO₂ capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream, Non-intubated, single patient use			
		CO₂ Nasal Cannula	
		Shelf Life[^]	
* 704-0010-00	Adult	36 months	10/Box
* 704-0012-00	Pediatric	36 months	10/Box
		CO₂ Nasal Cannula w/O₂ Delivery	
* 704-0014-00	Adult	36 months	10/Box
CE 704-0014-10	Adult, (for sale into EU/UK only)	36 months	10/Box
* 704-0016-00	Pediatric	36 months	10/Box
* 704-0015-00	Infant/Neonatal	36 months	10/Box
		CO₂ Oral/Nasal Cannula	
* 704-0017-00	Adult	36 months	10/Box
		CO₂ Oral/Nasal Cannula w/O₂ Delivery	
* 704-0020-00	Adult	36 months	10/Box
* 704-0022-00	Pediatric	36 months	10/Box
Sidestream, Intubated, single patient use			
Sidestream - short term up to 12 hours			
* 704-0023-00	Airway adapter + sampling line, ET tube > 4.0 mm	36 months	10/Box
CE 704-0023-10	Airway adapter + sampling line, ET tube > 4.0 mm	36 months	10/Box
* 704-0024-00	Airway adapter + sampling line, Pediatric/Infant sized	36 months	10/Box
* 704-0027-00	CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box
CE 704-0027-10	CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box

CAPNOGRAPHY Spacelabs Module 91517

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream - long term up to 120 hours		Shelf Life[^]	
* 704-0025-00	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	36 months	10/Box
CE 704-0025-10	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	36 months	10/Box
* 704-0026-00	Airway adapter + sampling line with Nafion, Infant	36 months	10/Box
* 704-0028-00	CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box
CE 704-0028-10	CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box

Mainstream, for Intubated Patients Only



Mainstream Airway Adapters, for use with CO₂ Sensor 704-0001-00			
* 704-0005-00	CO ₂ airway adapter, single use, Adult	36 months	24/Box
* 704-0002-00	CO ₂ airway adapter, reusable, Adult	36 months	1 Each

Indicates discontinued part, available while supplies last

[^] - Shelf life is from date of manufacture, printed on product label as an expiration date

* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

MULTIGAS ANALYZER Spacelabs Modules 90518 & 91518

Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).

These two legacy multi-gas modules are obsolete. However, Spacelabs continues to provide supplies that are available.



(90518)



(91518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Multigas Supplies & Accessories

Anesthetic Agent Samples Lines



015-0128-00	SAMPLE LINE,MALE TO MALE,LUER,SPU,7'X.04",25/BX		25/Box
015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line:0.048 inches ID tubing (90518/91518) <i>36 month shelf life from date of manufacture</i>	305 cm/10 ft.	50/Box
162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 0.048 in ID tubing (90518/91518) <i>36 month shelf life from date of manufacture</i>	305 cm/10 ft.	50/Box

Water Traps, Reusable



717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)		10/Pack
717-0448-00	D-FEND water trap cartridge, Green. Replace at least once every 24 hours or when occluded (91518)		10/Pack

Note: manufacturer (GE) has discontinued these parts. Available while supplies last.