








Invasive Blood Pressure



Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description	Length/Size	Qty
Integrated Adapter Cables (connects directly to transducer)			
Edward Lifesciences			
 015-0727-00	Single pressure integrated cable	400 cm/13 ft	1 Each
015-0728-00	Dual pressure integrated cable	400 cm/13 ft	1 Each
Abbott Laboratories/Hospira/ICU Medical			
 015-0729-00	Single pressure integrated cable	400 cm/13 ft	1 Each
015-0730-00	Dual pressure integrated cable	400 cm/13 ft	1 Each
Interface Cable & Adapters			
 Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)			
690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
Interface Adapter- Transducer interface cable required from transducer company			
 700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
Transducer Kits & Accessories			
Pressure Transducer Kits, Single use			
 690-0008-00	Transducer Kit, Patient Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>	152.5 cm/60 in	5/Box
 690-0009-00	Transducer Kit, Pole Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>	152.5 cm/60 in	5/Box

Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code Description Length/Size Qty

