










## Pulse Oximetry



**TruLink SpO2 Supplies & Accessories**  
**Spacelabs TruLink Technology**

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO2 supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can




Order Code	Description	Length/Size	Qty
<b>Spacelabs TruLink Technology</b>			
<b>Cables</b>			
 700-0030-00	TruLink SpO <sub>2</sub> cable	305 cm/10 ft.	1 Each
<b>Reusable sensors</b>			
 015-0660-00	TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
 015-0661-00	TruLink SpO <sub>2</sub> Y multi-site sensor <u>Includes:</u> (1) Sensor (1) ear clip (10) adult finger/neonate feet/hand wraps (1) pediatric finger or infant toe wraps	Universal Adult>30 kg Pediatric 10-50kg Infant 3-20kg Neonatal <3kg	1 Each
<b>Disposable Sensors</b>			
 015-0662-00	TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
 015-0663-00	TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
 015-0664-00	TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
 015-0665-00	TruLink SpO <sub>2</sub> vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
 015-0666-00	TruLink SpO <sub>2</sub> cloth sensor	Infant 3-20 kg/7-44 lbs.	24/Box
 015-0667-00	TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
<b>Dixtal Cables</b>			

**TruLink SpO2 Supplies & Accessories**  
**Spacelabs TruLink Technology**

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO2 supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can






	<b>Order Code</b>	<b>Description</b>	<b>Length/Size</b>	<b>Qty</b>
	700-0029-00	TruLink SpO <sub>2</sub> cable compatible with Novamatrix® sensors	30 cm/12 in.	1 Each

## Masimo® SpO2 Supplies & Accessories

### LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M). Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.











	Order Code	Description	Length/Size	Shelf Life*	Qty
<b>LNCS Masimo SET Technology</b>					
	<b>LNCS Adapter Cable</b>				
	700-0906-01	TruLink SpO <sub>2</sub> cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft		1 Each
	<b>LNCS Reusable Sensors</b>				
	690-0247-00	Masimo LNCS DC-I, SpO2 digit sensor	Adult > 30 kg		1 Each
	690-0248-00	Masimo LNCS DC-IP, SpO2 digit sensor	Pediatric 10-50 kg		1 Each
	690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO2 ear sensor	Adult > 30 kg		1 Each
	690-0291-00	Masimo LNCS Y-I multisite sensor	> 1 kg	36 months	1 Each
<b>LNCS Masimo SET Technology</b>					
	<b>LNCS Disposable Sensors</b>				
	690-0251-00	Masimo LNCS Amtx, adhesive SpO2 sensor, 18 in.	Adult > 30 kg	36 months	20/Box
	690-0252-00	Masimo LCNS Pmtx, adhesive SpO2 sensor, 18 in.	Pediatric 10 - 50 kg	36 months	20/Box
	690-0258-00	Masimo LNCS Inf-L, adhesive SpO2 sensor, 3 ft.	Infant 3 kg > 20 kg	36 months	20/Box
	690-0259-00	Masimo LNCS Neo-L, adhesive SpO2 sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	36 months	20/Box

# PULSE OXIMETRY Masimo RD SET™ Technology

## Masimo® SpO<sub>2</sub> Supplies & Accessories RD SET Technology

For use with Ultraview SL™ Command Module (91496-M). Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.








Order Code	Description	Length/Size	Life*	Qty
<b>RD SET Masimo Technology</b>				
<b>RD SET Adapter Cable</b>				
 010-1649-00	Masimo RD SET SpO <sub>2</sub> Cable	12 ft.		1 Each
<b>RD SET Reusable Sensors</b>				
 690-0267-00	Masimo RD SET DCI reusable SpO <sub>2</sub> sensor	Adult > 30 kg		1 Each
 690-0270-00	Masimo RD SET DCI-P reusable SpO <sub>2</sub> sensor	Pediatric 10-50 kg		1 Each
 690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO <sub>2</sub> sensor	Adult > 30 kg		1 Each
 690-0281-00	Masimo RD SET Y-I multi-site reusable SpO <sub>2</sub> sensor	> 1 kg	36 months	1 Each
<b>RD SET Masimo Technology</b>				
<b>RD SET Adhesive Disposable Sensors</b>				
 RD SET Adt 690-0283-00	Masimo RD SET Adt, adhesive SpO <sub>2</sub> sensor	Adult, > 30 kg finger or toe	36 months	20/Box
 RD SET Inf 690-0285-00	Masimo RD SET Inf, adhesive SpO <sub>2</sub> sensor	Infant, 3—20 kg, thumb or great toe	36 months	20/Box
 RD SET Neo 690-0286-00	Masimo RD SET Neo, adhesive SpO <sub>2</sub> sensor	Neonatal <3 kg, hand or foot Adult > 40 kg finger or toe	36 months	20/Box

**Nellcor® SpO<sub>2</sub> Supplies & Accessories**  
**Nellcor™ OxiMax® Technology**

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor. Cables and connectors are color-coded purple. **Note, use only genuine Nellcor sensors.**



Order Code	Description	Length/Size	Qty
<b>Nellcor OxiMax Technology</b>			
<b>Cable</b>			
	700-0792-00 TruLink SpO <sub>2</sub> cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
<b>Reusable Sensors</b>			
<b>Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien), now part of Medtronic.</b>			
	690-0003-01 Nellcor DS-100A, finger SpO <sub>2</sub> sensor	Adult	1 Each
	690-0004-01 Nellcor OXI-A/N, SpO <sub>2</sub> sensor	Adult/Neonatal	1 Each
	690-0039-01 Nellcor OXI-P/I, SpO <sub>2</sub> sensor	Pediatric/Infant	1 Each
	690-0215-01 Nellcor DS-Y, multi-site SpO <sub>2</sub> sensor	Universal	1 Each