


Temperature










TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Monitor Adapter Cables			
Adapters			
 700-0031-00	TruLink Dual Temperature Cable. *Required for all reusable probes and as directed for single-patient-use disposable system	30 cm/12 in.	1 Each
Reusable Probes/Sensors			
Skin Surface			
 690-0295-00	YSI 400 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0296-00	YSI 700 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0297-00	YSI 700 series compatible skin surface sensor temperature probe, pediatric size	305 cm/10 ft	1 Each
Esophageal / Rectal			
 690-0294-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
 690-0298-00	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
 690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each












Indicates discontinued part, available while supplies last

TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.




Order Code	Description	Length/Size	Qty
Single-Patient use temperature system			
Extension cables for single use probes			
 700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors, adapter required.	305 cm/10 ft	1 Each
 700-0401-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors. Direct connection to Command module, no adapter required.	305 cm/10 ft	1 Each
Skin			
 701-4000-00	Skin Temperature Sensor	Sterile Yes	Shelf Life* 36 months 20/Box
Esophageal/ Stethoscope			
 701-4009-00	9 French, YSI 400 series compat	Yes	36 months 20/Box
 701-4012-00	12 French, YSI 400 series compa	Yes	36 months 20/Box
 701-4018-00	18 French, YSI 400 series compa	Yes	36 months 20/Box
Esophageal/ Rectal			
 701-4109-00	9 French, YSI 400 series compat	Yes	36 months 20/Box
 701-4112-00	12 French, YSI 400 series compa	Yes	36 months 20/Box
Foley catheter with temperature probe			
 701-4214-00	14 French, YSI 400 series compa	Yes	36 months 20/Box
Tympanic temperature sensor			
 701-4300-00	Adult, YSI 400 series compatible	Yes	36 months 20/Box
* - Shelf life is from date of manufacture, printed on product label as an expiration date			
Exergen			
For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) & DM4 (91331) Monitors			
 010-2157-00	Temporal Artery Thermometer		1 Each

TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Reusable Temperature Probes



Single-Patient Use Temperature Probes

