



October 2024

Supplies & Accessories Product Catalog

An OSI Systems Company


© 2024 Spacelabs Healthcare Specifications subject to change without notice.

Spacelabs Healthcare reserves the right to make changes to the product at any time without notice. All brands and product names and trademarks of their respective owners. Product may not be available for sale in all locations. Check with your Spacelabs representative.

SPACELABS
HEALTHCARE

www.spacelabshealthcare.com

Non-Invasive Blood Pressure Cuffs



TruLink® Nylon Reusable

TruLink reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. TruLink cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, TruLink cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. TruLink cuffs are naturally latex and PVC-free.



**Range/
Length Qty**

Order Code Description Size

TruLink Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

715-1111-10	Single tube, HP	Infant	8 - 14 cm	1 Each
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each



Bayonet

Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



BLOOD PRESSURE CUFFS

TruLink® Soft

TruLink® Soft Disposable

TruLink fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

TruLink Soft Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

715-1113-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Bayonet

Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



TruLink® Vinyl Disposable

TruLink vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

TruLink Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

715-1112-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1122-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Bayonet

Single Tube Hose

714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



Blood Pressure Cuff Accessories

Spacelabs non-invasive blood pressure cuff accessories are compatible with all Spacelabs and TruLink branded non-invasive blood pressure cuffs. Choose from a variety of hoses, adapters and conversion accessories which offer ease-of-use and flexibility.

Order Code	Description	Size	Length	Qty
Non-Invasive Blood Pressure Hose				
Single Tube Hose				
<i>For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors</i>				
714-0018-02	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-02	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0019-02	90496/91496	Neonatal	275 cm/9 ft	1 Each



Double Tube to Single Tube Hose

For use with Command Modules 90496/91496

714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each
-------------	-------------	----------	-------------	--------



Hose Connectors & Adapters

Slip Lock Fitting

Fitting for Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.

714-0038-00	3mm ID Adult NIBP hoses			1 Each
714-0047-00	4mm ID Adult NIBP hoses			1 Each
714-0046-00	5mm ID Adult NIBP hoses			1 Each
714-0081-00	DUAL FEMALE SCREW TO SINGLE MALE BAYONET Y ADAPT <i>Use dual tube GE/Critikon cuffs with Spacelabs monitors</i>			1 Each
714-0083-00	MQML TO SINGLE FEMALE BAYONET Y ADAPTER			1 Each
714-0084-00	DUAL MQ TO SINGLE FEMALE BAYONET Y ADAPTER			1 Each







Indicates discontinued part, available while supplies last

BLOOD PRESSURE CUFFS UltraCheck® Nylon

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.






	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
Single Tube					
 Plastic Bayonet	US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
 Metal Bayonet	US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
	US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
	US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
	US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
	US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
	US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
	US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
	US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag
 Threaded Screw	US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
	US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
	US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
	US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
	US3544DM-05	Single tube, DM	Large Adult Long	32 - 42 cm	5/Bag
 Female Subminiature Twist Lock	US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
	US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
	US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
	US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
	US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
	Indicates discontinued part, available while supplies last				

BLOOD PRESSURE CUFFS UltraCheck® Nylon

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.








	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
 Twist Lock	Single Tube				
	US0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
	US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
 No Fitting	US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
	US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
 Manual Bulb	US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
	US1320MB	Single tube, MB w/gauge	Child	13 - 20 cm	1 Each
	US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
	US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each
	Indicates discontinued part, available while supplies last				

BLOOD PRESSURE CUFFS UltraCheck® Nylon

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
Double Tube					
	UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
Threaded Screw	UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
	UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
	UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
	UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
	UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
Subminiature Twist Lock	UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
	UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
	UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag
	UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
	UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
	Male & Female Mated				
	UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
	UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
	No Fitting				
	UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
	UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
Manual Bulb					
Indicates discontinued part, available while supplies last					

BLOOD PRESSURE CUFFS SoftCheck® Fabric

SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



	Order Code	Description	Size	Range/ Length	Qty
	Single Tube				
	ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
	ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
	ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
	ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
	ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
	ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
	ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
	ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
	ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box
	ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
	ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
	ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
	ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box



Indicates discontinued part, available while supplies last

BLOOD PRESSURE CUFFS SoftCheck® Fabric

SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

SoftCheck Fabric Single Patient Multiple Use BP Cuffs



No Fitting

Single Tube				
ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag



Threaded Screw

Double Tube				
DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box



Twist Lock

DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box
-------------	-----------------	-------	------------	--------



No fitting

DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
-----------	-------------------------	------------	------------	-------

Indicates discontinued part, available while supplies last

BLOOD PRESSURE CUFFS

SoftCheck® Vinyl

SoftCheck® Vinyl Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

SoftCheck Vinyl Disposables BP Cuffs

Single Tube



Plastic Bayonet

VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag



Male
Twist Lock

VST0814MQML5-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1826MQML5-20	Single tube, MQ	Small Adult	18 - 26 cm	20/Box
VST2635MQML-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
VST2635MQML5-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
VST4250MQML5-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box

Indicates discontinued part, available while supplies last

BLOOD PRESSURE CUFFS

UltraCheck® & SoftCheck® Curve

UltraCheck® Curve Nylon Reusable SoftCheck® Curve Fabric Disposable

Your patients come in all shapes and sizes. That's why we offer the Curve cuff for extra-large adults. Its conical-shaped design provides a better fit for the contours of your patient's arm. This can increase patient comfort while providing a more accurate reading.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

UltraCheck Curve Nylon Reusable BP Cuff

Single Tube				
BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MB	Single tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each
Double Tube				
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MB	Double tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

SoftCheck Curve Fabric Disposable BP Cuff

Single Tube				
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Box
BRST3854HP-05	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	5/Box

Indicates discontinued part, available while supplies last



Plastic Bayonet (HP)

Metal Bayonet (HP)

Female Screw Threaded (DM)

Female Locking (MQ)

Male Locking (MQML)

Mated Locking (MQMF)

Manual Bulb w/Gauge (MB)

BLOOD PRESSURE CUFFS SoftCheck® Fabric Neonatal

SoftCheck® Fabric Neonatal Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

SoftCheck Fabric Neonatal Disposable BP Cuffs

Single Tube						
	Micro-taper slip-fit	NN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
		NN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
		NN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
		NN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
		NN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
	Male Luer	NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
		NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
		NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
		NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
		NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
	Plastic Bayonet	NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
		NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
		NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box
<i>Indicates discontinued part, available while supplies last</i>						
Double Tube						
		NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box







BLOOD PRESSURE CUFFS

SoftCheck® Vinyl Neonatal

SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Vinyl Neonatal Disposable BP Cuffs					
Single Tube					
 Micro-taper slip-fit	VNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
 SnapQuik	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 Male Luer	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Plastic Bayonet	VNN1ST-HP-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN3ST-HP-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN5ST-HP-10	Single tube, MQ	Size 4	8-15 cm	10/Box
 Female Twist Lock	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box
Double Tube					
 Male Luer	VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
	VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
	VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
	VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
	VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

Pressure Infusors



Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Disposable Pressure Infusor

Single Patient Use Infusor with aneroid gauge - sold as one piece/complete unit



803FGA-24	Unifusor with thumbwheel valve & aneroid gauge	500mL	24/Box
804FGA-24	Unifusor with thumbwheel valve & aneroid gauge	1000mL	24/Box
805FGA-10	Unifusor with thumbwheel valve & aneroid gauge	3000mL	10/Box



803SGA-24	Unifusor with stopcock valve & aneroid gauge	500mL	24/Box
804SGA-24	Unifusor with stopcock valve & aneroid gauge	1000mL	24/Box
805SGA-10	Unifusor with stopcock valve & aneroid gauge	3000mL	10/Box

Disposable Infusor with piston gauge - sold as one piece/complete unit



903TGA-24	Unifusor with thumbwheel valve & piston gauge	500mL	24/Box
904TGA-24	Unifusor with thumbwheel valve & piston gauge	1000mL	24/Box
905TGA-10	Unifusor with thumbwheel valve & piston gauge	3000mL	10/Box



903SGA-24	Unifusor with stopcock valve & piston gauge	500mL	24/Box
904SGA-24	Unifusor with stopcock valve & piston gauge	1000mL	24/Box
905SGA-10	Unifusor with stopcock valve & piston gauge	3000mL	10/Box

Indicates discontinued part, available while supplies last

Unifusor II Disposable Pressure Infusor with Reusable Inflation System

Disposable Infusor - white Disposable bag with reusable inflation system



2803SGA-24	Unifusor II with aneroid gauge & stopcock valve	500mL	24/Box
2804SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000mL	24/Box
2805SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000mL	10/Box



2903SGA-24	Unifusor II with piston gauge & stopcock valve	500mL	24/Box
2904SGA-24	Unifusor II with piston gauge & stopcock valve	1000mL	24/Box
2905SGA-10	Unifusor II with piston gauge & stopcock valve	3000mL	10/Box



Disposable Replacement Parts

2803X-24	Unifusor II Disposable white infusion bag	500mL	24/Box
----------	---	-------	--------


Unifusor® Pressure Infusors

Unifusor Pressure Infusors are available in 3 Models: Unifusor Disposable, Unifusor II & Unifusor Plus Reusable. All models are available in 500, 1000 and 3,000 mL sizes. The Unifusors are 100% tested for leakage during the manufacturing process. Unifusor pressure Infusors include an IV pole loop, a fluid bag hood and a reinforced translucent sleeve which enables easy viewing of fluid type and level. This unit is designed for rapid IV infusion, invasive pressure monitoring and auto-transfusion procedures.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with white hose

1103SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	500mL	24/Box
1104SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	1000mL	24/Box
1113SGA-24	Unifusor Plus with piston gauge & stopcock valve	500mL	24/Box

Reusable replacement parts



1103XB-05	Unifusor Plus blue reusable bag with white hose	500mL	5/Box
1104XB-05	Unifusor Plus blue reusable bag with white hose	1000mL	5/Box
1105XB-05	Unifusor Plus blue reusable bag with white hose	3000mL	5/Box

Indicates discontinued part, available while supplies last

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with black hose

1105	Unifusor Plus with trigger / gauge assembly	3000mL	1 Each
------	---	--------	--------

Reusable replacement parts

1103X	Unifusor Plus blue reusable bag with black hose	500mL	1 Each
1104X	Unifusor Plus blue reusable bag with black hose	1000mL	1 Each
1105X	Unifusor Plus blue reusable bag with black hose	3000mL	1 Each

Electrocardiography


AAMI COLOR CODING



ECG Combiner Lead Wires & Cables

ECG COMBINER LEAD WIRE SYSTEM AAMI Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Disposable, DIN Shielded, Pinch, AAMI

700-0006-50	3-Lead	Pinch	127cm /50 in	25 Sets
* 700-0006-32	5 Lead	Pinch	127cm /50 in	25 Sets
* 700-0006-33	5 Lead (V-Lead)	Pinch	127cm /50 in	25 Sets



Disposable, DIN Shielded, Snap, AAMI

700-0006-61	5 Lead	Snap	127cm /50 in	25 Sets
-------------	--------	------	--------------	---------



Reusable, DIN Shielded, Pinch, AAMI

* 700-0006-35	3 Lead	Pinch	74cm /29 in	1 Set
* 700-0006-37	5 Lead	Pinch	74cm /29 in	1 Set
* 700-0006-38	5 Lead (V-Lead)	Pinch	74cm /29 in	1 Set
700-0010-00	3 Lead	Pinch	130cm/51 in	1 Set
* 700-0010-01	5 Lead	Pinch	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Pinch	130cm/51 in	1 Set



Reusable, DIN Shielded, Snap, AAMI

* 700-0007-35	3 Lead	Snap	74cm /29 in	1 Set
* 700-0007-37	5 Lead	Snap	74cm /29 in	1 Set

* Product not available for sale into CE Mark Countries (EU)

ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter




Radiolucent, Reusable, DIN Unshielded, Pinch, AAMI

700-0010-06	5 Lead, Radiolucent	Pinch	122cm/48 in	1 Set
700-0010-07	3 Lead, Radiolucent	Pinch	183cm/72 in	1 Set
700-0010-08	5 Lead, Radiolucent	Pinch	183cm/72 in	1 Set

ECG Combiner Lead Wires & Cables

ECG COMBINER LEAD WIRE SYSTEM AAMI Color Coded

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496					
* 700-0008-55	3 Lead	Spacelabs	305 cm/10 ft	1 Cable	
* 700-0008-57	5 Lead	Spacelabs	305 cm/10 ft	1 Cable	
700-0008-32	5 Lead	Spacelabs	457 cm/15 ft	1 Cable	
* 700-0008-56	10 Lead (5+5)	Spacelabs	305 cm/10 ft	1 Cable	
* Product not available for sale into CE Mark Countries (EU)					

ECG Combiner Adapters



Spacelabs Telemetry, AAMI			
700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)	1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)	1 Adapter

ECG Combiner Cables for Standardization




Other manufacturers					
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable	
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable	

Indicates discontinued part, available while supplies last

ECG LEAD WIRES & CABLES AAMI Color Coded - Shielded

ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

Snap Sets

700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set



Pinch

Pinch Sets

700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set

ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI




24" Snap Individual Color

012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

ECG LEAD WIRES & CABLES AAMI Color Coded - Shielded

ECG Lead Wire & Cables - Shielded

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Cables, Shielded, AAMI

Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496



700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

ECG Accessories



Electrodes

717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 pouches <i>24 month shelf life from date of manufacture</i>	Adult	50/Box
-------------	--	-------	--------

Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable

700-0008-12	3 Lead	Neonatal	1 Cable
-------------	--------	----------	---------


DIN Unshielded, Neonatal Electrodes.

717-0080-00	Pre-wired neonatal electrodes, Pre-gelled, silver/silver chloride, single use <i>24 month shelf life from date of manufacture</i>	3-Lead	100 Sets/Box
-------------	--	--------	--------------

DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap Sets

015-0732-00	3 Lead	90 cm/35 in	1 Set
015-0736-00	5 Lead	90 cm/35 in	1 Set



Pinch Sets

015-0731-00	3 Lead	90 cm/35 in	1 Set
012-0284-11	3 Lead	61 cm/24 in	1 set
015-0734-00	5 Lead	90 cm/35 in	1 Set

Indicates discontinued part, available while supplies last

Telemetry Transmitter Adapters



For use with DIN Unshielded Lead Wires

706-0191-00	Aria Telemetry Grouper (96281)	5/Bag
-------------	--------------------------------	-------

DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Telemetry Accessories



Protective Solutions and Electrodes

010-2350-00	Aria Tele Bumper		5/Box
-------------	------------------	--	-------



015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
-------------	-----------------------------------	--	---------



717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches <i>24 month shelf life from date of manufacture</i>	Adult	50/Box
-------------	--	-------	--------

DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY AAMI Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

SpO₂ Cable & Sensors

Note: Only the below genuine Spacelabs and Nellcor branded SpO₂ sensors are compatible with Spacelabs Telemetry systems

SpO₂ Cable



700-0014-00	TruLink SpO ₂ Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO ₂ sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

Reusable SpO₂ Sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO ₂ sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

Disposable SpO₂ Sensors



015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last

Electrocardiography


IEC COLOR CODING



ECG Combiner Lead Wires & Cables

Our ECG Combiner System of lead wires, cables and adapters allows you to standardize ECG supplies across monitoring brands. The system features a universal connection between cables and lead wires which allows a quick connect/disconnect between monitors and telemetry, and single patient multiple use lead wire options. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.

ECG COMBINER LEAD WIRE SYSTEM IEC Color Coded



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Multi-pinch

Reusable, DIN Shielded, Multi-Pinch, IEC

700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set
700-0006-43	5 Lead	Multi-Pinch	74 cm /29 in	1 Set
700-0006-44	5 Lead	Multi-Pinch	127cm /50 in	1 Set



Snap

Reusable, DIN Shielded, Snap, IEC

700-0007-41	3 Lead	Snap	74cm /29 in	1 Set
-------------	--------	------	-------------	-------

Indicates discontinued part, available while supplies last

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)

700-0008-81	3 Lead	IEC	305 cm/10 ft	1 Cable
700-0008-83	5 Lead	IEC	305 cm/10 ft	1 Cable
700-0008-84	10 Lead	IEC	305 cm/10 ft	1 Cable

ECG LEAD WIRES & CABLES IEC Color Coded - Shielded

ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, Shielded, Reusable, IEC



Snap

Snap Sets

700-0007-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0007-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0007-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0007-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0007-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0007-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set



Pinch

Pinch Sets

700-0006-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0006-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0006-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0006-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0006-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0006-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set

ECG Cables, Shielded, IEC



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496

700-0008-09	3 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-07	5 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-01	10 Lead, shielded	305 cm/10 ft	1 Cable

Indicates discontinued part, available while supplies last

ECG Accessories



Electrodes



717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches <i>24 month shelf life from date of manufacture</i>	Adult	50/Box
-------------	--	-------	--------

ECG LEAD WIRES & CABLES IEC Color Coded - Shielded

ECG Lead Wires & Cables - Shielded

IEC color coded shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.




Order Code	Description	Length	Qty
Neonatal ECG Supplies & Accessories			
	DIN Unshielded, Neonatal ECG Cable		
700-0008-13	3 Lead		1 Cable
	DIN Unshielded, Neonatal Electrodes		
717-0080-00	Pre-wired neonatal electrodes, pre-gelled, silver/silver chloride, CE <i>24 month shelf life from date of manufacture</i>	3-Lead	100 Sets/Box

DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY IEC Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC



Snap

Snap Sets			
015-0733-00	3 Lead	90 cm/35 in	1 Set
012-0328-10	5 Lead	61 cm/24 in	1 Set
015-0737-00	5 Lead	90 cm/35 in	1 Set



Pinch Sets			
015-0735-00	5 Lead	90 cm/35 in	1 Set

For use with DIN Unshielded Lead Wires			
706-0190-00	Aria Telemetry Grouper (96281), IEC		5/Bag

Telemetry Accessories



Disposable Pouches and Electrodes			
010-2350-00	Aria Tele Bumper		5/box
015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches <i>24 month shelf life from date of manufacture</i>	Adult	50/Box


Bold font indicates a new part number

Indicates discontinued part, available while supplies last

DIN Unshielded Lead Wires & Telemetry

TruLink ECG Din unshielded reusable lead wires are validated for use with Spacelabs telemetry systems. They are available in 3 and 5-lead sets, or available to purchase individually in color-coded multi-packs. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.

ECG LEAD WIRES & TELEMETRY IEC Color Coded - DIN Unshielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

SpO₂ Cable & Sensors

SpO₂ Cable



700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

Note: Only the below genuine Spacelabs and Nellcor branded SpO₂ sensors are compatible with Spacelabs Telemetry systems

Reusable SpO₂ Sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------



690-0003-01	Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each
-------------	--	-------------------------	--------

Disposable SpO₂ Sensors



015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

Indicates discontinued part, available while supplies last

Capnography & Multigas



NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

NomoLine LH Nasal CO2 Cannula		Shelf Life*		
-------------------------------	--	-------------	--	--



704-3817-00	Adult	24 months	2 m/7 ft	25/Box
704-3819-00	Infant	24 months	2 m/7 ft	25/Box
704-3820-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4624-00	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3821-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box



NomoLine LH Nasal/Oral CO2 Cannula				
------------------------------------	--	--	--	--

704-3824-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4628-00	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3825-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box



NomoLine LH Single Nasal Prong CO2 Cannula				
--	--	--	--	--

704-3826-00	Adult	24 months	2 m/7 ft	25/Box
-------------	-------	-----------	----------	--------

NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine LH Airway Adapter Set		Shelf Life*		
--------------------------------	--	-------------	--	--

704-3814-00	Adult/Pediatric	36 months	2 m/7 ft	25/Box
704-4620-00	Adult/Pediatric	36 months	4 m/13 ft	25/Box

NomoLine LH (Low Humidity) Procedural Oxygen Masks



NomoLine LH Nasal CO2 Cannula		Shelf Life*		
-------------------------------	--	-------------	--	--

704-5213-00 ^{USA}	Adult, Medium O2 Flow Rate	24 months	3 m/10 ft	30/Box
704-5214-00 ^{USA}	Pediatric, Medium O2 Flow	24 months	3 m/10 ft	30/Box
704-5215-00 ^{USA}	Adult, High O2 Flow Rate	24 months	3 m/10 ft	30/Box
704-5216-00 ^{USA}	Pediatric, High O2 Flow Rate	24 months	3 m/10 ft	30/Box

* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA

NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, with functionality in any orientation, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients



NomoLine HH Nasal CO2 Cannula		Shelf Life*		
704-3830-00	Adult	24 months	2 m/7 ft	25/Box
704-3832-00	Infant	24 months	2 m/7 ft	25/Box
704-3833-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-3834-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box
704-4933-00 ^{USA}	Infant w/O2	24 months	4 m/13 ft	25/Box




NomoLine HH Nasal/Oral CO2 Cannula				
704-3837-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4968-00 ^{USA}	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3838-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box



NomoLine HH Single Nasal Prong CO2 Cannula				
704-3839-00	Adult	24 months	2 m/7 ft	25/Box

NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine HH Airway Adapter Set		Shelf Life*		
704-3827-00	Adult/Pediatric	36 months	2 m/7 ft	25/Box
704-3828-00	Adult/Pediatric	36 months	3 m/10 ft	25/Box
704-4989-00 ^{USA}	Adult/Pediatric	36 months	4 m/13 ft	25/Box
704-4367-00	Infant/Neonatal	36 months	2 m/7 ft	25/Box

* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Sidestream Monitoring Supplies

Please note: the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.



CO2 Nasal Cannula with luer		Shelf Life*		
* 704-0176-00	Adult	60 months	2 m/7 ft	25/Box
CO2 Nasal Cannula w/O2 and luer				
* 704-0181-00	Adult	60 months	2 m/7 ft	25/Box
CE 704-0181-01	Adult	36 months	2 m/7 ft	25/Box
* 704-0180-00	Pediatric	60 months	2 m/7 ft	25/Box
** 704-0179-00	Infant	60 months	2 m/7 ft	25/Box
** 704-0186-00	Adult	60 months	3 m/10 ft	25/Box
CO2 Oral/Nasal Cannula w/O2 and luer				
* 704-0183-00	Adult	60 months	2 m/7 ft	25/Box
CE 704-0183-01	Adult	36 months	2 m/7 ft	25/Box
** 704-0187-00	Adult	60 months	3 m/10 ft	25/Box

Indicates discontinued part, available while supplies last
 * - Shelf life is from date of manufacture, printed on product label as an expiration date

Sidestream Monitoring Supplies - Intubated patients



Airway adapters and lines		Shelf Life*		
** 704-0191-00	Adult Line w/ T-Adapter	60 months	2 m/7 ft	25/Box
** 704-0192-00	Adult Line w/ T-Adapter	60 months	3 m/10 ft	25/Box
* 166-7121-00	sampling line	60 months	3 m/10 ft	50/Box
* 704-0175-00	T-Adapter	60 months		25/Box

* Product not available for sale into CE Mark Countries (EU)
 CE product is available for sale into CE Mark Countries (EU)
 ** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients.
Featuring Masimo IRMA mainstream technology



Sensor & Airway Adapters - required for mainstream monitoring

010-1980-00	Masimo IRMA mainstream sensor, reusable		1 Each
704-0173-00	Masimo IRMA mainstream airway adapter, single use <i>36 month shelf life from date of manufacture</i>	Adult/Pediatric	25/Box
704-0174-00	Masimo IRMA mainstream airway adapter, single use <i>36 month shelf life from date of manufacture</i>	Infant	10/Box

CAPNOGRAPHY





Spacelabs Module 91517

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCO₂ capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream, Non-intubated, single patient use			
		CO₂ Nasal Cannula	
		Shelf Life[^]	
* 704-0010-00	Adult	36 months	10/Box
* 704-0012-00	Pediatric	36 months	10/Box
		CO₂ Nasal Cannula w/O₂ Delivery	
* 704-0014-00	Adult	36 months	10/Box
CE 704-0014-10	Adult, (for sale into EU/UK only)	36 months	10/Box
* 704-0016-00	Pediatric	36 months	10/Box
* 704-0015-00	Infant/Neonatal	36 months	10/Box
		CO₂ Oral/Nasal Cannula	
* 704-0017-00	Adult	36 months	10/Box
		CO₂ Oral/Nasal Cannula w/O₂ Delivery	
* 704-0020-00	Adult	36 months	10/Box
* 704-0022-00	Pediatric	36 months	10/Box
Sidestream, Intubated, single patient use			
Sidestream - short term up to 12 hours			
* 704-0023-00	Airway adapter + sampling line, ET tube > 4.0 mm	36 months	10/Box
CE 704-0023-10	Airway adapter + sampling line, ET tube > 4.0 mm	36 months	10/Box
* 704-0024-00	Airway adapter + sampling line, Pediatric/Infant sized	36 months	10/Box
* 704-0027-00	CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box
CE 704-0027-10	CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box

CAPNOGRAPHY Spacelabs Module 91517

Capnography Accessories

The dual capnography module 91517 provides both Sidestream and Mainstream EtCo2 capabilities in a single module. Spacelabs carries validated TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream - long term up to 120 hours		Shelf Life[^]	
* 704-0025-00	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	36 months	10/Box
CE 704-0025-10	Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	36 months	10/Box
* 704-0026-00	Airway adapter + sampling line with Nafion, Infant	36 months	10/Box
* 704-0028-00	CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box
CE 704-0028-10	CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box

Mainstream, for Intubated Patients Only



Mainstream Airway Adapters, for use with CO₂ Sensor 704-0001-00			
* 704-0005-00	CO ₂ airway adapter, single use, Adult	36 months	24/Box
* 704-0002-00	CO ₂ airway adapter, reusable, Adult	36 months	1 Each

Indicates discontinued part, available while supplies last

[^] - Shelf life is from date of manufacture, printed on product label as an expiration date

* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

** Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

MULTIGAS ANALYZER Spacelabs Modules 90518 & 91518

Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).

These two legacy multi-gas modules are obsolete. However, Spacelabs continues to provide supplies that are available.



(90518)



(91518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Multigas Supplies & Accessories

Anesthetic Agent Samples Lines



015-0128-00	SAMPLE LINE,MALE TO MALE,LUER,SPU,7'X.04",25/BX		25/Box
015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line:0.048 inches ID tubing (90518/91518) <i>36 month shelf life from date of manufacture</i>	305 cm/10 ft.	50/Box
162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 0.048 in ID tubing (90518/91518) <i>36 month shelf life from date of manufacture</i>	305 cm/10 ft.	50/Box



Water Traps, Reusable

717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)		10/Pack
717-0448-00	D-FEND water trap cartridge, Green. Replace at least once every 24 hours or when occluded (91518)		10/Pack

Note: manufacturer (GE) has discontinued these parts. Available while supplies last.










Pulse Oximetry



TruLink SpO2 Supplies & Accessories
Spacelabs TruLink Technology

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO2 supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can




Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology			
Cables			
 700-0030-00	TruLink SpO ₂ cable	305 cm/10 ft.	1 Each
Reusable sensors			
 015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs.	1 Each
 015-0661-00	TruLink SpO ₂ Y multi-site sensor <u>Includes:</u> (1) Sensor (1) ear clip (10) adult finger/neonate feet/hand wraps (1) pediatric finger or infant toe wraps	Universal Adult>30 kg Pediatric 10-50kg Infant 3-20kg Neonatal <3kg	1 Each
Disposable Sensors			
 015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs.	24/Box
 015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs.	24/Box
 015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
 015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
 015-0666-00	TruLink SpO ₂ cloth sensor	Infant 3-20 kg/7-44 lbs.	24/Box
 015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
Dixtal Cables			

TruLink SpO2 Supplies & Accessories
Spacelabs TruLink Technology

TruLink branded pulse oximetry supplies and accessories are available for Spacelabs command modules (90496-U/91496-U). These cables and sensors are designed specifically for the Spacelabs module and provide a high level of validated reliable and accurate data capture. Economical and high quality, TruLink SPO2 supplies are available in reusable and disposable sensors for all patient types, including a popular reusable multi-site sensor that depending on patient type can






	Order Code	Description	Length/Size	Qty
	700-0029-00	TruLink SpO ₂ cable compatible with Novamatrix® sensors	30 cm/12 in.	1 Each

Masimo® SpO2 Supplies & Accessories

LNCS® Masimo SET® Technology

For use with Ultraview SL™ Command Module (91496-M). Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.











	Order Code	Description	Length/Size	Shelf Life*	Qty
LNCS Masimo SET Technology					
	LNCS Adapter Cable				
	700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft		1 Each
	LNCS Reusable Sensors				
	690-0247-00	Masimo LNCS DC-I, SpO2 digit sensor	Adult > 30 kg		1 Each
	690-0248-00	Masimo LNCS DC-IP, SpO2 digit sensor	Pediatric 10-50 kg		1 Each
	690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO2 ear sensor	Adult > 30 kg		1 Each
	690-0291-00	Masimo LNCS Y-I multisite sensor	> 1 kg	36 months	1 Each
LNCS Masimo SET Technology					
	LNCS Disposable Sensors				
	690-0251-00	Masimo LNCS Amtx, adhesive SpO2 sensor, 18 in.	Adult > 30 kg	36 months	20/Box
	690-0252-00	Masimo LCNS Pmtx, adhesive SpO2 sensor, 18 in.	Pediatric 10 - 50 kg	36 months	20/Box
	690-0258-00	Masimo LNCS Inf-L, adhesive SpO2 sensor, 3 ft.	Infant 3 kg > 20 kg	36 months	20/Box
	690-0259-00	Masimo LNCS Neo-L, adhesive SpO2 sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	36 months	20/Box

PULSE OXIMETRY Masimo RD SET™ Technology

Masimo® SpO₂ Supplies & Accessories RD SET Technology

For use with Ultraview SL™ Command Module (91496-M). Cables and connectors are color-coded gray. Note: use only genuine Masimo branded sensors. Masimo sensors are Not compatible with the 91496-U, 91496-N or Aria Tele modules.








Order Code	Description	Length/Size	Life*	Qty
RD SET Masimo Technology				
RD SET Adapter Cable				
 010-1649-00	Masimo RD SET SpO ₂ Cable	12 ft.		1 Each
RD SET Reusable Sensors				
 690-0267-00	Masimo RD SET DCI reusable SpO ₂ sensor	Adult > 30 kg		1 Each
 690-0270-00	Masimo RD SET DCI-P reusable SpO ₂ sensor	Pediatric 10-50 kg		1 Each
 690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO ₂ sensor	Adult > 30 kg		1 Each
 690-0281-00	Masimo RD SET Y-I multi-site reusable SpO ₂ sensor	> 1 kg	36 months	1 Each
RD SET Masimo Technology				
RD SET Adhesive Disposable Sensors				
 RD SET Adt 690-0283-00	Masimo RD SET Adt, adhesive SpO ₂ sensor	Adult, > 30 kg finger or toe	36 months	20/Box
 RD SET Inf 690-0285-00	Masimo RD SET Inf, adhesive SpO ₂ sensor	Infant, 3—20 kg, thumb or great toe	36 months	20/Box
 RD SET Neo 690-0286-00	Masimo RD SET Neo, adhesive SpO ₂ sensor	Neonatal <3 kg, hand or foot Adult > 40 kg finger or toe	36 months	20/Box

Nellcor® SpO₂ Supplies & Accessories
Nellcor™ OxiMax® Technology

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor. Cables and connectors are color-coded purple. **Note, use only genuine Nellcor sensors.**



	Order Code	Description	Length/Size	Qty
Nellcor OxiMax Technology				
	Cable			
	700-0792-00	TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
Reusable Sensors				
Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien), now part of Medtronic.				
   	690-0003-01	Nellcor DS-100A, finger SpO ₂ sensor	Adult	1 Each
	690-0004-01	Nellcor OXI-A/N, SpO ₂ sensor	Adult/Neonatal	1 Each
	690-0039-01	Nellcor OXI-P/I, SpO ₂ sensor	Pediatric/Infant	1 Each
	690-0215-01	Nellcor DS-Y, multi-site SpO ₂ sensor	Universal	1 Each

Temperature










TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Monitor Adapter Cables			
Adapters			
 700-0031-00	TruLink Dual Temperature Cable. *Required for all reusable probes and as directed for single-patient-use disposable system	30 cm/12 in.	1 Each
Reusable Probes/Sensors			
Skin Surface			
 690-0295-00	YSI 400 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0296-00	YSI 700 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
 690-0297-00	YSI 700 series compatible skin surface sensor temperature probe, pediatric size	305 cm/10 ft	1 Each
Esophageal / Rectal			
 690-0294-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
 690-0298-00	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
 690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each












Indicates discontinued part, available while supplies last

TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.




Order Code	Description	Length/Size	Qty
Single-Patient use temperature system			
Extension cables for single use probes			
 700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors, adapter required.	305 cm/10 ft	1 Each
 700-0401-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors. Direct connection to Command module, no adapter required.	305 cm/10 ft	1 Each
Skin			
 701-4000-00	Skin Temperature Sensor	Yes	36 months 20/Box
Esophageal/ Stethoscope			
 701-4009-00	9 French, YSI 400 series compat	Yes	36 months 20/Box
 701-4012-00	12 French, YSI 400 series compa	Yes	36 months 20/Box
 701-4018-00	18 French, YSI 400 series compa	Yes	36 months 20/Box
Esophageal/ Rectal			
 701-4109-00	9 French, YSI 400 series compat	Yes	36 months 20/Box
 701-4112-00	12 French, YSI 400 series compa	Yes	36 months 20/Box
Foley catheter with temperature probe			
 701-4214-00	14 French, YSI 400 series compa	Yes	36 months 20/Box
Tympanic temperature sensor			
 701-4300-00	Adult, YSI 400 series compatible	Yes	36 months 20/Box
<i>* - Shelf life is from date of manufacture, printed on product label as an expiration date</i>			
Exergen			
For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) & DM4 (91331) Monitors			
 010-2157-00	Temporal Artery Thermometer		1 Each

TEMPERATURE

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.





Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----


Reusable Temperature Probes



Single-Patient Use Temperature Probes









Invasive Blood Pressure



Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.




Order Code	Description	Length/Size	Qty
Integrated Adapter Cables (connects directly to transducer)			
Edward Lifesciences			
 015-0727-00	Single pressure integrated cable	400 cm/13 ft	1 Each
015-0728-00	Dual pressure integrated cable	400 cm/13 ft	1 Each
Abbott Laboratories/Hospira/ICU Medical			
 015-0729-00	Single pressure integrated cable	400 cm/13 ft	1 Each
015-0730-00	Dual pressure integrated cable	400 cm/13 ft	1 Each
Interface Cable & Adapters			
 Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)			
690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
Interface Adapter- Transducer interface cable required from transducer company			
 700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
Transducer Kits & Accessories			
Pressure Transducer Kits, Single use			
 690-0008-00	Transducer Kit, Patient Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>	152.5 cm/60 in	5/Box
 690-0009-00	Transducer Kit, Pole Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>	152.5 cm/60 in	5/Box

Invasive Pressure Supplies


Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers. Transducer kits are available in patient or pole mount options, as well as a neonatal configuration.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----



Cardiac Output



CARDIAC OUTPUT

Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Cardiac Output Cables & Probes

Cables



700-0027-00	TruLink Cardiac Output Adapter Cable. Connects CO cable to 91496 command module.	30 cm/12 in.	1 Each
015-0739-00	Cardiac Output cable Edwards compatible connects with 690-0033-00 temperature probe. Replaces 306655-003. The 700-0027-00 adapter not required.		1 Each
306655-004	Cardiac Output cable Datex/Ohmeda/Abbott compatible connects with 690-0023-00 or 690-0024-00 temperature probe. 700-0027-00 adapter required.		1 Each

Bold font indicates a new part number

Indicates discontinued part, available while supplies last.


Injectate Probes



690-0023-00	In-line injectate temperature probe. Datex/Ohmeda/Abbott connector style. Connects to 306655-004 cable and 690-0293-00 injectate kit integrated thermistor.		1 Each
690-0024-00	Ice bath temperature probe. Datex/Ohmeda/Abbott connector style. Connects to 306655-004 cable and needle end is submerged in an ice bath.		1 Each
690-0033-00	In-line injectate temperature probe. Edwards connector style. Connects to 306655-003 cable and 690-0293-00 injectate kit integrated thermistor.		1 Each

Cardiac Output Injectate Kits

Closed Loop Injectate Kits

690-0293-00		SP5045 Critikit In-Line Sensor, Control Syringe closed loop injectate kit <i>Packaged sterile with 36 month shelf life</i>	5/Box
-------------	---	---	-------

CARDIAC OUTPUT


Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

91496-C Module, Cardiac Output Connection guide



Recording & Printing



Paper

High quality Z-fold paper for use with the Bedside Printer Module (90449) and Central Printer (90469).


Roll paper for use with the Standalone printer (91449) and Spacelabs Portable Monitors.



(91449)
Printer

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Thermal Paper



Z-fold Paper for use with 90449 and 90469

006-0196-00	TruLink Thermal Paper, Z-Fold Pads, 2-channel printing, blank paper, 230 Sheets/pad	27.4 m/90 ft	10 Pads/Box
-------------	---	--------------	-------------



Roll Paper for use with 91449, 91367, 91369, 91370 & 91390

307438-001	TruLink Thermal Paper, Rolls, single- or dual-channel printing, 50 mm/1.97 inch wide blank paper	30.1 m/100 ft	10 Rolls/Box
------------	--	---------------	--------------

Recording Paper

20257-004	Thermal recorder paper, black	HearTwave II	each
010-2246-00	PAPR,THERM,Z-FOLD, 210X295MM X 180 SHTS/10RM,IEC	Cardioline	each

Bispectral Index



Bispectral Index Supplies

For use with Spacelabs BISx™ Module (91482).



Order Code	Description	Shelf Life*	Qty
BISx Cables & Sensors			
BISx Cables			
719-0004-01	BISx Patient Interface cable		1 Each
719-0006-20	BISx Pod		1 Each
719-0007-01	BISx host monitor cable		1 Each

BISx Sensors

Note: Sales of BISx sensors are limited to one box with the purchase of capital equipment. All future purchases should be made directly through Covidien, now part of Medtronic.

719-0003-10	BISx Extend Sensor, single patient use	18 months	25/Box
719-0002-10	BISx Quatro Sensor, single patient use	18 months	25/Box

* - Shelf life is from date of manufacture, printed on product label as an expiration date

Ultraview DM3

DUAL MODE VITAL SIGNS MONITOR



VITAL SIGNS MONITOR Ultraview DM3

Vital Sign Monitoring

Ultraview DM3 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

See previous sections for NIBP and SpO2 supplies:

- The DM3 uses the TruLink NIBP cuffs (single tube, bayonet fitting) and the standard air hose 714-0018-00.
- The DM3 used the same cables and sensors, by brand, for pulse oximetry.




DM3 91330

Order Code	Description	Size	Length/Size	Qty
Supplies unique to the DM3 monitor:				
019-0539-00	Replacement Temperature			1 Each
Indicates discontinued part, available while supplies last				

DM4

VITAL SIGNS MONITOR



VITAL SIGNS MONITOR Ultraview DM4










91331

Vital Sign Monitoring


Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.

The DM4 uses TruLink branded NIBP cuffs (single tube, bayonet fitting) and air hose (714-0018-00)

Order Code	Description	Length/Size	Qty
Supplies unique to the DM4 monitor			
RD SET Adapter Cable			
 010-2193-00	DM4 Masimo RD SET SpO2 cable compatible with Masimo RD SET sensors, see section on RD SET sensors	305 cm/10 ft	1 Each
Nellcor OxiMax Cable			
 010-2192-00	DM4 Nellcor OxiMax SpO2 cable compatible with Nellcor OxiMax sensors, see section on Nellcor sensors	305 cm/10 ft	1 Each
Filac			
 010-2174-00	Filac dock and insulation chamber		1 Each
 010-2175-00	Filac temperature probe		1 Each
 010-2228-00	Filac temperature oral probe covers, 20/bx		20/box 25/boxes
Exergen			
 010-2157-00	Temporal Artery Thermometer		1 Each
NeoNatal NIBP supplies			
 714-0087-00	NIBP,HOSE ASSY,NEONATE,DM4 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	1 Each

Elance

VITAL SIGNS MONITOR



Elace Vital Signs Monitor

Elace Vital Signs Monitor







The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add Microstream end tidal CO2.

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below
- TruLink ECG cables and leadwire systems are compatible



Elace Vital Signs Monitor
93300

Order Code	Description	Qty
Supplies unique to the Elance monitor		
ECG Lead Wires & Cables (can also use standard TruLink ECG cables and leadwire systems)		
ECG Integrated Cables and Lead Wires		
	700-0008-23 5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
	700-0008-24 3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)		
	SpO2 Cable	
700-0287-00	SpO2 adapter cable	305 cm/10 ft 1 Each
LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)		
	LNCS Adapter Cable	
690-0256-00	SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft 1 Each
Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)		
	Cable	
700-0290-00	SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft 1 Each
Temperature		
	Temperature Probe	
690-0261-01	Skin/Surface Temp probe, reusable, YSI-400	305cm/10ft 1 Each
Invasive Blood Pressure (See IBP section for compatible sensors and transducers)		
	Interface Cable	
690-0021-00	Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	366 cm/12 ft 1 Each
	Indicates discontinued part, available while supplies last	

C50

VITAL SIGNS MONITOR












Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
ECG leadwire systems				
 040-000908-00	CABLE,INTEGRATED 3-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000909-00	CABLE,INTEGRATED 5-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000910-00	CABLE,INTEGRATED 10-WIRE CLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000911-00	CABLE,INTEGRATED 3-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000912-00	CABLE,INTEGRATED 5-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
 040-000913-00	CABLE,INTEGRATED 10-WIRE CLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
Pulse Oximetry Supplies				
 040-000869-00	CABLE,INTEGRATED ADT FINGER CLIP SENSOR,C50	3 meters	Reusable	1 Each
 040-000769-00	CABLE,SPO2,12PIN TO DB9F,C50 <i>Sensors sold seperately</i>	3 meters	Reusable	1 Each
 040-000730-00	SENSOR,SPO2,ADT/NEO/PED,WRAP, 3',C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each








Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Non-Invasive Blood Pressure supplies				
 040-000808-00	ADULT NIBP TUBE (WHITE)2M	2 meters	Reusable	1 Each
	Reusable or Disposable Cuffs	See NIBP section under TruLink Reusable or Disposable cuffs. These are compatible with the C50. <i>Requires use of the 040-00808-00 air hose (above).</i>		
 714-0086-00	NIBP,HOSE ASSY,NEONATE,C50 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	Reusable	1 Each
Capnography Supplies: requires external Phasein ISA module, 099-000007-00				
 040-000400-00	CABLE,CO2/AG,C50 Connects the Phasein ISA module to the C50	1 meter	Reusable	1 Each
	For airway adapter, cannula and sampling line sets, please see the capnography supply section under Nomoline SPU supplies. These are compatible with the C50.			
Temperature sensors				
 040-000246-00	PROBE,TEMP,ADT,SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000387-00	PROBE,TEMP,NEO SKIN,10K,C50	3 meters	Reusable	1 Each
 040-000385-00	PROBE,TEMP,ADT,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each
 040-000386-00	PROBE,TEMP,NEO,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each






Spacelabs C50 Patient Monitoring System

Supplies and accessories for the Spacelabs C50 monitoring system.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Invasive Blood Pressure Cables and Transducer kits				
 040-000454-00	INTERFACE CABLE,IBP,EDWARD CONNCTR,C50	 3 meters	Reusable	1 Each
 040-000360-00	INTERFACE CABLE,IBP,ABBOTT CONNCTR,C50	 3 meters	Reusable	1 Each
Recording Paper				
 040-000116-00	RECORDER PAPER,THERMAL,C50	50mm X 20 meters		1 roll Each

Anesthesia



Anesthesia

Anesthesia System For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseSirius Anesthesia Machines.



BleaseFocus



BleaseSirius



ARKON

Order Code	Description	Qty	
Anesthesia Supplies			
Patient Flow Sensors		Shelf Life[^]	
<i>For use with BleaseSirius and BleaseFocus Anesthesia Machines only</i>			
* 10110090	Adult Patient flow sensor		1 Each
717-0428-00	D-Lite reusable air flow sensor, may require luer caps		1 Each
050-0842-00	Luer caps for D-Lite air flow sensor <i>See image for assembly, available for sale to EU</i>		Pack of 5
10110218	2M Patient flow sensor tubing assembly	60 months	1 Each
10110102	3M Patient flow sensor tubing assembly	60 months	1 Each
10110103	6M Patient flow sensor tubing assembly	60 months	1 Each
70300001	O2 probe with tee		1 Each
376-0561-01	Pediatric flow sensor kit *EU/UK ONLY*		1 Each
FR101029	Pediatric flow sensor and tubing kit, 2.0 meters	Yes	1 Each
FR101030	Pediatric flow sensor and tubing kit, 3.0 meters	Yes	1 Each
Indicates discontinued part, available while supplies last			
* Product not available for sale into CE Mark Countries (EU)			
[^] - Shelf life is from date of manufacture, printed on product label as an expiration date			



Ambulatory Blood Pressure Monitor



Ambulatory Blood Pressure

AMBULATORY BLOOD PRESSURE

Ambulatory blood pressure is a leading diagnostic tool. Rather than relying on one isolated reading, it compiles a series of measurements resulting in increased accuracy. Spacelabs offers a line of reusable and disposable cuffs, covers and pouches for use with the Spacelabs 90227, 90207 and 90217A APB Monitors.



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Reusable ABP Cuffs




UltraCheck Reusable Nylon ABP Cuffs

016-1320-02	Statcorp reusable ABP cuff	child	13 - 18 cm	1 Each
016-1826-02	Statcorp reusable ABP cuff	small adult	18 - 26 cm	1 Each
016-2635-02	Statcorp reusable ABP cuff	adult	26 - 35 cm	1 Each
016-3242-02	Statcorp reusable ABP cuff	large adult	32 - 42 cm	1 Each



Spacelabs Reusable Cloth ABP Cuffs (includes bladder, wrap, hose, fittings)

015-0118-05Q	Spacelabs reusable ABP cuff kit	child	12 - 20 cm	1 Each
015-0067-05Q	Spacelabs reusable ABP cuff kit	small adult	17 - 26 cm	1 Each
015-0068-06Q	Spacelabs reusable ABP cuff kit	adult	24 - 32 cm	1 Each
016-0077-06Q	Spacelabs reusable ABP cuff kit	large adult	32 - 42 cm	1 Each
016-0109-05Q	Spacelabs reusable ABP cuff kit	XL adult	38 - 50 cm	1 Each



Disposable Covers



Disposable Cuff Covers (for use with Spacelabs reusable ABP cuffs above)

706-0182-00	Universal cuff cover	Standard	50/Box
706-0183-00	Universal cuff cover	Large	50/Box

Disposable ABP Cuffs



Spacelabs Soft Fabric Disposable ABP Cuffs (no fitting)

015-1826-00	Spacelabs disposable ABP cuff	small adult	18 - 26 cm	20/Box
015-2635-00	Spacelabs disposable ABP cuff	adult	26 - 35 cm	20/Box
015-3242-00	Spacelabs disposable ABP cuff	large adult	32 - 42 cm	20/Box
015-4250-00	Spacelabs disposable ABP cuff	XL adult	42 - 50 cm	20/Box

Ambulatory Blood Pressure

AMBULATORY BLOOD PRESSURE

Ambulatory blood pressure is a leading diagnostic tool. Rather than relying on one isolated reading, it compiles a series of measurements resulting in increased accuracy. Spacelabs offers a line of reusable and disposable cuffs, covers and pouches for use with the Spacelabs 90227, 90207 and 90217A APB Monitors.



Order Code	Description	Size	Range	Qty
ZFT-QC8220	Quick release metal fittings, female bayonet (HP), 1/8 in. barb, For use with above single tube, disposable cuffs with no fitting.			each

ABP Accessories



Pouches - For use with ABP Monitor 90207

016-0410-01	Black polypropylene reusable pouch	1 Each
-------------	------------------------------------	--------



Pouches - For use with ABP Monitor 90217

016-0340-01	Black polypropylene reusable pouch	1 Each
-------------	------------------------------------	--------



Pouches - For use with ABP Monitor OnTrak 90227

015-0715-01	Black polypropylene reusable pouch	1 Each
-------------	------------------------------------	--------

Support

PK0001	1 inch with slide release buckle nylon webbing adjustable to max 59" inch length. For use as waist belt or shoulder strap.	1 Each
--------	--	--------



Miscellaneous ABP Accessories

712-0773-00	Adapter, converts ABP Luer Monitors (90207 and 90217) to quick disconnect	1 Each
-------------	---	--------



712-0060-00	QUICK CONNECTOR,1/8 TUBE,NIBP,CUFF	1 Each
-------------	------------------------------------	--------



Indicates discontinued part, available while supplies last

Holter and Event Monitoring



Holter and Event Monitoring Supplies


Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF



Evo



Aria



CardioCall

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Aria/evo Hook-Up Kits



50-0049	4 electrode hook-up kit with disposable evo chest pouch Includes: 4 clear tape electrodes, 4 prep pads, 1 prep scrub, 1 evo disposable pouch, 6 tape strips, 1 prep razor, 1 patient diary and pencil	Aria/evo	Each
---------	--	----------	------

Lead Wires



25452-525	4-Lead wire set, snap, 3 - 14" leads (orange/brown/red) and 1 - 7" lead (white)	Aria/evo	4LW/set
25452-535	4-Lead wire set, snap, 14" leads (red/white/orange/brown)	Aria/evo	4LW/set
25452-529	4-Lead wire set, snap, 40" leads (red/white/orange/brown)	Aria/evo	4LW/set
25452-519	7-Lead wire set, snap, 18" leads (white/black/red/green/orange/blue/brown)	423, 423A, 459, 461, 461A, 464, 483, 485	7 LW/set
19-7509	6-Lead wire set, snap, 24 in, AHA (requires 700-0316-00)	Lifecard CF	6LW/Set
23-0059	Patient cable, 1 channel, 2 electrodes (black/white), US	CardioCall	Each
23-0051	Patient cable, 2 channel, 3 electrodes (red/yellow/green)	CardioCall	Each


Electrodes and Prep supplies



392196-001	Adhesive tape remover pads, 100/box		100/Box
717-0083-00	Electrode, Transparent Tape, 5 per pouch X 10 Pouches		50/Box

Holter and Event Monitoring Supplies

Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF

Evo

Aria

CardioCall

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Lifecard CF Patient Cables



700-0316-00	3-channel, 6-lead, 22 in, AHA (requires 19-7509)	Lifecard CF	Each
700-0317-00	3-channel, 3-lead, snap, long	Lifecard CF	Each
700-0318-00	3-channel, 3-lead, snap, extra-long	Lifecard CF	Each
700-0319-00	3-channel, 3-lead, snap, short	Lifecard CF	Each
700-0320-00	CABLE,PX,LCF,2SNAP,CLIP,RD/WH.ROHS	Lifecard CF	Each
700-0321-00	CABLE,PX,LCF,2SNAP,SHORT,WO/CLIP,RD/WH.ROHS	Lifecard CF	Each
700-0322-00	CABLE,PX,LCF,7RECEPTACLE,AHA,WO/LEADWIRES,ROHS	Lifecard CF	Each
700-0323-00	2-channel, 4-lead, snap, long	Lifecard CF	Each
700-0324-00	2-channel, 4-lead, snap, short	Lifecard CF	Each

Lifecard 12 Patient Cables



700-0328-00	12-lead patient cable, EU	Lifecard 12	Each
700-0330-00	12-lead patient cable, long, EU	Lifecard 12	Each
700-0331-00	12-lead patient cable, long, US	Lifecard 12	Each

Pouches, Shoulder Straps & Belts



27249-7-401	Disposable pouch	Aria	50/Pack
50-0191	HOOK UP KIT WITH SINGLEUSE POUCH, CASE OF 20	Evo	20/Case
50-0042	Belt holster	Evo	Each
PK0001	Adjustable 1" nylon webbing w/ buckle, max 56" length	Evo	Each
46-0005	Disposable cord (use with 706-0189-00)	Evo	50m roll
50-0050	Disposable pouch	Evo	50/Pack
46-0408	Fabric pouch set, 2 pouches with laces, Cordlock	Lifecard CF	Each
46-0512	Disposable pouch	Lifecard CF	100/Pack



Holter and Event Monitoring Supplies


Spacelabs Healthcare combine state-of-the-art Holter recorders, fast accurate analyzers and 7-day ECG screening solutions. With our powerful information management platform Sentinel, clinicians the world over count on us for reliable, scalable and accurate Holter reports.



Lifecard CF



Evo



Aria



CardioCall

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Recorder Memory and Download



CFC90MB	Compact flash card memory, 90MByte	Lifecard CF	Each
CFC256MB	Compact flash card memory, 256MB	LC12	Each
019-1013-00	Compact Flashcard Reader	Lifecard CF	Each
019-1014-00	USB C to USB A adapter	Lifecard CF	Each

Bold font indicates new part numbers

Patient Diaries



073-0267-00	7 Patient Diary, English CE0123	Aria, Evo, Lifecard CF	Each
-------------	---------------------------------	---------------------------	------

Eclipse Mini and Eclipse Pro



Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.

Eclipse Pro and Mini Holter



Eclipse Pro



Eclipse Mini

Order Code	Description	Length/Style	Qty																																														
Eclipse Pro																																																	
Leadwire Sets																																																	
700-0378-00	3 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>M</td> <td>Yellow</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C5R</td> <td>Green</td> <td>205</td> <td>8</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	205	8	510	20	865	34	M	Yellow	255	10	560	22	915	36	C5R	Green	205	8	510	20	865	34	Each								
Snap	Color	Short			Medium		Long																																										
		mm	Inch	mm	Inch	mm	Inch																																										
C5	Red	205	8	510	20	865	34																																										
M	Yellow	255	10	560	22	915	36																																										
C5R	Green	205	8	510	20	865	34																																										
700-0391-00	3 Lead snap on with yoke Medium		Each																																														
700-0390-00	3 Lead snap on with yoke Long		Each																																														
700-0379-00	4 Lead snap on with yoke Short	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th rowspan="2">Color</th> <th colspan="2">Short</th> <th colspan="2">Medium</th> <th colspan="2">Long</th> </tr> <tr> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>C5</td> <td>Red</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> <tr> <td>C3</td> <td>Orange</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>M</td> <td>White</td> <td>255</td> <td>10</td> <td>560</td> <td>22</td> <td>915</td> <td>36</td> </tr> <tr> <td>C1</td> <td>Brown</td> <td>230</td> <td>9</td> <td>510</td> <td>20</td> <td>865</td> <td>34</td> </tr> </tbody> </table>	Snap	Color	Short		Medium		Long		mm	Inch	mm	Inch	mm	Inch	C5	Red	230	9	510	20	865	34	C3	Orange	255	10	560	22	915	36	M	White	255	10	560	22	915	36	C1	Brown	230	9	510	20	865	34	Each
Snap	Color	Short			Medium		Long																																										
		mm	Inch	mm	Inch	mm	Inch																																										
C5	Red	230	9	510	20	865	34																																										
C3	Orange	255	10	560	22	915	36																																										
M	White	255	10	560	22	915	36																																										
C1	Brown	230	9	510	20	865	34																																										
700-0394-00	4 Lead snap on with yoke Medium		Each																																														
700-0393-00	4 Lead snap on with yoke Long		Each																																														
700-0396-00	10 Lead snap on with yoke <i>IEC color coding generally used in European Union</i>	<table border="1"> <thead> <tr> <th rowspan="2">Snap</th> <th colspan="2">Length</th> </tr> <tr> <th>mm</th> <th>Inch</th> </tr> </thead> <tbody> <tr> <td>L</td> <td>610</td> <td>24</td> </tr> <tr> <td>F</td> <td>585</td> <td>23</td> </tr> <tr> <td>C6</td> <td>635</td> <td>25</td> </tr> <tr> <td>C5</td> <td>535</td> <td>21</td> </tr> </tbody> </table>	Snap	Length		mm	Inch	L	610	24	F	585	23	C6	635	25	C5	535	21	Each																													
Snap	Length																																																
	mm	Inch																																															
L	610	24																																															
F	585	23																																															
C6	635	25																																															
C5	535	21																																															
700-0377-00	10 Lead snap on with yoke <i>AHA color coding generally used in USA, outside of EU</i>	<table border="1"> <tbody> <tr> <td>C4</td> <td>460</td> <td>18</td> </tr> <tr> <td>C3</td> <td>355</td> <td>14</td> </tr> <tr> <td>C2</td> <td>330</td> <td>13</td> </tr> <tr> <td>C1</td> <td>330</td> <td>13</td> </tr> <tr> <td>N</td> <td>380</td> <td>15</td> </tr> <tr> <td>R</td> <td>610</td> <td>24</td> </tr> </tbody> </table>	C4	460	18	C3	355	14	C2	330	13	C1	330	13	N	380	15	R	610	24	Each																												
C4	460	18																																															
C3	355	14																																															
C2	330	13																																															
C1	330	13																																															
N	380	15																																															
R	610	24																																															
712-1052-01	ECLIPSE PRO ECG SENSORS <i>Not available for sale in USA</i> <i>24 month shelf life from date of manufacture</i>		Box of 10																																														
Accessories																																																	
706-0207-00	Eclipse Pro belt clip		Each																																														
073-0267-00	7 Patient Diary, English CE0123		Each																																														
392196-001	Adhesive tape remover pads, 100/box		100/Box																																														

Eclipse Pro is your versatile choice for screening, arrhythmia, pacemaker detection. A comfortable hook-up means your patient is more likely to be able to tolerate the recording for the desired length of time, and also avoid artifact and lead-off challenges. Eclipse Pro is available with 3, 4 or 10 lead options as well as a wearable patch. The Eclipse Mini is patch only.

Note: Patches may not be available in all regions at this time.


Eclipse Pro and Mini Holter




Eclipse Pro



Eclipse Mini

	Order Code	Description	Length/Style	Qty
Eclipse Mini				
	Eclipse Mini Supplies			
	712-1053-02	ECLIPSE MINI ECG SENSORS, DEFIB SAFE <i>24 month shelf life from date of manufacture</i>		Box of 10

ECG Stress Testing



Brand Name Stress Test Devices

ECG Stress Testing

Spacelabs offers supplies for the Suntech Tango M2 ABP system. Spacelabs has discontinued sale of these systems but continues to offer some accessories, while supplies last. Note: not available for sale into USA but is available for sale in the European Union and United Kingdom.



Suntech Tango M2

Order Code	Description	Dimensions	Qty
------------	-------------	------------	-----



Suntech Tango NIBP Supplies for Stress Test

010-2238-00	CUFF, ORBIT-K, SMALL ADULT		Each
010-2208-00	ECG PATIENT CABLE, IEC,10 LEAD CABLE, BANANA		Each
010-2210-00	SNAP / TAB ELECTRODE ADAPTERS, PK 10		Bag of 10
015-0711-00	Orbit cuff, Adult Plus, 33-44 cm		Each
015-0712-00	Orbit cuff, Large Adult, 39-46 cm		Each
98-0061-02	Orbit-K cuff, adult, 23-35 cm		Each
98-0002-50	Microphone pads		5/Pack


Indicates discontinued part, available while supplies last

CardioPulse ECG



CardioPulse Resting ECG










Compatible supplies and accessories for the portfolio of CardioPulse ECG 12-lead systems. Note: CardioPulse is only available for sale in the USA.



CardioPulse Prime



CardioPulse Go

	Order Code	Description	For Use With	Qty
Paper				
	010-2280-00	PAPER,THERM,Z-FOLD,210 MM,100SH	Prime	100 sheets
	010-2281-00	PAPER,THERMAL,Z-FOLD,114 MM,60SHX10PKS	Go	600 sheets
Lead Wires, Cables, Adapters				
	700-0386-00	CABLE,PATIENT,AHA,SNAP TYPE,CARDIOPULSE GO	Go	Each
	700-0387-00	CABLE,PATIENT,10-WIRE,AHA,BANANA PLUG,CP GO	Go	Each
	700-0388-00	CABLE,PATIENT,10-WIRE,AHA,CLIP,CP PRIME,1/BG	Prime	Each
	700-0389-00	CABLE,PATIENT,10-WIRE,AHA,BANANAPLUG,CP PRIME	Prime	Each
	700-0399-00	ADAPTER,D15,PATIENT CABLE COMPATIBLE,CP GO	Go	Each
Tab Electrodes, for use with electrode clips				
	58100	EasiClip, Skintact, tab electrode clips <i>For use with banana 4mm sized plugs</i>		Pack of 10
Barcode Scanner				
	010-2282-00	BARCODE SCANNER,LASER,INSTRUCTIONS,CARDIOPULSE		Each

Indicates discontinued part, available while supplies last

Resting ECG



CardioExpress resting ECG

CardioDirect resting ECG

CardioExpress SL6A, SL12A and SL18A provide fast and accurate 12-lead ECG acquisition. Advanced design and ergonomics enable the operator to quickly and easily control the unit to meet individual requirements. They provide storage direct to the units and connectivity via Sentinel, Cardiology Information Management System (CIMS).



CardioExpress



CardioDirect 12

Resting ECG systems
6, 12 and 18 leads

Order Code	Description	For Use With	Qty
ECG Patient Cables & Hook-up Sets			
MS1R-107584-A0	ECG cable, grabber connector, AHA	CardioExpress	Each
MS1R-107048-A0	ECG cable, 4MM, banana connector, AHA	CardioExpress	Each
010-2106-00	ECG cable, banana plugs, AHA	CardioExpress SL18A	Each
010-2107-00	ECG CABLE, SL18A, IEC, BANANA PLUGS	CardioExpress SL18A	Each
MS1R-106902-A0	ECG cable, 4MM, banana connector, IEC	CardioExpress	Each
25-0004	10 electrode patient cable, 4mm, EU	CD12	Each
25-0005	10 electrode patient cable, pinch, EU	CD12	Each
M15R-040172	CLIP/SNAP, BANANA SOCKET ADAPTERS, 10 PCS	CardioExpress	Each
M15R-040163	ADULT PRECORDIAL SUCTION ELECTRODES (6 PCS/SET)	CardioExpress	Each
M15R-040168	CHILD PRECORDIAL SUCTION ELECTRODES (6PCS/SET)	CardioExpress	Each
012-0907-00	ELECTRODE,ADULT LIMB CLAMP, QTY 4	CardioExpress	Each
58100	Easiclip, Skintact, 10/pk		10/Pack
AHA CardioExpress SL18A lead wires			
012-0964-00	LA ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0967-00	V3R ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0968-00	V1+V2+V3 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
012-0969-00	V4+V5+V6 ECG lead wire for SL-18A sampling box, AHA	CardioExpress SL18A	Each
ECG Patient Cables & Hook-up Sets			
IEC CardioExpress SL18A lead wires			
012-0975-00	R (RA) ECG Lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0976-00	L (LA) ECG LEADWIRE FORSL-18A SAMPLING BOX IEC	CardioExpress SL18A	Each
012-0982-00	C1+C2+C3 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0983-00	C4+C5+C6 ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0985-00	F (LL) ECG lead wire for SL-18A sampling box, IEC	CardioExpress SL18A	Each
012-0986-00	N (RL) ECG LEADWIRE FORSL-18A SAMPLING BOX IEC	CardioExpress SL18A	Each

Indicates discontinued part, available while supplies last

CardioExpress resting ECG

CardioDirect resting ECG

CardioExpress SL6A, SL12A and SL18A provide fast and accurate 12-lead ECG acquisition. Advanced design and ergonomics enable the operator to quickly and easily control the unit to meet individual requirements. They provide storage direct to the units and connectivity via Sentinel, Cardiology Information Management System (CIMS).

Resting ECG systems 6, 12 and 18 leads



CardioExpress



CardioDirect 12

Order Code	Description	For Use With	Qty
Miscellaneous ECG Accessories			
Recording Paper			
MS1R-107371-A0	Z-fold, 210 MM x 295 MM, 100 sheets	CardioExpress	Each
MS1R-32462	Z-fold, 210 MM x 295 MM, 200 sheets	CardioExpress	Each
MS1R-107560-A0	Z-fold, 110 MM x 140 MM, 144 sheets	CardioExpress	Each
MS1R-107451-A0	Z-fold, 215 MM x 280MM, 100 sheets	CardioExpress	Each
MS1R-32461	Roll, 210 MM x 30M	CardioExpress	Each
010-1897-00	RECORDING PAPER,ROLL,80MM X 20 M	CardioExpress	Each
Others			
010-2128-00	Carry case, SL12A	CardioExpress SL12A	Each
010-2123-00	USB Flash Drive	CardioExpress SL12A	Each
690-0275-00	Barcode scanner, one dimensional	CardioExpress	Each
690-0276-00	Barcode scanner, two dimensional	CardioExpress	Each

Indicates discontinued part, available while supplies last

Non-Medical



Non-Medical Accessories

Spacelabs carries a convenient assortment of Non-Medical supplies to help support your use of our monitoring systems.

Order Code Description With Qty



Note: Bumper is black. Yellow color used in image to highlight design and fit of the bumper.

Bumper Kits

010-2350-00	Aria Tele Bumpers	Aria Tele	5/Box
050-0816-00	Qube Mini Bumper Kit (black)	Qube Mini	Each



Cable Management Systems

016-1027-00	CORD CADDY, 4 PAIRS OF HOOKS	\$224.89	Each
016-1028-00	CORD CADDY, 6 PAIRS OF HOOKS	\$265.78	Each
016-1029-00	MOUNT,DOWNPOST,GCX,12 INCH (required to mount caddy)	\$65.52	Each

AV and Networking Cables

2287169	TRIPP LITE 10' HDMI M/M CABLE	\$19.14	Each
2195215	DP mini to DP cable, male/male, 10ft	\$26.14	Each
1973897	DP mini to DP cable, male/male, 6ft	\$20.22	Each
384343	USB 2.0 A to B Printer Cable, 10ft	\$11.84	Each
1561876	DP to DP cable, male/male, 10ft	\$28.82	Each
54403	DP to DP cable, male/male, 15ft	\$51.82	Each

Computer Supplies



5491141	NETGEAR GS308V3 8PRT SWITCH	\$43.10	Each
5615760	NETGEAR GS305 5PRT SWITCH	\$35.26	Each
3006282	LOGITECH MK270 WRLS KBD/MS SET	\$47.50	Each
010-2327-00	Poly Blackwire 5210 Monaural headset for use with SafeNSound	\$120.66	Each
6113257	USB-A port blockers (4) with reusable key (1)	\$56.46	Pack
6113260	USB-A port blockers (10)	\$66.80	Pack



Barcode Scanners

690-0275-00	Zebra LS2208 Barcode Scanner	\$250	Each
690-0276-00	Honeywell 1900GSR-2USB-C Barcode Scanner	\$400	Each

Non-Medical Accessories

Spacelabs carries a convenient assortment of Non-Medical supplies to help support your use of our monitoring systems.

Order Code Description With Qty



Displays and Mounts

210326	ERGOTRON DS100, DUAL DISP MOUNT	\$668.66	Each
6474205	NEC 43" LED-BACKLIT LCD DISPLAY	\$1,575.74	Each
6358610	NEC 50" LED-BACKLIT LCD DISPLAY	\$2,064.32	Each



Exergen TAT and Mounts: for use with Xprezzon, Qube, Qube Mini, DM4

010-2157-00	Exergen TAT 5000S, with coiled USB cable		Each
134201	Exergen TAT 5000S Wall mount, black		Each
134305	Exergen TAT 5000S locking wall mount, white		Each
134306	Exergen TAT 5000S locking wall mount with cap holder, white		Each

Universal Power Supplies and Network Switches



317698	APC,CS500,UPS,300WATT,500VA	\$273.40	Each
6064874	APC,BE850G2,UPS,450WATT,850VA	\$253.08	Each
5491141	NETGEAR GS308V3 8PRT SWITCH	\$43.10	Each
5615760	NETGEAR GS305 5PRT SWITCH	\$35.26	Each