

## Capnography & Multigas





(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

## NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.

**Order Code Description Length/Size Qty**



### NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

NomoLine LH Nasal CO2 Cannula		Shelf Life*		
704-3817-00	Adult	24 months	2 m/7 ft	25/Box
704-3819-00	Infant	24 months	2 m/7 ft	25/Box
704-3820-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4624-00	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3821-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box
NomoLine LH Nasal/Oral CO2 Cannula		Shelf Life*		
704-3824-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4628-00	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3825-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box
NomoLine LH Single Nasal Prong CO2 Cannula		Shelf Life*		
704-3826-00	Adult	24 months	2 m/7 ft	25/Box



### NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



NomoLine LH Airway Adapter Set		Shelf Life*		
704-3814-00	Adult/Pediatric	36 months	2 m/7 ft	25/Box
704-4620-00	Adult/Pediatric	36 months	4 m/13 ft	25/Box

### NomoLine LH (Low Humidity) Procedural Oxygen Masks



NomoLine LH Nasal CO2 Cannula		Shelf Life*		
704-5213-00 <sup>USA</sup>	Adult, Medium O2 Flow Rate	24 months	3 m/10 ft	30/Box
704-5214-00 <sup>USA</sup>	Pediatric, Medium O2 Flow	24 months	3 m/10 ft	30/Box
704-5215-00 <sup>USA</sup>	Adult, High O2 Flow Rate	24 months	3 m/10 ft	30/Box
704-5216-00 <sup>USA</sup>	Pediatric, High O2 Flow Rate	24 months	3 m/10 ft	30/Box

\* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

## NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.

**Order Code Description Length/Size Qty**



### NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients

NomoLine HH Nasal CO2 Cannula		Shelf Life*		
704-3830-00	Adult	24 months	2 m/7 ft	25/Box
704-3832-00	Infant	24 months	2 m/7 ft	25/Box
704-3833-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-3834-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box
704-4933-00 <sup>USA</sup>	Infant w/O2	24 months	4 m/13 ft	25/Box



NomoLine HH Nasal/Oral CO2 Cannula		Shelf Life*		
704-3837-00	Adult w/O2	24 months	2 m/7 ft	25/Box
704-4968-00 <sup>USA</sup>	Adult w/O2	24 months	4 m/13 ft	25/Box
704-3838-00	Pediatric w/O2	24 months	2 m/7 ft	25/Box




NomoLine HH Single Nasal Prong CO2 Cannula		Shelf Life*		
704-3839-00	Adult	24 months	2 m/7 ft	25/Box



### NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients

NomoLine HH Airway Adapter Set		Shelf Life*		
704-3827-00	Adult/Pediatric	36 months	2 m/7 ft	25/Box
704-3828-00	Adult/Pediatric	36 months	3 m/10 ft	25/Box
704-4989-00 <sup>USA</sup>	Adult/Pediatric	36 months	4 m/13 ft	25/Box
704-4367-00	Infant/Neonatal	36 months	2 m/7 ft	25/Box



\* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA

**Ready to order?**  **Click here to shop online**  
USA residents only

## CAPNOGRAPHY & MULTIGAS

### Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

#### Sidestream Monitoring Supplies

**Please note:** the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.



CO2 Nasal Cannula with luer		Shelf Life*		
-----------------------------	--	-------------	--	--

*	704-0176-00	Adult	60 months	2 m/7 ft	25/Box
---	-------------	-------	-----------	----------	--------



CO2 Nasal Cannula w/O2 and luer				
---------------------------------	--	--	--	--

*	704-0181-00	Adult	60 months	2 m/7 ft	25/Box
---	-------------	-------	-----------	----------	--------

CE	704-0181-01	Adult	36 months	2 m/7 ft	25/Box
----	-------------	-------	-----------	----------	--------

*	704-0180-00	Pediatric	60 months	2 m/7 ft	25/Box
---	-------------	-----------	-----------	----------	--------

**	704-0179-00	Infant	60 months	2 m/7 ft	25/Box
----	-------------	--------	-----------	----------	--------

**	704-0186-00	Adult	60 months	3 m/10 ft	25/Box
----	-------------	-------	-----------	-----------	--------



CO2 Oral/Nasal Cannula w/O2 and luer				
--------------------------------------	--	--	--	--

*	704-0183-00	Adult	60 months	2 m/7 ft	25/Box
---	-------------	-------	-----------	----------	--------

CE	704-0183-01	Adult	36 months	2 m/7 ft	25/Box
----	-------------	-------	-----------	----------	--------

**	704-0187-00	Adult	60 months	3 m/10 ft	25/Box
----	-------------	-------	-----------	-----------	--------

Indicates discontinued part, available while supplies last

\* - Shelf life is from date of manufacture, printed on product label as an expiration date

#### Sidestream Monitoring Supplies - Intubated patients



Airway adapters and lines		Shelf Life*		
---------------------------	--	-------------	--	--

**	704-0191-00	Adult Line w/ T-Adapter	60 months	2 m/7 ft	25/Box
----	-------------	-------------------------	-----------	----------	--------

**	704-0192-00	Adult Line w/ T-Adapter	60 months	3 m/10 ft	25/Box
----	-------------	-------------------------	-----------	-----------	--------



*	166-7121-00	sampling line	60 months	3 m/10 ft	50/Box
---	-------------	---------------	-----------	-----------	--------



*	704-0175-00	T-Adapter	60 months		25/Box
---	-------------	-----------	-----------	--	--------

\* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

\*\* Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

**Ready to order?  Click here to shop online**  
 USA residents only

## CAPNOGRAPHY & MULTIGAS

### Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



(92516) Capno Pod    (92517) Capnography Module    (92518) Multi-Gas Module

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

#### Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients.

Featuring Masimo IRMA mainstream technology



#### Sensor & Airway Adapters - required for mainstream monitoring








010-1980-00	Masimo IRMA mainstream sensor, reusable		1 Each
704-0173-00	Masimo IRMA mainstream airway adapter, single use <i>36 month shelf life from date of manufacture</i>	Adult/Pediatric	25/Box
704-0174-00	Masimo IRMA mainstream airway adapter, single use <i>36 month shelf life from date of manufacture</i>	Infant	10/Box

**Capnography Accessories**

The legacy dual capnography module 91517 is now obsolete. However, Spacelabs will continue to carry TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module, while these parts remain reasonably available from the manufacturer.



(91517)  
Dual Capnography Module


Order Code	Description	Size	Qty
<b>Sidestream, Non-intubated, single patient use</b>			
		<b>CO<sub>2</sub> Nasal Cannula</b>	
		<b>Shelf Life<sup>^</sup></b>	
* 704-0010-00	Adult	36 months	10/Box
* 704-0012-00	Pediatric	36 months	10/Box
		<b>CO<sub>2</sub> Nasal Cannula w/O<sub>2</sub> Delivery</b>	
* 704-0014-00	Adult	36 months	10/Box
CE 704-0014-10	Adult, (for sale into EU/UK only)	36 months	10/Box
* 704-0016-00	Pediatric	36 months	10/Box
* 704-0015-00	Infant/Neonatal	36 months	10/Box
		<b>CO<sub>2</sub> Oral/Nasal Cannula</b>	
* 704-0017-00	Adult	36 months	10/Box
		<b>CO<sub>2</sub> Oral/Nasal Cannula w/O<sub>2</sub> Delivery</b>	
* 704-0020-00	Adult	36 months	10/Box
* 704-0022-00	Pediatric	36 months	10/Box
<b>Sidestream, Intubated, single patient use</b>			
		<b>Sidestream - short term up to 12 hours</b>	
* 704-0023-00	Airway adapter + sampling line, ET tube > 4.0 mm	36 months	10/Box
CE 704-0023-10	Airway adapter + sampling line, ET tube > 4.0 mm	36 months	10/Box
* 704-0024-00	Airway adapter + sampling line, Pediatric/Infant sized	36 months	10/Box
* 704-0027-00	CO <sub>2</sub> sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box
CE 704-0027-10	CO <sub>2</sub> sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box

## Capnography Accessories

The legacy dual capnography module 91517 is now obsolete. However, Spacelabs will continue to carry TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module, while these parts remain reasonably available from the manufacturer.



(91517)  
Dual Capnography Module

Order Code	Description	Size	Qty
<b>Sidestream - long term up to 120 hours</b>		<b>Shelf Life<sup>^</sup></b>	
 *	704-0025-00 Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	36 months	10/Box
CE	704-0025-10 Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	36 months	10/Box
*	704-0026-00 Airway adapter + sampling line with Nafion, Infant	36 months	10/Box
*	704-0028-00 CO <sub>2</sub> sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box
CE	704-0028-10 CO <sub>2</sub> sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>	36 months	10/Box

### Mainstream, for Intubated Patients Only



#### Mainstream Airway Adapters, for use with CO<sub>2</sub> Sensor 704-0001-00



*	704-0005-00 CO <sub>2</sub> airway adapter, single use, Adult	36 months	24/Box
*	704-0002-00 CO <sub>2</sub> airway adapter, reusable, Adult	36 months	1 Each

**Indicates discontinued part, available while supplies last**  
<sup>^</sup> - Shelf life is from date of manufacture, printed on product label as an expiration date  
 \* Product not available for sale into CE Mark Countries (EU)  
**CE** product is available for sale into CE Mark Countries (EU)  
 \*\* Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

**Ready to order?**  **Click here to shop online**  
USA residents only

## MULTIGAS ANALYZER Modules 90518 & 91518

### Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers  
(90518 & 91518).

These two legacy multi-gas modules are obsolete. However, Spacelabs continues to provide supplies that are available.



(90518)



(91518)

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

#### Multigas Supplies & Accessories

##### Anesthetic Agent Samples Lines



015-0128-00	SAMPLE LINE,MALE TO MALE,LUER,SPU,7'X.04",25/BX		25/Box
015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line:0.048 inches ID tubing (90518/91518) <i>36 month shelf life from date of manufacture</i>	305 cm/10 ft.	50/Box



162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 0.048 in ID tubing (90518/91518) <i>36 month shelf life from date of manufacture</i>	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



##### Water Traps, Reusable

717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)		10/Pack
717-0448-00	D-FEND water trap cartridge, Green. Replace at least once every 24 hours or when occluded (91518)		10/Pack

*Note: manufacturer (GE) has discontinued these parts. Available while supplies last.*