


Cardiac Output



Ready to order?  **Click here to shop online**
USA residents only

Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Cardiac Output Cables & Probes

Cables



700-0027-00	TruLink Cardiac Output Adapter Cable. Connects CO cable 30 cm/12 in. to 91496 command module.		1 Each
015-0739-00	Cardiac Output cable Edwards compatible connects with 690-0033-00 temperature probe. Replaces 306655-003. The 700-0027-00 adapter not required.		1 Each
306655-004	Cardiac Output cable Datex/Ohmeda/Abbott compatible connects with 690-0023-00 or 690-0024-00 temperature probe. 700-0027-00 adapter required.		1 Each

Bold font indicates a new part number

Indicates discontinued part, available while supplies last.

Injectate Probes



690-0023-00	In-line injectate temperature probe. Datex/Ohmeda/Abbott connector style. Connects to 306655-004 cable and 690-0293-00 injectate kit integrated thermistor.		1 Each
690-0024-00	Ice bath temperature probe. Datex/Ohmeda/Abbott connector style. Connects to 306655-004 cable and needle end is submerged in an ice bath.		1 Each
690-0033-00	In-line injectate temperature probe. Edwards connector style. Connects to 306655-003 cable and 690-0293-00 injectate kit integrated thermistor.		1 Each

Cardiac Output Injectate Kits

Closed Loop Injectate Kits

690-0293-00



SP5045 Critikit In-Line Sensor, Control Syringe closed loop injectate kit
Packaged sterile with 36 month shelf life

5/Box

Cardiac Output Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for cardiac output. These supplies facilitate the thermodilution method. Spacelabs offers adapters, cables, reference and inline injectate probes and kits compatible with BD or Edwards Lifesciences injectate probes and most models of thermodilution catheters.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

91496-C Module, Cardiac Output Connection guide

