


# Electrocardiography

**AAMI COLOR CODING**



**ECG Combiner Lead Wires and Cables**

Our ECG Combiner System features a quick-yoke connection between cables and lead wires which allows a quickly connect/disconnect between monitors and telemetry. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

**ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter**



**Disposable, DIN Shielded, Pinch, AAMI**

700-0006-50	3-Lead	Pinch	127cm /50 in	25 Sets
* 700-0006-32	5 Lead	Pinch	127cm /50 in	25 Sets
* 700-0006-33	5 Lead (V-Lead)	Pinch	127cm /50 in	25 Sets



**Disposable, DIN Shielded, Snap, AAMI**

<b>700-0006-61</b>	<b>5 Lead</b>	<b>Snap</b>	<b>127cm /50 in</b>	<b>25 Sets</b>
--------------------	---------------	-------------	---------------------	----------------



**Reusable, DIN Shielded, Pinch, AAMI**

* 700-0006-35	3 Lead	Pinch	74cm /29 in	1 Set
* 700-0006-37	5 Lead	Pinch	74cm /29 in	1 Set
* 700-0006-38	5 Lead (V-Lead)	Pinch	74cm /29 in	1 Set
700-0010-00	3 Lead	Pinch	130cm/51 in	1 Set
* 700-0010-01	5 Lead	Pinch	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Pinch	130cm/51 in	1 Set




**Reusable, DIN Shielded, Snap, AAMI**

* 700-0007-35	3 Lead	Snap	74cm /29 in	1 Set
* 700-0007-37	5 Lead	Snap	74cm /29 in	1 Set

\* Product not available for sale into CE Mark Countries (EU)

## ECG Combiner Lead Wires and Cables

Our ECG Combiner System features a quick-yoke connection between cables and lead wires which allows a quickly connect/disconnect between monitors and telemetry. System is validated for use with Spacelabs command modules 90496/91496 and Aria/Ultraview Telemetry.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

### ECG Combiner Cables



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496				
* 700-0008-55	3 Lead	Spacelabs	305 cm/10 ft	1 Cable
* 700-0008-57	5 Lead	Spacelabs	305 cm/10 ft	1 Cable
700-0008-32	5 Lead	Spacelabs	457 cm/15 ft	1 Cable
* 700-0008-56	10 Lead (5+5)	Spacelabs	305 cm/10 ft	1 Cable

\* Product not available for sale into CE Mark Countries (EU)

### ECG Combiner Adapters



Spacelabs Telemetry, AAMI			
700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)	1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)	1 Adapter

### ECG Combiner Cables for Standardization




Other manufacturers				
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable

Indicates discontinued part, available while supplies last

**ECG Lead Wire & Cables - Shielded**

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



**Order Code      Description      Length/Size      Qty**

**ECG Lead Wire Sets, Shielded, Reusable, AAMI**



**Snap**

**Snap Sets**

700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set



**Pinch**

**Pinch Sets**

700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set

**ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI**



**24" Snap Individual Color**

012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

**ECG Lead Wire & Cables - Shielded**

TruLink ECG lead wires are designed and validated for use with Spacelabs monitoring systems. Ten foot reusable monitor cables are available in 3, 5, or 10 lead configurations. The reusable shielded TruLink lead wires feature individual color-coded leads with a 2-prong connection into the monitor cable. They are available in a variety of lengths with either a snap or pinch electrode interface. For use with Spacelabs command modules 90496/91496.



**Length/Size Qty**


**ECG Cables, Shielded, AAMI**



**Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496**

Order Code	Description	Length/Size	Qty
700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

**Neonatal ECG Supplies & Accessories**



**DIN Unshielded, Neonatal ECG Cable**


700-0008-12	3 Lead	Neonatal	1 Cable
-------------	--------	----------	---------

**DIN Unshielded, Neonatal Electrodes.**

717-0080-00	Pre-wired neonatal electrodes, Pre-gelled, silver/silver chloride, single use <i>24 month shelf life from date of manufacture</i>	3-Lead	100 Sets/Box
-------------	--	--------	--------------

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are available in 3 and 5-lead sets. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



#### Snap Sets

015-0732-00	3 Lead	90 cm/35 in	1 Set
015-0736-00	5 Lead	90 cm/35 in	1 Set



#### Pinch Sets

015-0731-00	3 Lead	90 cm/35 in	1 Set
<b>012-0284-11</b>	<b>3 Lead</b>	<b>61 cm/24 in</b>	<b>1 set</b>
015-0734-00	5 Lead	90 cm/35 in	1 Set

*Indicates discontinued part, available while supplies last*

### Telemetry Transmitter Adapters



#### For use with DIN Unshielded Lead Wires


706-0191-00	Aria Telemetry Grouper (96281)	5/Bag
-------------	--------------------------------	-------

**Ready to order?**  **Click here to shop online**  
 USA residents only

**ECG LEAD WIRES & TELEMETRY  
 AAMI - DIN Unshielded**

**DIN Unshielded Leadwire & Telemetry**

TruLink ECG Din unshielded reusable lead wires are available in 3 and 5-lead sets. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

**Telemetry Accessories**

**Protective Solutions and Electrodes**




010-2350-00	Aria Tele Bumper		5/Box
-------------	------------------	--	-------



015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
-------------	-----------------------------------	--	---------

## DIN Unshielded Leadwire & Telemetry

TruLink ECG Din unshielded reusable lead wires are available in 3 and 5-lead sets. They can be ordered with a snap or pinch electrode interface. A grouper adapter is also available to help manage leads for a quick connection to the transmitter. Validated for use with models 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### SpO<sub>2</sub> Cable & Sensors

**Note: Only the below genuine Spacelabs and Nellcor branded SpO<sub>2</sub> sensors are compatible with Spacelabs Telemetry systems**

#### SpO<sub>2</sub> Cable



700-0014-00	TruLink SpO2 Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO2 sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

#### Reusable SpO<sub>2</sub> Sensors



015-0660-00	TruLink SpO <sub>2</sub> finger sensor	Adult >30 kg/66 lbs.	1 Each
690-0003-01	Nellcor DS-100A, finger SpO2 sensor <i>Available for sale with monitor purchase</i>	Adult >30 kg/66 lbs.	1 Each

#### Disposable SpO<sub>2</sub> Sensors



015-0662-00	TruLink SpO <sub>2</sub> foam sensor	Adult > 30 kg/66 lbs.	24/Box
-------------	--------------------------------------	--------------------------	--------



015-0663-00	TruLink SpO <sub>2</sub> vinyl sensor	Adult >30 kg/66 lbs.	24/Box
-------------	---------------------------------------	-------------------------	--------



015-0664-00	TruLink SpO <sub>2</sub> foam sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	--------------------------------------	---------------------------------------	--------



015-0665-00	TruLink SpO <sub>2</sub> vinyl sensor	Pediatric 10-50 kg/ 22-110 lbs.	24/Box
-------------	---------------------------------------	---------------------------------------	--------



015-0667-00	TruLink SpO <sub>2</sub> cloth sensor	Neonatal <3 kg/7 lbs.	24/Box
-------------	---------------------------------------	--------------------------	--------

**Indicates discontinued part, available while supplies last**