

C50

VITAL SIGNS MONITOR





Spacelabs C50 Patient Monitoring System

Spacelabs will discontinue sale of the C50 monitor on January 31, 2026. Supplies and accessories are available while supplies last.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

	Order Code	Description	Dimensions	Usage	Qty
ECG leadwire systems					
	040-000908-00	CABLE,INTEGRATED 3-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000909-00	CABLE,INTEGRATED 5-WIRECLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000910-00	CABLE,INTEGRATED 10-WIRE CLIP,AHA,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000911-00	CABLE,INTEGRATED 3-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
	040-000912-00	CABLE,INTEGRATED 5-WIRECLIP,IEC,C50	3 meters with 4 meter limb leads	Reusable	1 Each
Pulse Oximetry Supplies					
	040-000869-00	CABLE,INTEGRATED ADT FINGER CLIP SENSOR,C50	3 meters	Reusable	1 Each
	040-000769-00	CABLE,SPO2,12PIN TO DB9F,C50 <i>Sensors sold seperately</i>	3 meters	Reusable	1 Each
	040-000730-00	SENSOR,SPO2,ADT/NEO/PED,WRAP, 3',C50 <i>Requires the 040-000769-00 cable</i>	1 meter	Reusable	1 Each








Spacelabs C50 Patient Monitoring System

Spacelabs will discontinue sale of the C50 monitor on January 31, 2026. Supplies and accessories are available while supplies last.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

	Order Code	Description	Dimensions	Usage	Qty
	Non-Invasive Blood Pressure supplies				
	040-000808-00	ADULT NIBP TUBE (WHITE)2M	2 meters	Reusable	1 Each
	Reusable or Disposable Cuffs	See NIBP section under TruLink Reusable or Disposable cuffs. These are compatible with the C50. Requires use of the 040-00808-00 air hose (above).			
	714-0086-00	NIBP,HOSE ASSY,NEONATE,C50 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	Reusable	1 Each
	Capnography Supplies: requires external Phasein ISA module, 099-000007-00				
	040-000400-00	CABLE,CO2/AG,C50 Connects the Phasein ISA module to the C50	1 meter	Reusable	1 Each
	For airway adapter, cannula and sampling line sets, please see the capnography supply section under Nomoline SPU supplies. These are compatible with the C50.				
   	Temperature sensors				
	040-000246-00	PROBE,TEMP,ADT,SKIN,10K,C50	3 meters	Reusable	1 Each
	040-000387-00	PROBE,TEMP,NEO SKIN,10K,C50	3 meters	Reusable	1 Each
	040-000385-00	PROBE,TEMP,ADT,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each
	040-000386-00	PROBE,TEMP,NEO,RECTAL/ESOPHAGEAL,10K,C50	3 meters	Reusable	1 Each




Spacelabs C50 Patient Monitoring System

Spacelabs will discontinue sale of the C50 monitor on January 31, 2026. Supplies and accessories are available while supplies last.

Note: the C50 system, including its supplies and accessories, may not be available for sale in all regions, including the USA.



C50

Order Code	Description	Dimensions	Usage	Qty
Invasive Blood Pressure Cables and Transducer kits				
 040-000454-00	INTERFACE CABLE,IBP,EDWARD CONNCTR,C50	 3 meters	Reusable	1 Each
Recording Paper				
 040-000116-00	RECORDER PAPER,THERMAL,C50	50mm X 20 meters		1 roll Each