

Elance

VITAL SIGNS MONITOR




Elace Vital Signs Monitor

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below
- TruLink ECG cables and leadwire systems are compatible



**Elace Vital Signs Monitor
93300**

Order Code	Description	Qty
Supplies unique to the Elance monitor		
ECG Lead Wires & Cables (can also use standard TruLink ECG cables and leadwire systems)		
ECG Integrated Cables and Lead Wires		
	700-0008-23 5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded 61cm/10in 305 cm/10 ft	1 Each
	700-0008-24 3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded 61cm/10in 305 cm/10 ft	1 Each
Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)		
SpO2 Cable		
	700-0287-00 SpO2 adapter cable 305 cm/10 ft	1 Each
LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)		
LNCS Adapter Cable		
	690-0256-00 SpO2 Adapter Cable, Masimo, LNCS 305 cm/10 ft	1 Each
Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)		
Cable		
	700-0290-00 SpO2 Adapter Cable, Nellcor OxiMax 305cm/10ft	1 Each
Temperature		
Temperature Probe		
	690-0261-01 Skin/Surface Temp probe, reusable, YSI-400 305cm/10ft	1 Each
Invasive Blood Pressure (See IBP section for compatible sensors and transducers)		
Interface Cable		
	690-0021-00 Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	1 Each
	Indicates discontinued part, available while supplies last	