


Non-Invasive Blood Pressure Cuffs



TruLink® Nylon Reusable

TruLink reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. TruLink cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, TruLink cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. TruLink cuffs are naturally latex and PVC-free.



**Range/
Length** **Qty**

Order Code	Description	Size	Range/ Length	Qty
TruLink Nylon Reusable BP Cuffs				
Single Tube				
<i>For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors</i>				
715-1111-10	Single tube, HP	Infant	8 - 14 cm	1 Each
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each
Single Tube Hose				
714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



TruLink® Soft Disposable

TruLink fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

TruLink Soft Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors

715-1113-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Single Tube Hose



714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each



TruLink® Vinyl Disposable




TruLink vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. TruLink cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, TruLink cuffs provide a wide range of options for correct readings.



	Order Code	Description	Size	Range/ Length	Qty
TruLink Vinyl Disposable BP Cuffs					
 Bayonet	Single Tube				
	<i>For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors</i>				
	715-1112-10	Single tube, HP	Infant	8 - 14 cm	5/Bag
	715-1122-10	Single tube, HP	Child	13 - 20 cm	5/Bag
	715-1132-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	715-1142-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
	715-1152-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	715-1162-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	715-1172-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
	715-1182-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag
	Single Tube Hose				
	714-0018-02	Single tube, quick connect	Adult	275 cm/9 ft	1 Each
	714-0060-02	Single tube, quick connect	Adult	366 cm/12 ft	1 Each

Blood Pressure Cuff Accessories





Spacelabs non-invasive blood pressure cuff accessories are compatible with all Spacelabs and TruLink branded non-invasive blood pressure cuffs. Choose from a variety of hoses, adapters and conversion accessories which offer ease-of-use and flexibility.

	Order Code	Description	Size	Length	Qty
Non-Invasive Blood Pressure Hose					
	Single Tube Hose				
	For use with Command Modules 90496/91496 and Ultraview DM3 (91330) and DM4 (91331) Monitors				
	714-0018-02	90496/91496/91330	Adult	275 cm/9 ft	1 Each
	714-0060-02	90496/91496/91330	Adult	366 cm/12 ft	1 Each
	714-0019-02	90496/91496	Neonatal	275 cm/9 ft	1 Each
	Double Tube to Single Tube Hose				
	For use with Command Modules 90496/91496				
	714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each
Hose Connectors & Adapters					
	Slip Lock Fitting				
	Fitting for Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.				
	714-0038-00	3mm ID Adult NIBP hoses			1 Each
	714-0046-00	5mm ID Adult NIBP hoses			1 Each
	714-0083-00	MQML TO SINGLE FEMALE BAYONET Y ADAPTER			1 Each
	714-0084-00	DUAL MQ TO SINGLE FEMALE BAYONET Y ADAPTER			1 Each
Indicates discontinued part, available while supplies last					

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.






	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
 Plastic Bayonet	Single Tube				
	US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
 Metal Bayonet	US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
	US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
	US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
	US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
	US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
	US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
	US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
	US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
 Threaded Screw	US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag
	US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
	US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
	US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
	US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
 Female Subminiature Twist Lock	US3544DM-05	Single tube, DM	Large Adult Long	32 - 42 cm	5/Bag
	US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
	US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
	US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
	US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
	US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag

Indicates discontinued part, available while supplies last

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.








	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
	Single Tube				
	US0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
	US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
	US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
Twist Lock					
	US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
No Fitting					
	US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
	US1320MB	Single tube, MB w/gauge	Child	13 - 20 cm	1 Each
	US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
	US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each
Manual Bulb	Indicates discontinued part, available while supplies last				

UltraCheck® Nylon Reusable

UltraCheck reusable nylon cuffs provide an economical and durable solution for your NIBP measurement needs. UltraCheck cuffs feature an anti-microbial coating, fully cleanable nylon construction, and are durable for multi-patient extended use. Available in 8 sizes, UltraCheck cuffs provide color-coded sizing and markings to assist you in selecting the right size for you patient. UltraCheck cuffs are naturally latex and PVC-free.



	Order Code	Description	Size	Length	Qty
UltraCheck Nylon Reusable BP Cuffs					
 Threaded Screw  Twist Lock  Male & Female Mated  No Fitting  Manual Bulb	Double Tube				
	UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
	UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
	UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
	UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
	UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
	UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
	UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag
	UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
	UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
	UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
	UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
	UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
	UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
	UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag

Indicates discontinued part, available while supplies last

SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



**Range/
Length**

Qty

Order Code

Description

Size

SoftCheck Fabric Disposable BP Cuffs

Single Tube

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box

Indicates discontinued part, available while supplies last



Plastic Bayonet







**Female
Subminiature**

SoftCheck® Fabric Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.





	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Fabric Single Patient Multiple Use BP Cuffs					
 No Fitting	Single Tube				
	ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
	ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
	ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
 Threaded Screw	Double Tube				
	DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box
 Twist Lock	DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box
 No fitting	DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag

Indicates discontinued part, available while supplies last

SoftCheck® Vinyl Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 8 sizes, SoftCheck cuffs provide a wide range of options for correct readings.



	Order Code	Description	Size	Range/ Length	Qty
 Plastic Bayonet	Single Tube				
	VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
	VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
	VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
	VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
	VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
	VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
	VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
	VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
 Male Twist Lock	VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
	VST0814MQML5-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
	VST1826MQML5-20	Single tube, MQ	Small Adult	18 - 26 cm	20/Box
	VST2635MQML-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
	VST2635MQML5-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
	VST4250MQML5-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box

Indicates discontinued part, available while supplies last

Ready to order?  **Click here to shop online**
USA residents only

BLOOD PRESSURE CUFFS UltraCheck® & SoftCheck® Curve

UltraCheck® Curve Nylon Reusable SoftCheck® Curve Fabric Disposable

Your patients come in all shapes and sizes. That's why we offer the Curve cuff for extra-large adults. Its conical-shaped design provides a better fit for the contours of your patient's arm. This can increase patient comfort while providing a more accurate reading.



Order Code	Description	Size	Range/ Length	Qty
------------	-------------	------	------------------	-----

UltraCheck Curve Nylon Reusable BP Cuff

Single Tube

BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MB	Single tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

Double Tube

BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MB	Double tube, manual bulb	Extra Large Adult	38 - 54 cm	1 Each

SoftCheck Curve Fabric Disposable BP Cuff

Single Tube

BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Box
BRST3854HP-05	Single tube, plastic HP	Extra Large Adult	38 - 54 cm	5/Box


Indicates discontinued part, available while supplies last



Plastic
Bayonet (HP)



Metal Bayonet
(HP)



Female Screw
Threaded (DM)



Female Locking
(MQ)



Male Locking
(MQML)



Mated Locking
(MQMF)






Manual Bulb
w/Gauge
(MB)

SoftCheck® Fabric Neonatal Disposable

SoftCheck fabric cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a super-soft Sontara fabric which is comfortable, helps reduce perspiration under the cuff, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.









	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Fabric Neonatal Disposable BP Cuffs					
 Micro-taper slip-fit	Single Tube				
	NN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	NN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	NN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	NN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	NN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
 Male Luer	NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Plastic Bayonet	NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
	NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
	NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Indicates discontinued part, available while supplies last

SoftCheck® Vinyl Neonatal Disposable

SoftCheck vinyl cuffs provide an economical and patient friendly solution for your NIBP measurement needs. SoftCheck cuffs feature a soft vinyl construction which is comfortable, cleanable, and is durable for long patient stays. The single patient - multiple use application offers a cost effective solution for infection control. Available in 5 different sizes for accurate readings.



	Order Code	Description	Size	Range/ Length	Qty
SoftCheck Vinyl Neonatal Disposable BP Cuffs					
 Micro-taper slip-fit	Single Tube				
	VNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box
 SnapQuik	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 Male Luer	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Plastic Bayonet	VNN1ST-HP-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN5ST-HP-10	Single tube, MQ	Size 4	8-15 cm	10/Box
 Female Twist Lock	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box
 Male Luer	Double Tube				
	VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
	VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
	VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
	VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
	VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box