

Elance

VITAL SIGNS MONITOR



Ready to order?  [Click here to shop online](#)
USA residents only

VITAL SIGNS MONITOR
Ultraview DM4



91331

Ultraview DM4 supports measurement of SpO2, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.
The DM4 uses TruLink branded NIBP cuffs (single tube, bayonet fitting) and air hose (714-0018-00)

Order Code	Description	Length/Size	Qty
Supplies unique to the DM4 monitor			
RD SET Adapter Cable			
 010-2193-00	DM4 Masimo RD SET SpO2 cable compatible with Masimo RD SET sensors, see section on RD SET sensors	305 cm/10 ft	1 Each
Filac			
 010-2174-00	Filac dock and insulation chamber		1 Each
Exergen			
 010-2157-00	Temporal Artery Thermometer		1 Each
NeoNatal NIBP supplies			
 714-0087-00	NIBP,HOSE ASSY,NEONATE,DM4 Use with cuffs NN1ST-10, NN2ST-10, NN3ST-10, NN4ST-10, NN5ST-10, VNN1ST-10, VNN2ST-10, VNN3ST-10, VNN4ST-10, VNN5ST-10	2.5 meters	1 Each

Elace Vital Signs Monitor

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below
- TruLink ECG cables and leadwire systems are compatible



**Elace Vital Signs Monitor
93300**

Order Code	Description	Qty
Supplies unique to the Elance monitor		
ECG Lead Wires & Cables (can also use standard TruLink ECG cables and leadwire systems)		
ECG Integrated Cables and Lead Wires		
	700-0008-23 5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
	700-0008-24 3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)		
	SpO2 Cable	
700-0287-00	SpO2 adapter cable	305 cm/10 ft 1 Each
LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)		
	LNCS Adapter Cable	
690-0256-00	SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft 1 Each
Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)		
	Cable	
700-0290-00	SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft 1 Each
Invasive Blood Pressure (See IBP section for compatible sensors and transducers)		
	Interface Cable	
690-0021-00	Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	366 cm/12 ft 1 Each
	Indicates discontinued part, available while supplies last	