








## Invasive Blood Pressure



## Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers.



	Order Code	Description		Length/Size	Qty
<b>Integrated Adapter Cables (connects directly to transducer)</b>					
	<b>Edward Lifesciences</b>				
	015-0727-00	Single pressure integrated cable		400 cm/13 ft	1 Each
	015-0728-00	Dual pressure integrated cable		400 cm/13 ft	1 Each
	<b>Abbott Laboratories/Hospira/ICU Medical</b>				
	015-0729-00	Single pressure integrated cable		400 cm/13 ft	1 Each
	015-0730-00	Dual pressure integrated cable		400 cm/13 ft	1 Each
<b>Interface Cable &amp; Adapters</b>					
	<b>Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)</b>				
	690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00		366 cm/12 ft	1 Each
	<b>Interface Adapter- Transducer interface cable required from transducer company</b>				
	700-0791-00	Single pressure cable		122 cm/4 ft	1 Each
	700-0028-00	Dual pressure cable		122 cm/4 ft	1 Each
<b>Transducer Kits &amp; Accessories</b>					
	<b>Pressure Transducer Kits, Single use</b>				
	690-0008-00	Transducer Kit, Patient Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>		152.5 cm/60 in	5/Box
	690-0009-00	Transducer Kit, Pole Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>		152.5 cm/60 in	5/Box

Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers.

