





## Temperature



## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.












	Order Code	Description	Length/Size	Qty
<b>Monitor Adapter Cables</b>				
	<b>Adapters</b>			
	700-0031-00	TruLink Dual Temperature Cable. *Required for all reusable probes and as directed for single-patient-use disposable system	30 cm/12 in.	1 Each
<b>Reusable Probes/Sensors</b>				
	<b>Skin Surface</b>			
	690-0295-00	YSI 400 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
	690-0296-00	YSI 700 series compatible skin surface sensor temperature probe, adult size	305 cm/10 ft	1 Each
	690-0297-00	YSI 700 series compatible skin surface sensor temperature probe, pediatric size	305 cm/10 ft	1 Each
	<b>Esophageal / Rectal</b>			
	690-0294-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
	690-0298-00	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
	690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
	118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
Indicates discontinued part, available while supplies last				

## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



	Order Code	Description	Length/Size	Qty	
Single-Patient use temperature system					
	Extension cables for single use probes				
	700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors, adapter required.	305 cm/10 ft	1 Each	
	700-0401-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors. Direct connection to Command module, no adapter required.	305 cm/10 ft	1 Each	
	Skin				
	701-4000-00	Skin Temperature Sensor	Sterile Yes	Shelf Life* 36 months	20/Box
	Esophageal/ Stethoscope				
	701-4009-00	9 French, YSI 400 series compat	Yes	36 months	20/Box
	701-4018-00	18 French, YSI 400 series compa	Yes	36 months	20/Box
	Esophageal/ Rectal				
	701-4109-00	9 French, YSI 400 series compat	Yes	36 months	20/Box
	701-4112-00	12 French, YSI 400 series compa	Yes	36 months	20/Box
	Tympanic temperature sensor				
	701-4300-00	Adult, YSI 400 series compatible	Yes	36 months	20/Box
* - Shelf life is from date of manufacture, printed on product label as an expiration date					
	Exergen				
	For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) & DM4 (91331) Monitors				
	010-2157-00	Temporal Artery Thermometer		1 Each	


## Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.








Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

### Reusable Temperature Probes

				
690-0294-00 YSI400 Adult Esophageal/ Rectal	690-0298-00 YSI700 Adult Esophageal/ Rectal	690-0295-00 YSI400 Adult Skin	690-0296-00 YSI700 Adult Skin	690-0297-00 YSI700 Pediatric Skin

### Single-Patient Use Temperature Probes

		
w/ Stethoscope 9Fr: 701-4009-00 12Fr: 701-4012-00 18Fr: 701-4018-00	Esophageal/Rectal 9Fr: 701-4109-00 12Fr: 701-4112-00	w/ Foley Catheter 14Fr: 701-4214-00 16Fr: 701-4216-00
		
Skin 701-4000-00	Tympanic 701-4300-00	

