

## Capnography & Multigas





## NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

**Order Code Description Length/Size Qty**



### NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients



#### NomoLine LH Nasal CO2 Cannula

#### Shelf Life\*

|             |                |           |           |        |
|-------------|----------------|-----------|-----------|--------|
| 704-3817-00 | Adult          | 24 months | 2 m/7 ft  | 25/Box |
| 704-3819-00 | Infant         | 24 months | 2 m/7 ft  | 25/Box |
| 704-3820-00 | Adult w/O2     | 24 months | 2 m/7 ft  | 25/Box |
| 704-4624-00 | Adult w/O2     | 24 months | 4 m/13 ft | 25/Box |
| 704-3821-00 | Pediatric w/O2 | 24 months | 2 m/7 ft  | 25/Box |



#### NomoLine LH Nasal/Oral CO2 Cannula


|             |                |           |           |        |
|-------------|----------------|-----------|-----------|--------|
| 704-3824-00 | Adult w/O2     | 24 months | 2 m/7 ft  | 25/Box |
| 704-4628-00 | Adult w/O2     | 24 months | 4 m/13 ft | 25/Box |
| 704-3825-00 | Pediatric w/O2 | 24 months | 2 m/7 ft  | 25/Box |



#### NomoLine LH Single Nasal Prong CO2 Cannula

|             |       |           |          |        |
|-------------|-------|-----------|----------|--------|
| 704-3826-00 | Adult | 24 months | 2 m/7 ft | 25/Box |
|-------------|-------|-----------|----------|--------|

### NomoLine LH (Low Humidity) Single Patient Use Sampling Lines - Intubated Patients



#### NomoLine LH Airway Adapter Set

#### Shelf Life\*

|             |                 |           |           |        |
|-------------|-----------------|-----------|-----------|--------|
| 704-3814-00 | Adult/Pediatric | 36 months | 2 m/7 ft  | 25/Box |
| 704-4620-00 | Adult/Pediatric | 36 months | 4 m/13 ft | 25/Box |

### NomoLine LH (Low Humidity) Procedural Oxygen Masks



#### NomoLine LH Nasal CO2 Cannula

#### Shelf Life\*

|                            |                            |           |           |        |
|----------------------------|----------------------------|-----------|-----------|--------|
| 704-5213-00 <sup>USA</sup> | Adult, Medium O2 Flow Rate | 24 months | 3 m/10 ft | 30/Box |
| 704-5214-00 <sup>USA</sup> | Pediatric, Medium O2 Flow  | 24 months | 3 m/10 ft | 30/Box |
| 704-5215-00 <sup>USA</sup> | Adult, High O2 Flow Rate   | 24 months | 3 m/10 ft | 30/Box |

\* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA



## NomoLine Sampling Lines

Masimo NomoLine technology eliminates common problems associated with conventional sidestream gas analysis. Designed for low-flow applications, NomoLine sampling lines can be used in a variety of clinical scenarios on both intubated and non-intubated adult, pediatric, infant, and neonatal patients, in both low- and high-humidity configurations. These are designed for use with Spacelabs Capnography and Multi-Gas Modules 92516, 92517 and 92518.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|



### NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Non-Intubated Patients



| NomoLine HH Nasal CO2 Cannula |                | Shelf Life* |           |        |
|-------------------------------|----------------|-------------|-----------|--------|
| 704-3830-00                   | Adult          | 24 months   | 2 m/7 ft  | 25/Box |
| 704-3832-00                   | Infant         | 24 months   | 2 m/7 ft  | 25/Box |
| 704-3833-00                   | Adult w/O2     | 24 months   | 2 m/7 ft  | 25/Box |
| 704-3834-00                   | Pediatric w/O2 | 24 months   | 2 m/7 ft  | 25/Box |
| 704-4933-00 <sup>USA</sup>    | Infant w/O2    | 24 months   | 4 m/13 ft | 25/Box |




| NomoLine HH Nasal/Oral CO2 Cannula |                | Shelf Life* |           |        |
|------------------------------------|----------------|-------------|-----------|--------|
| 704-3837-00                        | Adult w/O2     | 24 months   | 2 m/7 ft  | 25/Box |
| 704-4968-00 <sup>USA</sup>         | Adult w/O2     | 24 months   | 4 m/13 ft | 25/Box |
| 704-3838-00                        | Pediatric w/O2 | 24 months   | 2 m/7 ft  | 25/Box |



| NomoLine HH Single Nasal Prong CO2 Cannula |       | Shelf Life* |          |        |
|--------------------------------------------|-------|-------------|----------|--------|
| 704-3839-00                                | Adult | 24 months   | 2 m/7 ft | 25/Box |

### NomoLine HH (High Humidity) Single Patient Use Sampling Lines - Intubated Patients



| NomoLine HH Airway Adapter Set |                 | Shelf Life* |           |        |
|--------------------------------|-----------------|-------------|-----------|--------|
| 704-3827-00                    | Adult/Pediatric | 36 months   | 2 m/7 ft  | 25/Box |
| 704-3828-00                    | Adult/Pediatric | 36 months   | 3 m/10 ft | 25/Box |
| 704-4989-00 <sup>USA</sup>     | Adult/Pediatric | 36 months   | 4 m/13 ft | 25/Box |
| 704-4367-00                    | Infant/Neonatal | 36 months   | 2 m/7 ft  | 25/Box |

\* - Shelf life is from date of manufacture, printed on product label as an expiration date

USA - Parts only available for sale in the USA

## Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

### Sidestream Monitoring Supplies

**Please note:** the NomoLine multi-patient-use adapter has been discontinued. Without this adapter there is no means to connect the below sidestream supplies to our capnography modules. We do, however, make these supplies available for order for those customers who have remaining supply of the adapter, or use these for other applications.



#### CO2 Nasal Cannula w/O2 and luer

|    |             |           |           |           |        |
|----|-------------|-----------|-----------|-----------|--------|
| *  | 704-0181-00 | Adult     | 60 months | 2 m/7 ft  | 25/Box |
| *  | 704-0180-00 | Pediatric | 60 months | 2 m/7 ft  | 25/Box |
| ** | 704-0179-00 | Infant    | 60 months | 2 m/7 ft  | 25/Box |
| ** | 704-0186-00 | Adult     | 60 months | 3 m/10 ft | 25/Box |

#### CO2 Oral/Nasal Cannula w/O2 and luer

|    |             |       |           |           |        |
|----|-------------|-------|-----------|-----------|--------|
| *  | 704-0183-00 | Adult | 60 months | 2 m/7 ft  | 25/Box |
| CE | 704-0183-01 | Adult | 36 months | 2 m/7 ft  | 25/Box |
| ** | 704-0187-00 | Adult | 60 months | 3 m/10 ft | 25/Box |

Indicates discontinued part, available while supplies last

\* - Shelf life is from date of manufacture, printed on product label as an expiration date

### Sidestream Monitoring Supplies - Intubated patients



#### Airway adapters and lines

#### Shelf Life\*

|    |             |                         |           |           |        |
|----|-------------|-------------------------|-----------|-----------|--------|
| ** | 704-0191-00 | Adult Line w/ T-Adapter | 60 months | 2 m/7 ft  | 25/Box |
| *  | 166-7121-00 | sampling line           | 60 months | 3 m/10 ft | 50/Box |
| *  | 704-0175-00 | T-Adapter               | 60 months |           | 25/Box |

\* Product not available for sale into CE Mark Countries (EU)

CE product is available for sale into CE Mark Countries (EU)

\*\* Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

## Capnography & Multigas

Spacelabs carries a variety of TruLink and brand name capnography supplies that have been validated for use with Spacelabs modules. For Sidestream and Mainstream monitoring, we offer a variety of low-cost but high-function adapters, cannula, sampling lines, sensors, and airway T-adapters.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

**Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients.**

Featuring Masimo IRMA mainstream technology



### Sensor & Airway Adapters - required for mainstream monitoring

|             |                                                                                                          |                 |        |
|-------------|----------------------------------------------------------------------------------------------------------|-----------------|--------|
| 010-1980-00 | Masimo IRMA mainstream sensor, reusable                                                                  |                 | 1 Each |
| 704-0173-00 | Masimo IRMA mainstream airway adapter, single use<br><i>36 month shelf life from date of manufacture</i> | Adult/Pediatric | 25/Box |



(91517)  
Dual Capnography Module

## Capnography Accessories

The legacy dual capnography module 91517 is now obsolete. However, **Spacelabs will continue to carry TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module until December 31, 2026.** If you need a large last-time buy of these supplies, please have your orders in before this date.

|                                                                                     | Order Code                                                        | Description                                                                                        | Size                          | Qty    |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|-------------------------------|--------|
| <b>Sidestream, Non-intubated, single patient use</b>                                |                                                                   |                                                                                                    |                               |        |
|   | <b>CO<sub>2</sub> Nasal Cannula</b>                               |                                                                                                    | <b>Shelf Life<sup>^</sup></b> |        |
|                                                                                     | * 704-0010-00                                                     | Adult                                                                                              | 36 months                     | 10/Box |
|                                                                                     | * 704-0012-00                                                     | Pediatric                                                                                          | 36 months                     | 10/Box |
|                                                                                     | <b>CO<sub>2</sub> Nasal Cannula w/O<sub>2</sub> Delivery</b>      |                                                                                                    |                               |        |
|                                                                                     | * 704-0014-00                                                     | Adult                                                                                              | 36 months                     | 10/Box |
|                                                                                     | * 704-0016-00                                                     | Pediatric                                                                                          | 36 months                     | 10/Box |
|                                                                                     | * 704-0015-00                                                     | Infant/Neonatal                                                                                    | 36 months                     | 10/Box |
|                                                                                     | <b>CO<sub>2</sub> Oral/Nasal Cannula</b>                          |                                                                                                    |                               |        |
|                                                                                     | * 704-0017-00                                                     | Adult                                                                                              | 36 months                     | 10/Box |
|                                                                                     | <b>CO<sub>2</sub> Oral/Nasal Cannula w/O<sub>2</sub> Delivery</b> |                                                                                                    |                               |        |
|                                                                                     | * 704-0020-00                                                     | Adult                                                                                              | 36 months                     | 10/Box |
|                                                                                     | * 704-0022-00                                                     | Pediatric                                                                                          | 36 months                     | 10/Box |
| <b>Sidestream, Intubated, single patient use</b>                                    |                                                                   |                                                                                                    |                               |        |
|  | <b>Sidestream - short term up to 12 hours</b>                     |                                                                                                    |                               |        |
|                                                                                     | * 704-0023-00                                                     | Airway adapter + sampling line, ET tube > 4.0 mm                                                   | 36 months                     | 10/Box |
|                                                                                     | CE 704-0023-10                                                    | Airway adapter + sampling line, ET tube > 4.0 mm                                                   | 36 months                     | 10/Box |
|                                                                                     | * 704-0024-00                                                     | Airway adapter + sampling line, Pediatric/Infant sized                                             | 36 months                     | 10/Box |
|                                                                                     | * 704-0027-00                                                     | CO <sub>2</sub> sampling line w/male luer<br><i>Recommended for heat moisture exchange filters</i> | 36 months                     | 10/Box |
|                                                                                     | CE 704-0027-10                                                    | CO <sub>2</sub> sampling line w/male luer<br><i>Recommended for heat moisture exchange filters</i> | 36 months                     | 10/Box |

## Capnography Accessories

The legacy dual capnography module 91517 is now obsolete. However, **Spacelabs will continue to carry TruLink and brand name cannula, sampling lines, and airway T-adapters for use with this module until December 31, 2026.** If you need a large last-time buy of these supplies, please have your orders in before this date.



(91517)  
Dual Capnography Module

| Order Code                                    | Description                                                                                                   | Size                          | Qty    |
|-----------------------------------------------|---------------------------------------------------------------------------------------------------------------|-------------------------------|--------|
| <b>Sidestream - long term up to 120 hours</b> |                                                                                                               | <b>Shelf Life<sup>^</sup></b> |        |
| * 704-0025-00                                 | Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm                                                 | 36 months                     | 10/Box |
| CE 704-0025-10                                | Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm                                                 | 36 months                     | 10/Box |
| * 704-0026-00                                 | Airway adapter + sampling line with Nafion, Infant                                                            | 36 months                     | 10/Box |
| * 704-0028-00                                 | CO <sub>2</sub> sampling line w/male luer and Nafion<br><i>Recommended for heat moisture exchange filters</i> | 36 months                     | 10/Box |
| CE 704-0028-10                                | CO <sub>2</sub> sampling line w/male luer and Nafion<br><i>Recommended for heat moisture exchange filters</i> | 36 months                     | 10/Box |

## Mainstream, for Intubated Patients Only



### Mainstream Airway Adapters, for use with CO<sub>2</sub> Sensor 704-0001-00



|               |                                                   |           |        |
|---------------|---------------------------------------------------|-----------|--------|
| * 704-0005-00 | CO <sub>2</sub> airway adapter, single use, Adult | 36 months | 24/Box |
| * 704-0002-00 | CO <sub>2</sub> airway adapter, reusable, Adult   | 36 months | 1 Each |

**Indicates discontinued part, available while supplies last**

<sup>^</sup> - Shelf life is from date of manufacture, printed on product label as an expiration date

\* Product not available for sale into CE Mark Countries (EU)

**CE** product is available for sale into CE Mark Countries (EU)

\*\* Product not available for sale into CE Mark Countries (EU). NO alternative product currently available.

Ready to order?  Click here to shop online  
USA residents only

MULTIGAS ANALYZER  
Modules 90518 & 91518

Multigas Analyzer Supplies

For use with Ultraview SL™ Multigas Analyzers  
(90518 & 91518).

These two legacy multi-gas modules are obsolete. There is limited quantities of supplies available.



(90518)



(91518)

| Order Code | Description | Length/Size | Qty |
|------------|-------------|-------------|-----|
|------------|-------------|-------------|-----|

Multigas Supplies & Accessories



Water Traps, Reusable

|             |                                                                                                        |         |
|-------------|--------------------------------------------------------------------------------------------------------|---------|
| 717-0447-00 | D-FEND water trap cartridge, Black.<br>Replace at least once every two months or when occluded (91518) | 10/Pack |
| 717-0448-00 | D-FEND water trap cartridge, Green.<br>Replace at least once every 24 hours or when occluded (91518)   | 10/Pack |

*Note: Discontinued part, available while supplies last.*