

Elance

VITAL SIGNS MONITOR



Elace Vital Signs Monitor

Most TruLink supplies work with the Elance monitor:

- TruLink NIBP cuffs and hoses are compatible
- SpO2 sensors are compatible with same brand cable below
- TruLink ECG cables and leadwire systems are compatible

The below are unique to the Elance monitor.

Note, the Elance monitor had a last shipment date of March 31, 2020. Spacelabs will continue to provide the below supplies, as reasonably available, until March 31, 2027.



**Elace Vital Signs Monitor
93300**

Order Code	Description	Qty
Supplies unique to the Elance monitor		
ECG Lead Wires & Cables (can also use standard TruLink ECG cables and leadwire systems)		
ECG Integrated Cables and Lead Wires		
 700-0008-23	5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
700-0008-24	3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft 1 Each
Spacelabs TruLink SpO2 Technology (use cable below with TruLink branded sensors only)		
SpO2 Cable		
 700-0287-00	SpO2 adapter cable	305 cm/10 ft 1 Each
LNCS Masimo SET SpO2 Technology (use cable below with Masimo LNCS sensors only)		
LNCS Adapter Cable		
 690-0256-00	SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft 1 Each
Nellcor OxiMax SpO2 Technology (use cable below with Nellcor OxiMax sensors only)		
Cable		
 700-0290-00	SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft 1 Each
Invasive Blood Pressure (See IBP section for compatible sensors and transducers)		
Interface Cable		
 690-0021-00	Compatible with Datex Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. See IBP section for other compatible supplies.	366 cm/12 ft 1 Each
Indicates discontinued part, available while supplies last		