



Invasive Blood Pressure



Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers.



	Order Code	Description		Length/Size	Qty
Integrated Adapter Cables (connects directly to transducer)					
	Edward Lifesciences				
	015-0727-00	Single pressure integrated cable		400 cm/13 ft	1 Each
	015-0728-00	Dual pressure integrated cable		400 cm/13 ft	1 Each
	Abbott Laboratories/Hospira/ICU Medical				
	015-0729-00	Single pressure integrated cable		400 cm/13 ft	1 Each
015-0730-00	Dual pressure integrated cable	400 cm/13 ft		1 Each	
Interface Cable & Adapters					
	Interface Cable (Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers)				
	690-0021-00	Requires adapter cable 700-0091-00 or 700-0028-00		366 cm/12 ft	1 Each
	Interface Adapter- Transducer interface cable required from transducer company				
	700-0791-00	Single pressure cable		122 cm/4 ft	1 Each
	700-0028-00	Dual pressure cable		122 cm/4 ft	1 Each
Transducer Kits & Accessories					
	Pressure Transducer Kits, Single use				
	690-0008-00	Transducer Kit, Patient Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>		152.5 cm/60 in	5/Box
	690-0009-00	Transducer Kit, Pole Mount, DTX/Plus <i>Packaged sterile with 36 month shelf life</i>		152.5 cm/60 in	5/Box

Invasive Pressure Supplies

Supplies are available through Spacelabs for 90496/91496 modules configured for invasive pressure. Customers can choose from Edwards Lifesciences or Abbott Laboratories integrated single or dual pressure cables that connect directly to the transducer(s). For BD models, Spacelabs offers single or dual adapters and interface cables to allow use of compatible transducers.

